

MaximBio

Plant Disease Diagnoses

Viruses

Fungi

Growth Hormone

Viroid

Phytoplasma

Insect ID

Bacteria

GMO

Nematodes



ELISA Assay Kit

Pathosreen Kit

ELISA assay에 필요한 모든 항체와 시약이 제공됩니다.

Antibody-coated 96-well microtiter plates

Enzyme conjugated antibody

PNP substrate tablets, 5 mg each

PNP substrate buffer, 5X concentrate

ECl buffer, 5X concentrate

Positive control (if available)

Reagent Set

ELISA에 필요한 antibody 가 제공됩니다.

Coating antibody & conjugated antibody

Buffer pack

ELISA에 필요한 시약이 제공됩니다.

Carbonate Coating Buffer

PBST wash buffer (powder)

Extraction buffer, ECl buffer

Tween-20, 3M sodium hydroxide


PNP Substrate buffer, PNP Substrate tablets



ImmunoStrips


식물병 진단시약 - 진단막대에 식물 extract solution만 묻히면 바로 식물병 감염 여부를 알 수 있습니다.

- Sample extraction buffer는 따로 판매하고 있습니다.



Acidovorax avenae subsp. *citrulli*
Calibrachoa Mottle Virus
Clavibacter michiganensis subsp. *michiganensis*
Cucumber Mosaic Virus
Impatiens Necrotic Spot Virus
Pepino Mosaic Virus
Ralstonia solanacearum
Squash Mosaic Virus
Tobacco Mosaic Virus
Tomato Spotted Wilt Virus
Xanthomonas campestris pv. *pelargonii*
Bt-Cry1Ab/1Ac
Bt-Cry2A
Bt-Cry3A
Bt-Cry9C
Roundup Ready

Enzyme Linked ImmunoSorbent Assay (ELISA)




DAS ELISA

(Direct double antibody sandwich)

Capture by: polyclonal antibody


Detection by: polyclonal antibody



Compound direct ELISA

Capture by: polyclonal antibody


Detection by: monoclonal antibody mixture



Indirect ELISA (PTA ELISA)

Capture by: polystyrene plate

Detection by: monoclonal antibody mixture




TAS ELISA

(Indirect triple antibody)

Capture by: polyclonal antibody

Detection by: monoclonal antibody



IF

(immunofluorescence assay)

Capture by: polystyrene plate

Detection by: polyclonal antibody

Pab : polyclonal antibody
 Mab : monoclonal antibody
 Ag : antigen
 Ab : antibody

[A] 72

Abutilon Mosaic Virus
 African Cassava Mosaic Virus
 Alfalfa Mosaic Virus
 Alstroemeria Carla Virus
 Alstroemeria Flower Banding Virus
 Alstroemeria Mosaic Virus
 American Plum Line Pattern Virus
 Andean Potato Latent Virus
 Andean Potato Mottle Virus
 Apple Chlorotic Leaf Spot Virus
 Apple Mosaic Virus
 Apple Stem Grooving Virus
 Arabis Mosaic Virus
 Arracacha Virus B - OCA strain
 Asparagus Virus 2

[B] 74

Banana Bract Mosaic Virus
 Banana Bunchy Top Virus
 Banana Streak Virus
 Barley Mild Mosaic Virus
 Barley Stripe Mosaic Virus
 Barley Yellow Dwarf Virus
 Barley Yellow Mosaic Virus
 Bean Common Mosaic Necrosis Virus
 Bean Common Mosaic Virus
 Bean Golden Mosaic Virus
 Bean Leaf Roll Virus
 Bean Pod Mottle Virus
 Bean Yellow Mosaic Virus
 Beet Mild Yellowing Virus
 Beet Mosaic Virus
 Beet Necrotic Yellow Vein Virus
 Beet Soil Born Virus
 Beet Western Yellows Virus
 Beet Yellows Virus
 Bell Pepper Mottle Virus
 Belladonna Mottle Virus
 Blackeye Cowpea Mosaic Virus
 Blackgram Mottle Virus
 Blueberry Leaf Mottle Virus
 Blueberry Necrotic Shock Virus
 Blueberry Scorch Virus
 Blueberry Shock Virus
 Blueberry Shoestring Virus
 Broad Bean Stain Virus
 Broad Bean True Mosaic Virus
 Broad Bean Wilt Virus
 Brome Mosaic Virus
 Brome Streak Mosaic Virus

[C] 80

Cabbage Blackring Virus
 Cactus Virus X
 Calibrachoa Mottle Virus
 Carnation Etched Ring Virus
 Carnation Italian Ringspot Virus
 Carnation Latent Virus
 Carnation Mottle Virus
 Carnation Necrotic Fleck Virus
 Carnation Ringspot Virus
 Carnation Vein Mottle Virus
 Cauliflower Mosaic Virus
 Celery Mosaic Virus
 Cereal Yellow Dwarf Virus - rpv
 Cherry Leaf Roll Virus
 Chilli Veinal Mottle Virus
 Chrysanthemum Aspermy Virus
 Chrysanthemum Stem Necrosis Virus
 Chrysanthemum Virus B
 Citrus Tatter Leaf Virus
 Citrus Tristeza Virus
 Clover Yellow Mosaic Virus
 Clover Yellow Vein Virus
 Cocksfoot Mild Mosaic Virus
 Cowpea Aphid-Borne Mosaic Virus
 Cowpea Chlorotic Mottle Virus
 Cowpea Mild Mottle Virus
 Cowpea Mottle Virus
 Cowpea Yellow Mosaic Virus
 Cowpea Mosaic Virus
 Cowpea Severe Mosaic Virus
 Cucumber Green Mottle Mosaic Virus
 Cucumber Mosaic Virus
 Cucumber Necrosis Virus
 Cymbidium Mosaic Virus
 Cymbidium Ringspot Virus

[D] 87

Dahlia Mosaic Virus
 Dasheen Mosaic Virus

[E] 87

Eggplant Mosaic Virus
 Eggplant Mottled Dwarf Virus
 Elderberry Latent Virus

[F] 87

Flower - Break Virus
 Freesia Mosaic Virus

[G] 87

Garlic Common Latent Virus
 Garlic Virus D

Gemini Virus Screen	
Grapevine Algerian Latent Virus	
Grapevine Court-Noue Disease	
Grapevine Fanleaf Virus	
Grapevine Fleck Virus	
Grapevine Leaf Roll associated Virus	
Grapevine Virus A	
Grapevine Vitivirus B	
Groundnut Bud Necrosis Virus	
Groundnut Ringspot Virus	
[H]	89
Helenium Virus S	
Hibiscus Chlorotic Ringspot Virus	
Hop American Latent Virus	
Hop Latent Virus	
Hop Mosaic Virus	
Hosta Virus X	
Hydrangea Ringspot Virus	
[I]	90
Impatiens Necrotic Spot Virus	
Indian Cassava Mosaic Virus	
Iris Fulva Mosaic Virus	
Iris Mild Mosaic Virus	
Iris Severe Mosaic Virus	
Iris Yellow Spot Virus	
[J]	91
Jhonson Grass Chlorotic Strip Mosaic Virus	
Jhonson Grass Mosaic Virus	
[K]	92
Kalanchoe Latent Virus	
Kyuri Green Mottle Mosaic Virus	
[L]	92
Leek Yellow Stripe Virus	
Lettuce Mosaic Virus	
Lilac Chlorotic Leafspot Virus	
Lily Mottle Virus	
Lily Symptomless Virus	
Lily Virus X	
Line Pattern Virus	
[M]	93
Maize Chlorotic Dwarf Virus	
Maize Chlorotic Mottle Virus	
Maize Dwarf Mosaic Virus	
Maize Mosaic Virus	
Maize Rayado Fino Virus	
Maize Rough Dwarf Virus	
Maize Streak Virus	
Maize Stripe Virus	
Maize White Line Mosaic Virus	
Maracuja Mosaic Virus	
Moroccan Watermelon Mosaic Virus	
Melon Necrotic Spot Virus	
Mirafiori Lettuce Big-Vein Virus	
Mirafiori Lettuce Mosaic Virus	
Mirafiori Lettuce Virus	
Mushroom Virus 1	
[N]	95
Nandina Mosaic Virus	
Narcissus Latent Virus	
Narcissus Late Season Yellows Virus	
Narcissus Mosaic Virus	
Narcissus Q Virus	
Narcissus Tip Necrosis Virus	
Narcissus Yellow Stripe Virus	
Narcissus Virus Set	
Neckar River Virus	
Nerine Latent Virus	
Nerine Virus X	
Nerine Virus Y	
Nerine Yellow Stripe Virus	
[O]	96
Odontoglossum Ringspot Virus	
Okra Mosaic Virus	
Onion Yellow Dwarf Virus	
[P]	97
Papaya Mosaic Virus	
Papaya Ringspot Virus	
Paprika Mild Mottle Virus	
Passiflora Latent Virus	
Pea Early-Browning Virus	
Pea Enation Mosaic Virus	
Pea Mosaic Virus	
Pea Seed-Borne Mosaic Virus	
Peach Rosette Mosaic Virus	
Peanut Mottle Virus	
Peanut Stunt Virus	
Pelargonium Flower Break Virus	
Pelargonium Leaf Curl Virus	
Pelargonium Line Pattern Virus	
Pelargonium Ringspot Virus	
Pelargonium Zonate Spot Virus	
Pepino Mosaic Virus	
Pepino Mosaic Virus	
Pepper Mild Mottle Virus	
Pepper Mild Mottle Virus	
Pepper Mottle Virus	
Pepper Veinal Mottle Virus	
Peru Tomato Mosaic Virus	
Petunia Asteroid Mosaic Virus	
Pineapple Mealybug Wilt associated Virus - 1	
Pineapple Mealybug Wilt associated Virus - 2	
Plum Pox Virus	
Poinsettia Cryptic Virus	
Poinsettia Mosaic Virus	
Poplar Mosaic Virus	

Potato Andean Mottle Virus
 Potato Aucuba Mosaic Virus
 Potato Black Ringspot Virus
 Potato Latent Virus
 Potato Leaf Roll Virus
 Potato Mop Top Virus
 Potato Virus A
 Potato Virus B - OCA strain
 Potato Virus M
 Potato Virus S
 Potato Virus T
 Potato Virus V
 Potato Virus X
 Potato Virus Y
 Potato Yellow Dwarf Virus
 Potato Yellowing Virus
 Potyvirus Group
 Prune Dwarf Virus
 Prunus Necrotic Ringspot Virus

[R] **105**

Radish Mosaic Virus
 Raspberry Bushy Dwarf Virus
 Raspberry Ringspot Virus
 Red Clover Cryptic Virus 2
 Red Clover Mottle Virus
 Red Clover Vein Mosaic Virus
 Red LaSoda Virus
 Ribgrass Mosaic Virus
 Rice Dwarf Virus
 Rice Ragged Stunt Virus
 Rice Stripe Virus
 Rice Yellow Mottle Virus
 Ryegrass Cryptic Virus
 Ryegrass Mosaic Virus

[S] **107**

Saguaro Cactus Virus
 Scrophularia Mottle Virus
 Satsuma Dwarf Virus
 Shallot Latent Virus
 Shallot Virus X
 Shallot Yellow Stripe Virus
 Sieg River Virus
 Soil-borne Cereal Mosaic Virus
 Soil-borne Wheat Mosaic Virus
 Sonchus Yellow Net Virus
 Sorghum Mosaic Virus
 Southern Bean Mosaic Virus
 Soybean Mosaic Virus
 Squash Leaf Curl Virus
 Squash Mosaic Virus
 Strawberry Latent Ringspot Virus
 Strawberry Mild Yellow Edge Virus
 Strawberry Mild Yellow Edge Potex Virus
 Streptocarpus Flower-Break Virus

Sugarcane Bacilliform Virus
 Sugarcane Mosaic Virus
 Sweetpotato NCM- ELISA Kits

[T] **111**

Tobacco Etch Virus
 Tobacco Mild Green Mosaic Virus
 Tobacco Mosaic Virus
 Tobacco Necrosis Virus
 Tobacco Rattle Virus
 Tobacco Ringspot Virus
 Tobacco Streak Virus
 Tobacco Vein Banding Mosaic Virus
 Tobacco Vein Mottling Virus
 Tomato Aspermy Virus
 Tomato Black Ring Virus
 Tomato Bushy Stunt Virus
 Tomato Chlorotic Spot Virus
 Tomato Mosaic Virus
 Tomato Mottle Virus
 Tomato Ringspot Virus
 Tomato Spotted Wilt Virus
 Tomato Yellow Leaf Curl Virus
 Gemini Virus Screen
 (Indian, European, African strains of TYLCV)
 Tospo Virus
 Tulip Breaking Virus
 Tulip Crinkle Virus
 Tulip Virus X
 Turnip Mosaic Virus
 Turnip Yellow Mosaic Virus
 Turnip Yellows Virus

[V] **118**

Viola Mottle Virus

[W] **118**

Watermelon Fruit Blotch Virus
 Watermelon Mosaic Virus 1
 Watermelon Mosaic Virus 2
 Watermelon Sliver Mottle Virus
 Wheat Dwarf Virus
 Wheat Soil-Borne Mosaic Virus
 Wheat Spindle Streak Mosaic Virus
 Wheat Streak Mosaic Virus
 Wheat Yellow Mosaic Virus
 Whitefly-Transmitted (Begomoviruses)
 Geminiviruses
 White Clover Cryptic Virus 1
 White Lupine Mosaic Virus
 Wild Potato Mosaic Virus

[Y] **119**

Yam Mosaic Virus

[Z] **119**

Zucchini Green Mottle Mosaic Virus
 Zucchini Yellow Fleck Virus
 Zucchini Yellow Mosaic Virus

Viroid **120**

Chrysanthemum Chlorotic Mottle Viroid
 Chrysanthemum Stunt Viroid
 Potato Spindle Tuber Viroid

Bacteria **120**

Acidovorax avenae subsp. citrulli
 Agrobacterium radiobacter subsp. tumefaciens
 Agrobacterium tumefaciens (strain-specific)
 Banana Blood Disease
 Burkholderia gladioli pv. allicola
 Burkholderia glumae
 Burkholderia cepacia
 Cladosporium spp.
 Clavibacter michiganensis subsp. insidiosum
 Clavibacter michiganensis subsp. michiganensis
 Clavibacter michiganensis subsp. sepedonicus
 Clavibacter michiganensis subsp. sepedonicus-Cy3
 Clavibacter michiganensis subsp. sepedonicus-FITC
 Clavibacter michiganensis subsp. tessellarius
 Curtobacterium flaccumfaciens pv. flaccumfaciens
 Curtobacterium flaccumfaciens pv. oortii
 Erwinia amylovora
 Erwinia carotovora
 Erwinia carotovora subsp. atroseptica
 Erwinia carotovora subsp. carotovora
 Erwinia chrysanthemi
 Erwinia chrysanthemi pv. dianthicola
 Erwinia herbicola
 Erwinia salicis
 Erwinia salicilis
 Erwinia stewartii
 Pantoea stewartii
 Pantoea stewartii subsp. stewartii
 Pectobacterium chrysanthemi
 Pierce's Disease
 Pseudomonas avenae
 Pseudomonas andropogonis
 Pseudomonas cichorii
 Pseudomonas fuscovaginae
 Pseudomonas gladioli pv. allicola
 Pseudomonas glumae
 Pseudomonas marginalis pv. marginalis

Pseudomonas savastanoi pv. phaseolicola
 Pseudomonas solanacearum - Ps1
 Pseudomonas solanacearum
 Pseudomonas syringae pv. avellanae
 Pseudomonas syringae pv. glycinea
 Pseudomonas syringae pv. lachrymans
 Pseudomonas syringae pv. mors-prunorum
 Pseudomonas syringae pv. phaseolicola
 Pseudomonas syringae pv. pisi
 Pseudomonas syringae pv. porri
 Pseudomonas syringae pv. syringae
 Pseudomonas syringae pv. tomato
 Pseudomonas syringae subsp. savastanoi
 Pseudomonas viridiflava
 Ralstonia solanacearum
 Rhodococcus fascians
 Spongospora subterranea f. sp. subterranea
 Spiroplasma citri
 Xanthomonas
 Xanthomonas albilineans
 Xanthomonas arboricola pv. pruni
 Xanthomonas axonopodis pv. arboricola
 Xanthomonas axonopodis pv. begoniae
 Xanthomonas axonopodis pv. citri
 Xanthomonas campestris
 Xanthomonas campestris pv. armoraciae
 Xanthomonas campestris pv. begoniae
 Xanthomonas campestris pv. campestris
 Xanthomonas campestris pv. carotae
 Xanthomonas campestris pv. citri
 Xanthomonas campestris pv. dieffenbachiae
 Xanthomonas campestris pv. hyacinthi
 Xanthomonas campestris pv. jacinti
 Xanthomonas campestris pv. oryzae
 Xanthomonas campestris pv. pelargonii
 Xanthomonas campestris pv. phaseoli
 Xanthomonas campestris pv. pruni
 Xanthomonas campestris pv. vesicatoria
 Xanthomonas fragariae
 Xanthomonas hortorum pv. pelargonii
 Xanthomonas maltophilia
 Xylella fastidiosa
 Xylophilus ampelinus (bacterial necrosis)

Fungi **113**

Ascochyta pisi
 Aspergillus & penicillium spp.
 Botrytis allii
 Botrytis cinerea
 Cladosporium spp.
 Colletotrichum acutatum
 Colletotrichum spp.
 Ergosterol (Total Fungus detection)
 Humicola lanuginosa
 Laccase

Monilinia fructicola
 Monilinia spp.
 Mouldetect
 Mycosphasrella pinodes
 Nectria galligena
 Phoma medicaginis pinodella
 Phytophthora
 Phytophthora infestans
 Phythium species.
 Phytophthora spp.
 Pseudocercospora (Eyespot)
 Pythium ultimum
 Pythium spp.
 Rhizoctonia solani I
 Rhizoctonia solani II
 Rhizoctonia spp.
 Rhizoctonia species.
 Phytophthora+Pythium+Rhizoctonia species.
 Sclerotinia sclerotiorum
 Septoria nodorum+ Septoria tritici
 Septoria nodorum
 Septoria tritici
 Tapesia spp. (Eyespot)
 Verticillium dahliae
 Verticillium spec. rape isolate
 Verticillium spec. strawberry isolate
 Verticillium spp.

Phytoplasma 137

Apple Proliferation Phytoplasma
 Aster Yellows Phytoplasma
 Corn Stunt Spiroplasma
 Elm Yellows Phytoplasma
 Flavescence doree Phytoplasma
 Stolbur Phytoplasma
 Universal Phytoplasma
 Universal Phytoplasma nested

GMO 138

Bt-Cry1Ab/1Ac
 Bt-Cry1C protein
 Bt-Cry1F protein
 Bt-Cry2A protein
 Bt-Cry3A protein
 Bt-Cry3Bb1 protein
 Bt-Cry9C protein
 Bacillus thuringiensis subsp. kurstaki
 Neomycin Phosphotransferase II
 Roundup Ready (CP4 EPSPS)

Insect ID 140

Helicoverpa/Heliothis - ID

Nematodes 140

Globodera rostochiensis/pallida

Plant Growth Hormones 141

Abscisic Acid
 Dihydrozeatin Riboside
 Indole-3-acetic Acid
 Isopentenyladenosine
 trans-Zeatin Riboside

Accessories 142

CHEMICALS AND BUFFERS

3M sodium hydroxide
 3M sulfuric acid
 Ames buffer
 APLPV EC buffer
 BCIP substrate
 Bovine serum albumin
 BP substrate solution, alkaline phosphatase
 Buffer pack
 Carbonate coating buffer
 Chicken egg albumin
 Diethyldithiocarbamate
 ECI buffer
 ECG buffer
 Extraction buffer
 MRS component
 Nonfat dried milk
 OPD substrate buffer / Sticks
 PBST wash buffer
 PCR quality water
 PNP substrate buffer
 PNP substrate tablets
 Polyvinylpyrrolidone, PVP-40T
 PZSV EC buffer
 TMB peroxidase substrate solution
 Tween 20

LABORATORY EQUIPMENT AND SUPPLIES

Canister Racks X
KLECO 4200
KLECO 8200
Microcentrifuge tubes / Holder
Micro dilution tubes
Micro dilution tubes, racked
Nylon membrane
Pipetter
Pipet Tips
Plate
Plate reader, automated
Plate Washer Bottle
Precoated strip plate
Roller Press
Sample bag
Sample cup holder tray
Sample cups
Screwcap vial
Sealed straws
Tissue Homogenizer

PCR primer for viral group 145

바이러스검정 RCR-Mix Kit 146

Antibody Service 147

Abutilon Mosaic Virus

ABMV	DAS ELISA	
KS 1001 - 0500R	Reagent Set	500 test
KS 1001 - 1000R	Reagent Set	1000 test
KS 1001 - 5000R	Reagent Set	5000 test
KS 1001 - P	Positive Control	
KS 1001 - N	Negative Control	

African Cassava Mosaic Virus

ACMV	DAS ELISA	Poly/Poly
CAB 73600/0500	Coating Ab	500 test
CAB 73600/1000	Coating Ab	1000 test
CAB 73600/5000	Coating Ab	5000 test
ECA 73600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 73600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 73600/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 73600/0500	Reagent Set,Alk.Phos.	500 test
SRA 73600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 73600/5000	Reagent Set,Alk.Phos.	5000 test

African Cassava Mosaic Virus

ACMV	TAS ELISA	
KS 1378 - 0500R	Reagent Set	500 test
KS 1378 - 1000R	Reagent Set	1000 test
KS 1378 - 5000R	Reagent Set	5000 test
KS 1378 - P	Positive Control	
KS 1378 - N	Negative Control	

Alfalfa Mosaic Virus

AMV	DAS ELISA	Poly/Mono
CAB 87600/0500	Coating Ab	500 test
CAB 87600/1000	Coating Ab	1000 test
CAB 87600/5000	Coating Ab	5000 test
ECA 87600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 87600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 87600/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 87600/0500	Conjugated Ab,Peroxidase	500 test
ECP 87600/1000	Conjugated Ab,Peroxidase	1000 test
ECP 87600/5000	Conjugated Ab,Peroxidase	5000 test
PSA 87600/0500	PathoScreen Kit,Alk.Phos.	96 test
PSA 87600/1000	PathoScreen Kit,Alk.Phos.	288 test
PSA 87600/5000	PathoScreen Kit,Alk.Phos.	480 test
PSP 87600/0500	PathoScreen Kit, Peroxidase	96 test
PSP 87600/1000	PathoScreen Kit, Peroxidase	288 test
PSP 87600/5000	PathoScreen Kit, Peroxidase	480 test
SRA 87600/0500	Reagent Set,Alk.Phos.	500 test
SRA 87600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 87600/5000	Reagent Set,Alk.Phos.	5000 test
SRP 87600/0500	Reagent Set,Peroxidase	500 test

0/1000	Reagent Set,Peroxidase	1000 test
0/5000	Reagent Set,Peroxidase	5000 test
0	Negative Control	
0	Positive Control	

emia Carla Virus

DAS ELISA		
0500R	Reagent Set	500 test
1000R	Reagent Set	1000 test
5000R	Reagent Set	5000 test
P	Positive Control	
N	Negative Control	

emia Flower Banding Virus

DAS ELISA		
0500R	Reagent Set	500 test
1000R	Reagent Set	1000 test
5000R	Reagent Set	5000 test
P	Positive Control	
N	Negative Control	

emia Mosaic Virus

	DAS ELISA	Poly/Poly
0/0500	Coating Ab	500 test
0/1000	Coating Ab	1000 test
0/5000	Coating Ab	5000 test
0/0500	Conjugated Ab,Alk.Phos.	500 test
0/1000	Conjugated Ab,Alk.Phos.	1000 test
0/5000	Conjugated Ab,Alk.Phos.	5000 test
0/0500	Reagent Set,Alk.Phos.	500 test
0/1000	Reagent Set,Alk.Phos.	1000 test
0/5000	Reagent Set,Alk.Phos.	5000 test

can Plum Line Pattern Virus

	DAS ELISA	Poly/Poly
0/0500	Coating Ab	500 test
0/1000	Coating Ab	1000 test
0/5000	Coating Ab	5000 test
0/0500	Conjugated Ab,Alk.Phos.	500 test
0/1000	Conjugated Ab,Alk.Phos.	1000 test
0/5000	Conjugated Ab,Alk.Phos.	5000 test
0/0500	Reagent Set,Alk.Phos.	500 test
0/1000	Reagent Set,Alk.Phos.	1000 test
0/5000	Reagent Set,Alk.Phos.	5000 test
0	Negative Control	
0	Positive Control	

n Potato Latent Virus

	DAS ELISA	Poly/Poly
0/0500	Coating Ab	500 test

CAB 40200/1000	Coating Ab	1000 test
CAB 40200/5000	Coating Ab	5000 test
ECA 40200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 40200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 40200/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 40200/0500	Reagent Set,Alk.Phos.	500 test
SRA 40200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 40200/5000	Reagent Set,Alk.Phos.	5000 test
LNC 40200	Negative Control	

LNC 30700	Negative Control	
-----------	------------------	--

Apple Stem Grooving Virus

ASGV	DAS ELISA	
KS 1012 - 0500R	Reagent Set	500 test
KS 1012 - 1000R	Reagent Set	1000 test
KS 1012 - 5000R	Reagent Set	5000 test
KS 1012 - P	Positive Control	
KS 1012 - N	Negative Control	

Arabis Mosaic Virus

ArMV	DAS ELISA	Poly/Poly
CAB 23201/0500	Coating Ab	500 test
CAB 23201/1000	Coating Ab	1000 test
CAB 23201/5000	Coating Ab	5000 test
ECA 23201/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 23201/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 23201/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 23201/0500	Conjugated Ab,Peroxidase	500 test
ECP 23201/1000	Conjugated Ab,Peroxidase	1000 test
ECP 23201/5000	Conjugated Ab,Peroxidase	5000 test
PSA 23201/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 23201/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 23201/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 23201/0096	PathoScreen Kit,Peroxidase	96 test
PSP 23201/0288	PathoScreen Kit,Peroxidase	288 test
PSP 23201/0480	PathoScreen Kit,Peroxidase	480 test
SRA 23201/0500	Reagent Set,Alk.Phos.	500 test
SRA 23201/1000	Reagent Set,Alk.Phos.	1000 test
SRA 23201/5000	Reagent Set,Alk.Phos.	5000 test
SRP 23201/0500	Reagent Set,Peroxidase	500 test
SRP 23201/1000	Reagent Set,Peroxidase	1000 test
SRP 23201/5000	Reagent Set,Peroxidase	5000 test
LNC 23201	Negative Control	
LPC 23201	Positive Control	

Arabis Mosaic Virus in grape

ArMV	TAS ELISA	Poly/Poly
XAB 23202/0500	Coating Ab	500 test
XAB 23202/1000	Coating Ab	1000 test
XAB 23202/5000	Coating Ab	5000 test
XCA 23202/0500	Conjugated Ab,Alk.Phos.	500 test
XCA 23202/1000	Conjugated Ab,Alk.Phos.	1000 test
XCA 23202/5000	Conjugated Ab,Alk.Phos.	5000 test
XRA 23202/0500	Conjugated Ab,Peroxidase	500 test
XRA 23202/1000	Conjugated Ab,Peroxidase	1000 test
XRA 23202/5000	Conjugated Ab,Peroxidase	5000 test
LNC 23202	Negative Control	

Andean Potato Mottle Virus

APMoV	DAS ELISA	
KS 1009 - 0500R	Reagent Set	500 test
KS 1009 - 1000R	Reagent Set	1000 test
KS 1009 - 5000R	Reagent Set	5000 test
KS 1009 - P	Positive Control	
KS 1009 - N	Negative Control	

Apple Chlorotic Leaf Spot Virus

ACLSV	DAS ELISA	
KS 1010 - 0500R	Reagent Set	500 test
KS 1010 - 1000R	Reagent Set	1000 test
KS 1010 - 5000R	Reagent Set	5000 test
KS 1010 - P	Positive Control	
KS 1010 - N	Negative Control	

Apple Mosaic Virus

ApMV	DAS ELISA	Poly/Poly
CAB 30700/0500	Coating Ab	500 test
CAB 30700/1000	Coating Ab	1000 test
CAB 30700/5000	Coating Ab	5000 test
ECA 30700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 30700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 30700/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 30700/0500	Conjugated Ab,Peroxidase	500 test
ECP 30700/1000	Conjugated Ab,Peroxidase	1000 test
ECP 30700/5000	Conjugated Ab,Peroxidase	5000 test
PSA 30700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 30700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 30700/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 30700/0096	PathoScreen Kit,Peroxidase	96 test
PSP 30700/0288	PathoScreen Kit,Peroxidase	288 test
PSP 30700/0480	PathoScreen Kit,Peroxidase	480 test
SRA 30700/0500	Reagent Set,Alk.Phos.	500 test
SRA 30700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 30700/5000	Reagent Set,Alk.Phos.	5000 test
SRP 30700/0500	Reagent Set,Peroxidase	500 test
SRP 30700/1000	Reagent Set,Peroxidase	1000 test
SRP 30700/5000	Reagent Set,Peroxidase	5000 test

Arabid Mosaic Virus+Grapevine Fan Leaf Viruses

ArMV+GFLV	DAS ELISA	Poly/Poly
KS 1014 - 0500R	Reagent Set	500 test
KS 1014 - 1000R	Reagent Set	1000 test
KS 1014 - 5000R	Reagent Set	5000 test
KS 1014 - P	Positive Control	
KS 1014 - N	Negative Control	

Arracacha Virus B - OCA strain

AVBO	DAS ELISA	
KS 1015 - 0500R	Reagent Set	500 test
KS 1015 - 1000R	Reagent Set	1000 test
KS 1015 - 5000R	Reagent Set	5000 test
KS 1015 - P	Positive Control	
KS 1015 - N	Negative Control	

Asparagus Virus 2

AV2	Indirect ELISA	- /Poly
CAB 71000/0500	Coating Ab	500 test
CAB 71000/1000	Coating Ab	1000 test
CAB 71000/5000	Coating Ab	5000 test
ECA 71000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 71000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 71000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 71000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 71000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 71000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 71000/0500	Reagent Set,Alk.Phos.	500 test
SRA 71000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 71000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 71000	Negative Control	
LPC 71000	Positive Control	

Banana Bract Mosaic Virus

BBrMV	DAS ELISA	Poly/Poly
CAB 24600/0500	Coating Ab	500 test
CAB 24600/1000	Coating Ab	1000 test
CAB 24600/5000	Coating Ab	5000 test
ECA 24600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 24600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 24600/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 24600/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 24600/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 24600/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 24600/0500	Reagent Set,Alk.Phos.	500 test
SRA 24600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 24600/5000	Reagent Set,Alk.Phos.	5000 test
LNC 24600	Negative Control	
LPC 24600	Positive Control	

Banana Bunchy Top Virus

BBTV	Compound Direct ELISA	Poly/Mono
CAB 24700/0500	Coating Ab	500 test
CAB 24700/1000	Coating Ab	1000 test
CAB 24700/5000	Coating Ab	5000 test
ECA 24700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 24700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 24700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 24700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 24700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 24700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 24700/0500	Reagent Set,Alk.Phos.	500 test
SRA 24700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 24700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 24700	Negative Control	
LPC 24700	Positive Control	

Banana Streak Virus

BSV	DAS ELISA	
KS 1019 - 0500R	Reagent Set	500 test
KS 1019 - 1000R	Reagent Set	1000 test
KS 1019 - 5000R	Reagent Set	5000 test
KS 1019 - P	Positive Control	
KS 1019 - N	Negative Control	

Barley Mild Mosaic Virus

BaMMV	DAS ELISA	
KS 1020 - 0500R	Reagent Set	500 test
KS 1020 - 1000R	Reagent Set	1000 test
KS 1020 - 5000R	Reagent Set	5000 test
KS 1020 - P	Positive Control	
KS 1020 - N	Negative Control	

Barley Stripe Mosaic Virus

BSMV	DAS ELISA	Poly/Poly
CAB 19500/0500	Coating Ab	500 test
CAB 19500/1000	Coating Ab	1000 test
CAB 19500/5000	Coating Ab	5000 test
ECA 19500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 19500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 19500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 19500/0500	Conjugated Ab,Peroxidase	500 test
ECP 19500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 19500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 19500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 19500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 19500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 19500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 19500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 19500/0480	PathoScreen Kit,Peroxidase	480 test

SRA 19500/0500	Reagent Set,Alk.Phos.	500 test
SRA 19500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 19500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 19500/0500	Reagent Set,Peroxidase	500 test
SRP 19500/1000	Reagent Set,Peroxidase	1000 test
SRP 19500/5000	Reagent Set,Peroxidase	5000 test
LNC 19500	Negative Control	
LPC 19500	Positive Control	



Barley Yellow Dwarf Virus - mav

BYDV-mav		TAS ELISA
KS 1023 - 0500R	Reagent Set	500 test
KS 1023 - 1000R	Reagent Set	1000 test
KS 1023 - 5000R	Reagent Set	5000 test
KS 1023 - P	Positive Control	
KS 1023 - N	Negative Control	

Barley Yellow Dwarf Virus - pav

BYDV-pav		DAS ELISA
KS 1026 - 0500R	Reagent Set	500 test
KS 1026 - 1000R	Reagent Set	1000 test
KS 1026 - 5000R	Reagent Set	5000 test
KS 1026 - P	Positive Control	
KS 1026 - N	Negative Control	

Barley Yellow Dwarf Virus - pav

BYDV-pav		TAS ELISA
KS 1025 - 0500R	Reagent Set	500 test
KS 1025 - 1000R	Reagent Set	1000 test
KS 1025 - 5000R	Reagent Set	5000 test
KS 1025 - P	Positive Control	
KS 1025 - N	Negative Control	

Barley Yellow Dwarf Virus - pav

BYDV-pav		Indirect sandwich ELISA	Poly/Poly
CAB 27500/0500	Coating Ab		500 test
CAB 27500/1000	Coating Ab		1000 test
CAB 27500/5000	Coating Ab		5000 test
ECA 27500/0500	Conjugated Ab,Alk.Phos.		500 test
ECA 27500/1000	Conjugated Ab,Alk.Phos.		1000 test
ECA 27500/5000	Conjugated Ab,Alk.Phos.		5000 test
PSA 27500/0096	PathoScreen Kit,Alk.Phos.		96 test
PSA 27500/0288	PathoScreen Kit,Alk.Phos.		288 test
PSA 27500/0480	PathoScreen Kit,Alk.Phos.		480 test
SRA 27500/0500	Reagent Set,Alk.Phos.		500 test
SRA 27500/1000	Reagent Set,Alk.Phos.		1000 test
SRA 27500/5000	Reagent Set,Alk.Phos.		5000 test
LNC 27500	Negative Control		
LPC 27500	Positive Control		

Barley Yellow Dwarf Virus - mav+pav

BYDV-mav+pav		DAS ELISA
KS 1027 - 0500R	Reagent Set	500 test
KS 1027 - 1000R	Reagent Set	1000 test
KS 1027 - 5000R	Reagent Set	5000 test
KS 1027 - P	Positive Control	
KS 1027 - N	Negative Control	

Barley Yellow Dwarf Virus - B (See Barley Yellow Dwarf Virus - pav)

Barley Yellow Dwarf Virus - F (See Barley Yellow Dwarf Virus - mav)

Barley Yellow Dwarf Virus - mav

BYDV-mav	DAS ELISA	Mono/Poly
CAB 26500/0500	Coating Ab	500 test
CAB 26500/1000	Coating Ab	1000 test
CAB 26500/5000	Coating Ab	5000 test
ECA 26500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 26500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 26500/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 26500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 26500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 26500/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 26500/0500	Reagent Set,Alk.Phos.	500 test
SRA 26500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 26500/5000	Reagent Set,Alk.Phos.	5000 test
LNC 26500	Negative Control	
LPC 26500	Positive Control	



Barley Yellow Dwarf Virus - rpv

BYDV-rpv	DAS ELISA	
KS 1029 - 0500R	Reagent Set	500 test
KS 1029 - 1000R	Reagent Set	1000 test
KS 1029 - 5000R	Reagent Set	5000 test
KS 1029 - P	Positive Control	
KS 1029 - N	Negative Control	

Barley Yellow Dwarf Virus - rpv

BYDV-rpv	TAS ELISA	
KS 1030 - 0500R	Reagent Set	500 test
KS 1030 - 1000R	Reagent Set	1000 test
KS 1030 - 5000R	Reagent Set	5000 test
KS 1030 - P	Positive Control	
KS 1030 - N	Negative Control	

Barley Yellow Dwarf Virus - rpv (See Cereal Yellow Dwarf Virus - rpv)

Barley Yellow Dwarf Virus - sgv

BYDV-sgv	DAS ELISA	Poly/Poly
CAB 27000/0500	Coating Ab	500 test
CAB 27000/1000	Coating Ab	1000 test
CAB 27000/5000	Coating Ab	5000 test
ECA 27000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 27000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 27000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 27000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 27000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 27000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 27000/0500	Reagent Set,Alk.Phos.	500 test
SRA 27000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 27000/5000	Reagent Set,Alk.Phos.	5000 test

LNC 27000	Negative Control
LPC 27000	Positive Control

Barley Yellow Dwarf Virus screen (BYDV-mav + pav + rpv)

TAS ELISA

KS 1402 - 0500R	Reagent Set	500 test
KS 1402 - 1000R	Reagent Set	1000 test
KS 1402 - 5000R	Reagent Set	5000 test
KS 1402 - P	Positive Control	
KS 1402 - N	Negative Control	

Barley Yellow Mosaic Virus

BaYMV DAS ELISA

KS 1033 - 0500R	Reagent Set	500 test
KS 1033 - 1000R	Reagent Set	1000 test
KS 1033 - 5000R	Reagent Set	5000 test
KS 1033 - P	Positive Control	
KS 1033 - N	Negative Control	

Barley Yellow Mosaic Virus

BaYMV TAS ELISA

KS 1401 - 0500R	Reagent Set	500 test
KS 1401 - 1000R	Reagent Set	1000 test
KS 1401 - 5000R	Reagent Set	5000 test
KS 1401 - P	Positive Control	
KS 1401 - N	Negative Control	

Bean Common Mosaic Necrosis Virus

BCMNV DAS ELISA

KS 1034 - 0500R	Reagent Set	500 test
KS 1034 - 1000R	Reagent Set	1000 test
KS 1034 - 5000R	Reagent Set	5000 test
KS 1034 - P	Positive Control	
KS 1034 - N	Negative Control	

Bean Common Mosaic Necrosis Virus-Serotype A

BCMNV DAS ELISA

KS 1381 - 0500R	Reagent Set	500 test
KS 1381 - 1000R	Reagent Set	1000 test
KS 1381 - 5000R	Reagent Set	5000 test
KS 1381 - P	Positive Control	
KS 1381 - N	Negative Control	

Bean Common Mosaic Virus

BCMV DAS ELISA

KS 1035 - 0500R	Reagent Set	500 test
KS 1035 - 1000R	Reagent Set	1000 test
KS 1035 - 5000R	Reagent Set	5000 test
KS 1035 - P	Positive Control	

KS 1035 - N	Negative Control
-------------	------------------

Bean Common Mosaic Virus - Serotype B

BCMV	DAS ELISA	
KS 1382 - 0500R	Reagent Set	500 test
KS 1382 - 1000R	Reagent Set	1000 test
KS 1382 - 5000R	Reagent Set	5000 test
KS 1382 - P	Positive Control	
KS 1382 - N	Negative Control	

Bean Common Mosaic Virus

BCMV	Indirect ELISA	-- / Mono
CAB 46000/0500	Coating Ab	500 test
CAB 46000/1000	Coating Ab	1000 test
CAB 46000/5000	Coating Ab	5000 test
ECA 46000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 46000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 46000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 46000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 46000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 46000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 46000/0500	Reagent Set,Alk.Phos.	500 test
SRA 46000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 46000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 46000	Negative Control	
LPC 46000	Positive Control	

Bean Golden Mosaic Virus

BGMV	DAS ELISA	
KS 1321 - 0500R	Reagent Set	500 test
KS 1321 - 1000R	Reagent Set	1000 test
KS 1321 - 5000R	Reagent Set	5000 test
KS 1321 - P	Positive Control	
KS 1321 - N	Negative Control	

Bean Leaf Roll Virus

BLRV	DAS ELISA	
KS 1038 - 0500R	Reagent Set	500 test
KS 1038 - 1000R	Reagent Set	1000 test
KS 1038 - 5000R	Reagent Set	5000 test
KS 1038 - P	Positive Control	
KS 1038 - N	Negative Control	

Bean Pod Mottle Virus

BPMV	DAS ELISA	Poly/Poly
CAB 46400/5000	Coating Ab	5000 test
ECA 46400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 46400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 46400/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 46400/0096	PathoScreen Kit,Alk.Phos.	96 test

PSA 46400/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 46400/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 46400/0500	Reagent Set,Alk.Phos.	500 test
SRA 46400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 46400/5000	Reagent Set,Alk.Phos.	5000 test
LNC 46400	Negative Control	
LPC 46400	Positive Control	

Bean Yellow Mosaic Virus

BYMV	DAS ELISA	
KS 1040 - 0500R	Reagent Set	500 test
KS 1040 - 1000R	Reagent Set	1000 test
KS 1040 - 5000R	Reagent Set	5000 test
KS 1040 - P	Positive Control	
KS 1040 - N	Negative Control	

Bean Yellow Mosaic Virus - pea isolate

BYMV-pea	DAS ELISA	
KS 1403 - 0500R	Reagent Set	500 test
KS 1403 - 1000R	Reagent Set	1000 test
KS 1403 - 5000R	Reagent Set	5000 test
KS 1403 - P	Positive Control	
KS 1403 - N	Negative Control	

Beet Mild Yellowing Virus - rape isolate

BMVY-rape	DAS ELISA	
KS 1041 - 0500R	Reagent Set	500 test
KS 1041 - 1000R	Reagent Set	1000 test
KS 1041 - 5000R	Reagent Set	5000 test
KS 1041 - P	Positive Control	
KS 1041 - N	Negative Control	

Beet Mild Yellowing Virus

BMVY	TAS ELISA	
KS 1404 - 0500R	Reagent Set	500 test
KS 1404 - 1000R	Reagent Set	1000 test
KS 1404 - 5000R	Reagent Set	5000 test
KS 1404 - P	Positive Control	
KS 1404 - N	Negative Control	

Beet Mosaic Virus

BtMV	DAS ELISA	
KS 1042 - 0500R	Reagent Set	500 test
KS 1042 - 1000R	Reagent Set	1000 test
KS 1042 - 5000R	Reagent Set	5000 test
KS 1042 - P	Positive Control	
KS 1042 - N	Negative Control	

Beet Necrotic Yellow Vein Virus (Rhizomania)

BNYVV	DAS ELISA	
KS 1043 - 0500R	Reagent Set	500 test
KS 1043 - 1000R	Reagent Set	1000 test
KS 1043 - 5000R	Reagent Set	5000 test
KS 1043 - P	Positive Control	
KS 1043 - N	Negative Control	

Beet Necrotic Yellow Vein Virus

BNYVV	TAS ELISA	
KS 1405 - 0500R	Reagent Set	500 test
KS 1405 - 1000R	Reagent Set	1000 test
KS 1405 - 5000R	Reagent Set	5000 test
KS 1405 - P	Positive Control	
KS 1405 - N	Negative Control	

Beet Necrotic Yellow Vein Virus

BNYVV	Compound Direct ELISA	Poly/Mono
CAB 16301/0500	Coating Ab	500 test
CAB 16301/1000	Coating Ab	1000 test
CAB 16301/5000	Coating Ab	5000 test
ECA 16301/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 16301/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 16301/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 16301/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 16301/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 16301/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 16301/0500	Reagent Set,Alk.Phos.	500 test
SRA 16301/1000	Reagent Set,Alk.Phos.	1000 test
SRA 16301/5000	Reagent Set,Alk.Phos.	5000 test
LNC 16301	Negative Control	
LPC 16301	Positive Control	

Beet Soil Born Virus

BSBV	DAS ELISA	
KS 1322 - 0500R	Reagent Set	500 test
KS 1322 - 1000R	Reagent Set	1000 test
KS 1322 - 5000R	Reagent Set	5000 test
KS 1322 - P	Positive Control	
KS 1322 - N	Negative Control	

Beet Western Yellows Virus

BWYV	DAS ELISA	
KS 1045 - 0500R	Reagent Set	500 test
KS 1045 - 1000R	Reagent Set	1000 test
KS 1045 - 5000R	Reagent Set	5000 test
KS 1045 - P	Positive Control	
KS 1045 - N	Negative Control	

Beet Western Yellows Virus

BWYV	Compound Direct ELISA	Poly/Mono
CAB 32100/0500	Coating Ab	500 test
CAB 32100/1000	Coating Ab	1000 test
CAB 32100/5000	Coating Ab	5000 test
ECA 32100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 32100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 32100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 32100/0500	Reagent Set,Alk.Phos.	500 test
SRA 32100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 32100/5000	Reagent Set,Alk.Phos.	5000 test
LNC 32100	Negative Control	
LPC 32100	Positive Control	

Beet Yellows Virus

BYV	DAS ELISA	
KS 1047 - 0500R	Reagent Set	500 test
KS 1047 - 1000R	Reagent Set	1000 test
KS 1047 - 5000R	Reagent Set	5000 test
KS 1047 - P	Positive Control	
KS 1047 - N	Negative Control	

Beet Yellows Virus

BYV	TAS ELISA	
KS 1406 - 0500R	Reagent Set	500 test
KS 1406 - 1000R	Reagent Set	1000 test
KS 1406 - 5000R	Reagent Set	5000 test
KS 1406 - P	Positive Control	
KS 1406 - N	Negative Control	

Bell Pepper Mottle Virus

BPMV	DAS ELISA	
KS 1323 - 0500R	Reagent Set	500 test
KS 1323 - 1000R	Reagent Set	1000 test
KS 1323 - 5000R	Reagent Set	5000 test
KS 1323 - P	Positive Control	
KS 1323 - N	Negative Control	

Belladonna Mottle Virus

BeMV	DAS ELISA	
KS 1324 - 0500R	Reagent Set	500 test
KS 1324 - 1000R	Reagent Set	1000 test
KS 1324 - 5000R	Reagent Set	5000 test
KS 1324 - P	Positive Control	
KS 1324 - N	Negative Control	

Blackeye Cowpea Mosaic Virus

BICMV	DAS ELISA	
KS 1048 - 0500R	Reagent Set	500 test

KS 1048 - 1000R	Reagent Set	1000 test
KS 1048 - 5000R	Reagent Set	5000 test
KS 1048 - P	Positive Control	
KS 1048 - N	Negative Control	

Blackgram Mottle Virus

BMoV	DAS ELISA	
KS 1049 - 0500R	Reagent Set	500 test
KS 1049 - 1000R	Reagent Set	1000 test
KS 1049 - 5000R	Reagent Set	5000 test
KS 1049 - P	Positive Control	
KS 1049 - N	Negative Control	

Blueberry Leaf Mottle Virus

BLMV	Compound Direct ELISA	Poly/Poly
CAB 88001/0500	Coating Ab	500 test
CAB 88001/1000	Coating Ab	1000 test
CAB 88001/5000	Coating Ab	5000 test
ECA 88001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 88001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 88001/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 88001/0500	Reagent Set,Alk.Phos.	500 test
SRA 88001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 88001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 88001	Negative Control	
LPC 88001	Positive Control	

Blueberry Necrotic Shock Virus

BINShV	TAS ELISA	
KS 1052 - 0500R	Reagent Set	500 test
KS 1052 - 1000R	Reagent Set	1000 test
KS 1052 - 5000R	Reagent Set	5000 test
KS 1052 - P	Positive Control	
KS 1052 - N	Negative Control	

Blueberry Scorch Virus

BBScV	DAS ELISA	Poly/Poly
CAB 19100/0500	Coating Ab	500 test
CAB 19100/1000	Coating Ab	1000 test
CAB 19100/5000	Coating Ab	5000 test
ECA 19100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 19100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 19100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 19100/0500	Reagent Set,Alk.Phos.	500 test
SRA 19100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 19100/5000	Reagent Set,Alk.Phos.	5000 test
LPC 19100	Positive Control	
LPC 19100	Negative Control	

Blueberry Shock Virus

BIShV	Compound Direct ELISA	Poly/Mono
CAB 19200/0500	Coating Ab	500 test
CAB 19200/1000	Coating Ab	1000 test
CAB 19200/5000	Coating Ab	5000 test
ECA 19200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 19200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 19200/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 19200/0500	Reagent Set,Alk.Phos.	500 test
SRA 19200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 19200/5000	Reagent Set,Alk.Phos.	5000 test
LPC 19200	Positive Control	

Blueberry Shoestring Virus


BSSV	DAS ELISA	Poly/Poly
CAB 19000/0500	Coating Ab	500 test
CAB 19000/1000	Coating Ab	1000 test
CAB 19000/5000	Coating Ab	5000 test
ECP 19000/0500	Conjugated Ab, Peroxidase	500 test
ECP 19000/1000	Conjugated Ab, Peroxidase	1000 test
ECP 19000/5000	Conjugated Ab, Peroxidase	5000 test
SRP 19000/0500	Reagent Set, Peroxidase	500 test
SRP 19000/1000	Reagent Set, Peroxidase	1000 test
SRP 19000/5000	Reagent Set, Peroxidase	5000 test
LNC 19000	Negative Control	
LPC 19000	Positive Control	

Broad Bean Stain Virus

BBSV	DAS ELISA	
KS 1325 - 0500R	Reagent Set	500 test
KS 1325 - 1000R	Reagent Set	1000 test
KS 1325 - 5000R	Reagent Set	5000 test
KS 1325 - P	Positive Control	
KS 1325 - N	Negative Control	

Broad Bean True Mosaic Virus

BBTMV	DAS ELISA	
KS 1326 - 0500R	Reagent Set	500 test
KS 1326 - 1000R	Reagent Set	1000 test
KS 1326 - 5000R	Reagent Set	5000 test
KS 1326 - P	Positive Control	
KS 1326 - N	Negative Control	



Broad Bean Wilt Virus 1

BBWV1	DAS ELISA	
KS 1055 - 0500R	Reagent Set	500 test
KS 1055 - 1000R	Reagent Set	1000 test
KS 1055 - 5000R	Reagent Set	5000 test
KS 1055 - P	Positive Control	
KS 1055 - N	Negative Control	

Broad Bean Wilt Virus 2

BBWV2	DAS ELISA	
KS 1327 - 0500R	Reagent Set	500 test
KS 1327 - 1000R	Reagent Set	1000 test
KS 1327 - 5000R	Reagent Set	5000 test

Brome Mosaic Virus

BMV	DAS ELISA	Poly/Poly
CAB 3291800 /05500	Coating Ab	500 test
CAB 3291800 /05000	Coating Ab	1000 test
CAB 3291800 /05000	Coating Ab	5000 test
ECA 3291800 /05500	Conjugated Ab,Alk.Phos.	500 test
ECA 3291800 /05000	Conjugated Ab,Alk.Phos.	1000 test
ECA 3291800 /05000	Conjugated Ab,Alk.Phos.	5000 test
SRA 3291800 /05500	Reagent Set,Alk.Phos.	500 test
SRA 3291800 /05000	Reagent Set,Alk.Phos.	1000 test
SRA 3291800 /05000	Reagent Set,Alk.Phos.	5000 test
LNC 29300	Negative Control	
LPC 29300	Positive Control	

Brome Streak Mosaic Virus

BSMV	DAS ELISA	
KS 1057 - 0500R	Reagent Set	500 test
KS 1057 - 1000R	Reagent Set	1000 test
KS 1057 - 5000R	Reagent Set	5000 test
KS 1057 - P	Positive Control	
KS 1057 - N	Negative Control	

Cabbage Blackring Virus (See Turip Mosaic Virus)

is Virus X

DAS ELISA

- 0500R	Reagent Set	500 test
- 1000R	Reagent Set	1000 test
- 5000R	Reagent Set	5000 test
- P	Positive Control	
- N	Negative Control	

achoa Mottle Virus

DAS ELISA

Poly/Poly

00/0500	Coating Ab	500 test
00/1000	Coating Ab	1000 test
00/5000	Coating Ab	5000 test
00/0500	Conjugated Ab,Alk.Phos.	500 test
00/1000	Conjugated Ab,Alk.Phos.	1000 test
00/5000	Conjugated Ab,Alk.Phos.	5000 test
00/0500	Conjugated Ab,Peroxidase	500 test
00/1000	Conjugated Ab,Peroxidase	1000 test
00/5000	Conjugated Ab,Peroxidase	5000 test
00/0096	PathoScreen Kit,Alk.Phos.	96 test
00/0288	PathoScreen Kit,Alk.Phos.	288 test
00/0480	PathoScreen Kit,Alk.Phos.	480 test
00/0096	PathoScreen Kit, Peroxidase	96 test
00/0288	PathoScreen Kit, Peroxidase	288 test
00/0480	PathoScreen Kit, Peroxidase	480 test
00/0500	Reagent Set, Alk.Phos.	500 test
00/1000	Reagent Set, Alk.Phos.	1000 test
00/5000	Reagent Set, Alk.Phos.	5000 test
00/0500	Reagent Set, Peroxidase	96 test
00/1000	Reagent Set, Peroxidase	288 test
00/5000	Reagent Set, Peroxidase	480 test
000/0005	ImmunoStrip Test Kit	5 strips
000/0025	ImmunoStrip Test Kit	25 strips
000/0250	ImmunoStrip Test Kit	250 strips
000/0005	ImmunoStrip Test	5 strips
000/0025	ImmunoStrip Test	25 strips
000/0250	ImmunoStrip Test	250 strips
00	Negative Control	
00	Positive Control	

ImmunoStrip Test kit

50 strips and buffer filled mesh bags

ImmunoStrips Buffer

ImmunoStrips must be used with SEB1 buffer
is offered separately. Please see
pages for a listing of SEB1 buffer products.

Tomato Etched Ring Virus

DAS ELISA

- 0500R	Reagent Set	500 test
- 1000R	Reagent Set	1000 test

KS 1058 - 5000R	Reagent Set	5000 test
KS 1058 - P	Positive Control	
KS 1058 - N	Negative Control	

Carnation Etched Ring Virus

CERV	Direct ELISA	
KS 1407 - 0500R	Reagent Set	500 test
KS 1407 - 1000R	Reagent Set	1000 test
KS 1407 - 5000R	Reagent Set	5000 test
KS 1407 - P	Positive Control	
KS 1407 - N	Negative Control	

Carnation Etched Ring Virus

CERV	Compound Direct ELISA	Poly/Mono
CAB 34100/0500	Coating Ab	500 test
CAB 34100/1000	Coating Ab	1000 test
CAB 34100/5000	Coating Ab	5000 test
ECA 34100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 34100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 34100/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 34100/0500	Conjugated Ab,Peroxidase	500 test
ECP 34100/1000	Conjugated Ab,Peroxidase	1000 test
ECP 34100/5000	Conjugated Ab,Peroxidase	5000 test
PSA 34100/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 34100/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 34100/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 34100/0096	PathoScreen Kit,Peroxidase	96 test
PSP 34100/0288	PathoScreen Kit,Peroxidase	288 test
PSP 34100/0480	PathoScreen Kit,Peroxidase	480 test
SRA 34100/0500	Reagent Set,Alk.Phos.	500 test
SRA 34100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 34100/5000	Reagent Set,Alk.Phos.	5000 test
SRP 34100/0500	Reagent Set,Peroxidase	500 test
SRP 34100/1000	Reagent Set,Peroxidase	1000 test
SRP 34100/5000	Reagent Set,Peroxidase	5000 test
LNC 34100	Negative Control	
LPC 34100	Positive Control	

Carnation Italian Ringspot Virus

CIRSV	DAS ELISA	
KS 1383 - 0500R	Reagent Set	500 test
KS 1383 - 1000R	Reagent Set	1000 test
KS 1383 - 5000R	Reagent Set	5000 test
KS 1383 - P	Positive Control	
KS 1383 - N	Negative Control	

Carnation Latent Virus

CLV	DAS ELISA	
KS 1060 - 0500R	Reagent Set	500 test

KS 1060 - 1000R	Reagent Set	1000 test
KS 1060 - 5000R	Reagent Set	5000 test
KS 1060 - P	Positive Control	
KS 1060 - N	Negative Control	

Carnation Latent Virus

CLV	Compound Direct/DAS ELISA	Poly/Mono
CAB 34000/0500	Coating Ab	500 test
CAB 34000/1000	Coating Ab	1000 test
CAB 34000/5000	Coating Ab	5000 test
ECA 34000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 34000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 34000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 34000/0500	Conjugated Ab,Peroxidase	500 test
ECP 34000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 34000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 34000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 34000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 34000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 34000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 34000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 34000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 34000/0500	Reagent Set,Alk.Phos.	500 test
SRA 34000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 34000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 34000/0500	Reagent Set,Peroxidase	500 test
SRP 34000/1000	Reagent Set,Peroxidase	1000 test
SRP 34000/5000	Reagent Set,Peroxidase	5000 test
LNC 34000	Negative Control	
LPC 34000	Positive Control	

Carnation Mottle Virus

CarMV	DAS ELISA	Poly/Poly
CAB 68000/0500	Coating Ab	500 test
CAB 68000/1000	Coating Ab	1000 test
CAB 68000/5000	Coating Ab	5000 test
ECA 68000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 68000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 68000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 68000/0500	Conjugated Ab,Peroxidase	500 test
ECP 68000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 68000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 68000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 68000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 68000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 68000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 68000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 68000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 68000/0500	Reagent Set,Alk.Phos.	500 test
SRA 68000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 68000/5000	Reagent Set,Alk.Phos.	5000 test

식물병 Viruses

SRP 68000/0500	Reagent Set,Peroxidase	500 test
SRP 68000/1000	Reagent Set,Peroxidase	1000 test
SRP 68000/5000	Reagent Set,Peroxidase	5000 test
LNC 68000	Negative Control	
LPC 68000	Positive Control	

Carnation Necrotic Fleck Virus

CNFV	Indirect/DAS ELISA	Mono/Mono
CAB 68500/0500	Coating Ab	500 test
CAB 68500/1000	Coating Ab	1000 test
CAB 68500/5000	Coating Ab	5000 test
ECA 68500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 68500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 68500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 68500/0500	Conjugated Ab,Peroxidase	500 test
ECP 68500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 68500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 68500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 68500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 68500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 68500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 68500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 68500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 68500/0500	Reagent Set,Alk.Phos.	500 test
SRA 68500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 68500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 68500/0500	Reagent Set,Peroxidase	500 test
SRP 68500/1000	Reagent Set,Peroxidase	1000 test
SRP 68500/5000	Reagent Set,Peroxidase	5000 test
LNC 68500	Negative Control	
LPC 68500	Positive Control	

Carnation Ringspot Virus

CRSV	DAS ELISA	Poly/Poly
CAB 68300/0500	Coating Ab	500 test
CAB 68300/1000	Coating Ab	1000 test
CAB 68300/5000	Coating Ab	5000 test
ECA 68300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 68300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 68300/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 68300/0500	Conjugated Ab,Peroxidase	500 test
ECP 68300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 68300/5000	Conjugated Ab,Peroxidase	5000 test
PSA 68300/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 68300/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 68300/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 68300/0096	PathoScreen Kit,Peroxidase	96 test
PSP 68300/0288	PathoScreen Kit,Peroxidase	288 test
PSP 68300/0480	PathoScreen Kit,Peroxidase	480 test

SRA 68300/0500	Reagent Set,Alk.Phos.	500 test
SRA 68300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 68300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 68300/0500	Reagent Set,Peroxidase	500 test
SRP 68300/1000	Reagent Set,Peroxidase	1000 test
SRP 68300/5000	Reagent Set,Peroxidase	5000 test
LNC 68300	Negative Control	
LPC 68300	Positive Control	

Carnation Vein Mottle Virus

CVMV	DAS ELISA	
KS 1066 - 0500R	Reagent Set	500 test
KS 1066 - 1000R	Reagent Set	1000 test
KS 1066 - 5000R	Reagent Set	5000 test
KS 1066 - P	Positive Control	
KS 1066 - N	Negative Control	

Cauliflower Mosaic Virus

CaMV	DAS ELISA	Poly/Poly
CAB 45600/0500	Coating Ab	500 test
CAB 45600/1000	Coating Ab	1000 test
CAB 45600/5000	Coating Ab	5000 test
ECP 45600/0500	Conjugated Ab, Peroxidase	500 test
ECP 45600/1000	Conjugated Ab, Peroxidase	1000 test
ECP 45600/5000	Conjugated Ab, Peroxidase	5000 test
SRP 45600/0500	Reagent Set, Peroxidase	500 test
SRP 45600/1000	Reagent Set, Peroxidase	1000 test
SRP 45600/5000	Reagent Set, Peroxidase	5000 test
LPC 45600	Positive Control	

Celery Mosaic Virus

CeMV	DAS ELISA	
KS 1328 - 0500R	Reagent Set	500 test
KS 1328 - 1000R	Reagent Set	1000 test
KS 1328 - 5000R	Reagent Set	5000 test
KS 1328 - P	Positive Control	
KS 1328 - N	Negative Control	

Cereal Yellow Dwarf Virus - rpv

CYDV-rpv	Compound Direct ELISA	Poly/Mono
CAB 28001/0500	Coating Ab	500 test
CAB 28001/1000	Coating Ab	1000 test
CAB 28001/5000	Coating Ab	5000 test
ECA 28001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 28001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 28001/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 28001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 28001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 28001/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 28001/0500	Reagent Set,Alk.Phos.	500 test

SRA 28001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 28001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 28001	Negative Control	
LPC 28001	Positive Control	

Cherry Leaf Roll Virus

CLRV	DAS ELISA	
KS 1068 - 0500R	Reagent Set	500 test
KS 1068 - 1000R	Reagent Set	1000 test
KS 1068 - 5000R	Reagent Set	5000 test
KS 1068 - P	Positive Control	
KS 1068 - N	Negative Control	

Cherry Leaf Roll Virus - birch isolate

CLRV-birch	DAS ELISA	
KS 1069 - 0500R	Reagent Set	500 test
KS 1069 - 1000R	Reagent Set	1000 test
KS 1069 - 5000R	Reagent Set	5000 test
KS 1069 - P	Positive Control	
KS 1069 - N	Negative Control	

Cherry Leaf Roll Virus - cherry strain

CLRV-ch	DAS ELISA	
KS 1070 - 0500R	Reagent Set	500 test
KS 1070 - 1000R	Reagent Set	1000 test
KS 1070 - 5000R	Reagent Set	5000 test
KS 1070 - P	Positive Control	
KS 1070 - N	Negative Control	

Cherry Leaf Roll Virus - elderberry isolate

CLRV-elderberry	DAS ELISA	
KS 1071 - 0500R	Reagent Set	500 test
KS 1071 - 1000R	Reagent Set	1000 test
KS 1071 - 5000R	Reagent Set	5000 test
KS 1071 - P	Positive Control	
KS 1071 - N	Negative Control	

Chilli Veinal Mottle Virus

ChiVMV	DAS ELISA	
KS 1329 - 0500R	Reagent Set	500 test
KS 1329 - 1000R	Reagent Set	1000 test
KS 1329 - 5000R	Reagent Set	5000 test
KS 1329 - P	Positive Control	
KS 1329 - N	Negative Control	

**Chrysanthemum Aspermy Virus
(See Tomato Aspermy Virus)****Chrysanthemum Stem Necrosis Virus**

CSNV	DAS ELISA	
KS 1330 - 0500R	Reagent Set	500 test
KS 1330 - 1000R	Reagent Set	1000 test
KS 1330 - 5000R	Reagent Set	5000 test
KS 1330 - P	Positive Control	
KS 1330 - N	Negative Control	

Chrysanthemum Virus B

CVB	DAS ELISA	Poly/Poly
CAB 81000/0500	Coating Ab	500 test
CAB 81000/1000	Coating Ab	1000 test
CAB 81000/5000	Coating Ab	5000 test
ECA 81000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 81000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 81000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 81000/0500	Conjugated Ab,Peroxidase	500 test
ECP 81000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 81000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 81000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 81000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 81000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 81000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 81000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 81000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 81000/0500	Reagent Set,Alk.Phos.	500 test
SRA 81000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 81000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 81000/0500	Reagent Set,Peroxidase	500 test
SRP 81000/1000	Reagent Set,Peroxidase	1000 test
SRP 81000/5000	Reagent Set,Peroxidase	5000 test
LNC 81000	Negative Control	
LPC 81000	Positive Control	

**Citrus Tatter Leaf Virus
(See Apple Stem Grooving Virus)****Citrus Tristeza Virus**

CTV	DAS ELISA	
KS 1075 - 0500R	Reagent Set	500 test
KS 1075 - 1000R	Reagent Set	1000 test
KS 1075 - 5000R	Reagent Set	5000 test
KS 1075 - P	Positive Control	
KS 1075 - N	Negative Control	

Citrus Tristeza Virus

CTV	Compound Direct ELISA	Poly/Mono
CAB 78900/0500	Coating Ab	500 test
CAB 78900/1000	Coating Ab	1000 test
CAB 78900/5000	Coating Ab	5000 test

식물병 Viruses

ECA 78900/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 78900/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 78900/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 78900/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 78900/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 78900/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 78900/0500	Reagent Set,Alk.Phos.	500 test
SRA 78900/1000	Reagent Set,Alk.Phos.	1000 test
SRA 78900/5000	Reagent Set,Alk.Phos.	5000 test
LNC 78900	Negative Control	
LPC 78900	Positive Control	

Clover Yellow Mosaic Virus

CYMV	DAS ELISA	
KS 1331 - 0500R	Reagent Set	500 test
KS 1331 - 1000R	Reagent Set	1000 test
KS 1331 - 5000R	Reagent Set	5000 test
KS 1331 - P	Positive Control	
KS 1331 - N	Negative Control	

Clover Yellow Vein Virus

CYVV	DAS ELISA	
KS 1077 - 0500R	Reagent Set	500 test
KS 1077 - 1000R	Reagent Set	1000 test
KS 1077 - 5000R	Reagent Set	5000 test
KS 1077 - P	Positive Control	
KS 1077 - N	Negative Control	

Cocksfoot Mild Mosaic Virus

CMMV	DAS ELISA	
KS 1332 - 0500R	Reagent Set	500 test
KS 1332 - 1000R	Reagent Set	1000 test
KS 1332 - 5000R	Reagent Set	5000 test
KS 1332 - P	Positive Control	
KS 1332 - N	Negative Control	

Cowpea Aphid-Borne Mosaic Virus

CABMV	DAS ELISA	
KS 1333 - 0500R	Reagent Set	500 test
KS 1333 - 1000R	Reagent Set	1000 test
KS 1333 - 5000R	Reagent Set	5000 test
KS 1333 - P	Positive Control	
KS 1333 - N	Negative Control	

Cowpea Chlorotic Mottle Virus

CCMV	DAS ELISA	
KS 1078 - 0500R	Reagent Set	500 test
KS 1078 - 1000R	Reagent Set	1000 test
KS 1078 - 5000R	Reagent Set	5000 test
KS 1078 - P	Positive Control	

KS 1078 - N	Negative Control
-------------	------------------

Cowpea Mild Mottle Virus

CPMMV	DAS ELISA	
KS 1334 - 0500R	Reagent Set	500 test
KS 1334 - 1000R	Reagent Set	1000 test
KS 1334 - 5000R	Reagent Set	5000 test
KS 1334 - P	Positive Control	
KS 1334 - N	Negative Control	

Cowpea Mottle Virus

CPMoV	DAS ELISA	
KS 1335 - 0500R	Reagent Set	500 test
KS 1335 - 1000R	Reagent Set	1000 test
KS 1335 - 5000R	Reagent Set	5000 test
KS 1335 - P	Positive Control	
KS 1335 - N	Negative Control	

Cowpea Yellow Mosaic Virus

CPYMV	DAS ELISA	
KS 1336 - 0500R	Reagent Set	500 test
KS 1336 - 1000R	Reagent Set	1000 test
KS 1336 - 5000R	Reagent Set	5000 test
KS 1336 - P	Positive Control	
KS 1336 - N	Negative Control	

Cowpea Mosaic Virus

CPMV	DAS ELISA	Poly/Poly
CAB 24200/0500	Coating Ab	500 test
CAB 24200/1000	Coating Ab	1000 test
CAB 24200/5000	Coating Ab	5000 test
ECA 24200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 24200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 24200/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 24200/0500	Reagent Set,Alk.Phos.	500 test
SRA 24200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 24200/5000	Reagent Set,Alk.Phos.	5000 test
LNC 24200	Negative Control	
LPC 24200	Positive Control	

Cowpea Severe Mosaic Virus

CPSMV	DAS ELISA	
KS 1080 - 0500R	Reagent Set	500 test
KS 1080 - 1000R	Reagent Set	1000 test
KS 1080 - 5000R	Reagent Set	5000 test
KS 1080 - P	Positive Control	
KS 1080 - N	Negative Control	

Cucumber Green Mottle Mosaic Virus

CGMMV	DAS ELISA	Poly/Poly
CAB 45700/0500	Coating Ab	500 test
CAB 45700/1000	Coating Ab	1000 test
CAB 45700/5000	Coating Ab	5000 test
ECA 45700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 45700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 45700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 45700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 45700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 45700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 45700/0500	Reagent Set,Alk.Phos.	500 test
SRA 45700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 45700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 45700	Negative Control	
LPC 45700	Positive Control	

Cucumber Green Mottle Mosaic Virus - C
(See Kyuri Green Mottle Mosaic Virus)

Cucumber Green Mottle Mosaic Virus - Wa
(See Cucumber Green Mottle Mosaic Virus)

Cucumber Mosaic Virus

CMV	DAS ELISA	
KS 1337 - 0500R	Reagent Set	500 test
KS 1337 - 1000R	Reagent Set	1000 test
KS 1337 - 5000R	Reagent Set	5000 test
KS 1337 - P	Positive Control	
KS 1337 - N	Negative Control	

Cucumber Mosaic Virus

CMV	Indirect ELISA	
KS 1434 - 0500R	Reagent Set	500 test
KS 1434 - 1000R	Reagent Set	1000 test
KS 1434 - 5000R	Reagent Set	5000 test
KS 1434 - P	Positive Control	
KS 1434 - N	Negative Control	

Cucumber Mosaic Virus

CMV	Compound Direct ELISA	Poly/Mono
CAB 44501/0500	Coating Ab	500 test
CAB 44501/1000	Coating Ab	1000 test
CAB 44501/5000	Coating Ab	5000 test
ECA 44501/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44501/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44501/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 44501/0500	Conjugated Ab,Peroxidase	500 test
ECP 44501/1000	Conjugated Ab,Peroxidase	1000 test
ECP+44501/5000	Conjugated Ab,Peroxidase	5000 test

PSA 44501/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 44501/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 44501/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 44501/0096	PathoScreen Kit,Peroxidase	96 test
PSP 44501/0288	PathoScreen Kit,Peroxidase	288 test
PSP 44501/0480	PathoScreen Kit,Peroxidase	480 test
SRA 44501/0500	Reagent Set,Alk.Phos.	500 test
SRA 44501/1000	Reagent Set,Alk.Phos.	1000 test
SRA 44501/5000	Reagent Set,Alk.Phos.	5000 test
SRP 44501/0500	Reagent Set,Peroxidase	500 test
SRP 44501/1000	Reagent Set,Peroxidase	1000 test
SRP 44501/5000	Reagent Set,Peroxidase	5000 test
*ISK 44501/0005	ImmunoStrip Test Kit	5 strips
*ISK 44501/0025	ImmunoStrip Test Kit	25 strips
*ISK 44501/0250	ImmunoStrip Test Kit	250 strips
*STX 44501/0005	ImmunoStrip Test	5 strips
*STX 44501/0025	ImmunoStrip Test	25 strips
*STX 44501/0250	ImmunoStrip Test	250 strips
LNC 44501	Negative Control	
LPC 44501	Positive Control	

*** CMV ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** CMV ImmunoStrip Buffer**

CMV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

CMV, INSV, TSWV, TMV ImmunoStrip Comb

*ISK 37800/0024	ImmunoStrip Comb Kit, 24 combs and 25 buffer filled mesh bags	
*STX 37800/0024	ImmunoStrip Comb	24 combs

* ImmunoStrip Comb test system for CMV, INSV, TSWV, TMV are packaged in sets of 24.

ImmunoStrip combs must be used with SEB1 buffer.

Please see accessories for a listing of SEB1 buffer products. ImmunoStrip kits include sample extract bags.

Cucumber Mosaic Virus - Subgroup I

CMV I	Compound Direct ELISA	Poly/Mono
CAB 44700/0500	Coating Ab	500 test
CAB 44700/1000	Coating Ab	1000 test
CAB 44700/5000	Coating Ab	5000 test
ECA 44700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 44700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 44700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 44700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 44700/0500	Reagent Set,Alk.Phos.	500 test
SRA 44700/1000	Reagent Set,Alk.Phos.	1000 test

식물병 Viruses

SRA 44700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 44700	Negative Control	
LPC 44700	Positive Control	

Cucumber Mosaic Virus - Subgroup II

CMV II	Compound Direct ELISA	Poly
CAB 44800/0500	Coating Ab	500 test
CAB 44800/1000	Coating Ab	1000 test
CAB 44800/5000	Coating Ab	5000 test
ECA 44800/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44800/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44800/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 44800/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 44800/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 44800/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 44800/0500	Reagent Set,Alk.Phos.	500 test
SRA 44800/1000	Reagent Set,Alk.Phos.	1000 test
SRA 44800/5000	Reagent Set,Alk.Phos.	5000 test
LPC 44800	Positive Control	
LNC 44800	Negative Control	

Cucumber Mosaic Virus - DTL

(See Cucumber Mosaic Virus - Subgroup I)

Cucumber Mosaic Virus - ToRS

(See Cucumber Mosaic Virus - Subgroup II)

Cucumber Mosaic Virus - DTL/ToRS

CMV-DTL/ToRS	DAS ELISA	
KS 1089 - 0500R	Reagent Set	500 test
KS 1089 - 1000R	Reagent Set	1000 test
KS 1089 - 5000R	Reagent Set	5000 test
KS 1089 - P	Positive Control	
KS 1089 - N	Negative Control	

Cucumber Necrosis Virus

CNV	DAS ELISA	
KS 1338 - 0500R	Reagent Set	500 test
KS 1338 - 1000R	Reagent Set	1000 test
KS 1338 - 5000R	Reagent Set	5000 test
KS 1338 - P	Positive Control	
KS 1338 - N	Negative Control	

Cymbidium Mosaic Virus

CyMV	DAS ELISA	
KS 1090 - 0500R	Reagent Set	500 test
KS 1090 - 1000R	Reagent Set	1000 test
KS 1090 - 5000R	Reagent Set	5000 test
KS 1090 - P	Positive Control	
KS 1090 - N	Negative Control	

Cymbidium Mosaic Virus

CyMV	Compound Direct ELISA	Poly/Mono
CAB 89500/0500	Coating Ab	500 test
CAB 89500/1000	Coating Ab	1000 test
CAB 89500/5000	Coating Ab	5000 test
ECA 89500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 89500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 89500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 89500/0500	Conjugated Ab,Peroxidase	500 test
ECP 89500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 89500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 89500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 89500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 89500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 89500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 89500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 89500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 89500/0500	Reagent Set,Alk.Phos.	500 test
SRA 89500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 89500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 89500/0500	Reagent Set,Peroxidase	500 test
SRP 89500/1000	Reagent Set,Peroxidase	1000 test
SRP 89500/5000	Reagent Set,Peroxidase	5000 test
LNC 89500	Negative Control	
LPC 89500	Positive Control	

Cymbidium Ringspot Virus

CyRSV	DAS ELISA	Poly/Poly
CAB 33001/0500	Coating Ab	500 test
CAB 33001/1000	Coating Ab	1000 test
CAB 33001/5000	Coating Ab	5000 test
ECA 33001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 33001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 33001/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 33001/0500	Conjugated Ab,Peroxidase	500 test
ECP 33001/1000	Conjugated Ab,Peroxidase	1000 test
ECP 33001/5000	Conjugated Ab,Peroxidase	5000 test
PSA 33001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 33001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 33001/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 33001/0096	PathoScreen Kit,Peroxidase	96 test
PSP 33001/0288	PathoScreen Kit,Peroxidase	288 test
PSP 33001/0480	PathoScreen Kit,Peroxidase	480 test
SRA 33001/0500	Reagent Set,Alk.Phos.	500 test
SRA 33001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 33001/5000	Reagent Set,Alk.Phos.	5000 test
SRP 33001/0500	Reagent Set,Peroxidase	500 test
SRP 33001/1000	Reagent Set,Peroxidase	1000 test
SRP 33001/5000	Reagent Set,Peroxidase	5000 test
LNC 33001	Negative Control	

LPC 33001	Positive Control
-----------	------------------

Dahlia Mosaic Virus

DMV	DAS ELISA	
KS 1093 - 0500R	Reagent Set	500 test
KS 1093 - 1000R	Reagent Set	1000 test
KS 1093 - 5000R	Reagent Set	5000 test
KS 1093 - P	Positive Control	
KS 1093 - N	Negative Control	

Dasheen Mosaic Virus

DsMV	DAS ELISA	Poly/Poly
CAB 48000/0500	Coating Ab	500 test
CAB 48000/1000	Coating Ab	1000 test
CAB 48000/5000	Coating Ab	5000 test
ECA 48000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 48000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 48000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 48000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 48000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 48000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 48000/0500	Reagent Set,Alk.Phos.	500 test
SRA 48000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 48000/5000	Reagent Set,Alk.Phos.	5000 test
LPC 48000	Positive Control	

Eggplant Mosaic Virus

EMV	DAS ELISA	
KS 1095 - 0500R	Reagent Set	500 test
KS 1095 - 1000R	Reagent Set	1000 test
KS 1095 - 5000R	Reagent Set	5000 test
KS 1095 - P	Positive Control	
KS 1095 - N	Negative Control	

Eggplant Mottled Dwarf Virus

EMDV	DAS ELISA	
KS 1339 - 0500R	Reagent Set	500 test
KS 1339 - 1000R	Reagent Set	1000 test
KS 1339 - 5000R	Reagent Set	5000 test
KS 1339 - P	Positive Control	
KS 1339 - N	Negative Control	

Elderberry Latent Virus

ELV	DAS ELISA	
KS 1340 - 0500R	Reagent Set	500 test
KS 1340 - 1000R	Reagent Set	1000 test
KS 1340 - 5000R	Reagent Set	5000 test
KS 1340 - P	Positive Control	
KS 1340 - N	Negative Control	

Flower - Break Virus

FBV	DAS ELISA	
KS 1096 - 0500R	Reagent Set	500 test
KS 1096 - 1000R	Reagent Set	1000 test
KS 1096 - 5000R	Reagent Set	5000 test
KS 1096 - P	Positive Control	
KS 1096 - N	Negative Control	

Freesia Mosaic Virus

FreMV	DAS ELISA	
KS 1097 - 0500R	Reagent Set	500 test
KS 1097 - 1000R	Reagent Set	1000 test
KS 1097 - 5000R	Reagent Set	5000 test
KS 1097 - P	Positive Control	
KS 1097 - N	Negative Control	

Garlic Common Latent Virus

GCLV	DAS ELISA	Poly/Poly
CAB 20400/0500	Coating Ab	500 test
CAB 20400/1000	Coating Ab	1000 test
CAB 20400/5000	Coating Ab	5000 test
ECA 20400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20400/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 20400/0500	Reagent Set,Alk.Phos.	500 test
SRA 20400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20400/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20400	Negative Control	

Garlic Virus D

GarV-D	DAS ELISA	
KS 1384 - 0500R	Reagent Set	500 test
KS 1384 - 1000R	Reagent Set	1000 test
KS 1384 - 5000R	Reagent Set	5000 test
KS 1384 - P	Positive Control	
KS 1384 - N	Negative Control	

Gemini Virus Screen**(Indian, European, African strains of TYLCV)**

TYLCV	TAS ELISA	
KS 1396 - 0100R	Reagent Set	100 tests
KS 1396 - P	Positive Control	
KS 1396 - N	Negative Control	

Grapevine Algerian Latent Virus

GALV	DAS ELISA	
KS 1101 - 0500R	Reagent Set	500 test
KS 1101 - 1000R	Reagent Set	1000 test
KS 1101 - 5000R	Reagent Set	5000 test
KS 1101 - P	Positive Control	

KS 1101 - N	Negative Control
-------------	------------------

Grapevine Court-Noue Disease

GFLV+ArMV	DAS ELISA	
KS 1408 - 0500R	Reagent Set	500 test
KS 1408 - 1000R	Reagent Set	1000 test
KS 1408 - 5000R	Reagent Set	5000 test
KS 1408 - 0500S	Full Kit	500 test
KS 1408 - 1000S	Full Kit	1000 test
KS 1408 - 5000S	Full Kit	5000 test
KS 1408 - P	Positive Control	
KS 1408 - N	Negative Control	

Grapevine Fanleaf Virus

GFLV	DAS ELISA	Poly/Poly
CAB 27701/0500	Coating Ab	500 test
CAB 27701/1000	Coating Ab	1000 test
CAB 27701/5000	Coating Ab	5000 test
ECA 27701/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 27701/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 27701/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 27701/0500	Reagent Set,Alk.Phos.	500 test
SRA 27701/1000	Reagent Set,Alk.Phos.	1000 test
SRA 27701/5000	Reagent Set,Alk.Phos.	5000 test
LNC 27701	Negative Control	
LPC 27701	Positive Control	

Grapevine Fanleaf Virus in grape

GFLV	DAS ELISA	
XAB 27702/0500	Coating Ab	500 test
XAB 27702/1000	Coating Ab	1000 test
XAB 27702/5000	Coating Ab	5000 test
XCA 27702/0500	Conjugated Ab,Alk.Phos.	500 test
XCA 27702/1000	Conjugated Ab,Alk.Phos.	1000 test
XCA 27702/5000	Conjugated Ab,Alk.Phos.	5000 test
XRA 27702/0500	Reagent Set,Alk.Phos.	500 test
XRA 27702/1000	Reagent Set,Alk.Phos.	1000 test
XRA 27702/5000	Reagent Set,Alk.Phos.	5000 test
LNC 27702	Negative Control	
LPC 27702	Positive Control	

Grapevine Fleck Virus

GFkV	DAS ELISA	
KS 1103 - 0500R	Reagent Set	500 test
KS 1103 - 1000R	Reagent Set	1000 test
KS 1103 - 5000R	Reagent Set	5000 test
KS 1103 - P	Positive Control	
KS 1103 - N	Negative Control	

Grapevine Fleck Virus

GFkV	TAS ELISA	
KS 1409 - 0500R	Reagent Set	500 test
KS 1409 - 1000R	Reagent Set	1000 test
KS 1409 - 5000R	Reagent Set	5000 test
KS 1409 - 0500S	Full Kit	500 test
KS 1409 - 1000S	Full Kit	1000 test
KS 1409 - 5000S	Full Kit	5000 test
KS 1409 - P	Positive Control	
KS 1409 - N	Negative Control	

Grapevine Leaf Roll associated Virus 1

GLRaV 1	DAS ELISA	
KS 1104 - 0500R	Reagent Set	500 test
KS 1104 - 1000R	Reagent Set	1000 test
KS 1104 - 5000R	Reagent Set	5000 test
KS 1104 - 0500S	Full Kit	500 test
KS 1104 - 1000S	Full Kit	1000 test
KS 1104 - 5000S	Full Kit	5000 test
KS 1104 - P	Positive Control	
KS 1104 - N	Negative Control	

Grapevine Leaf Roll associated Virus (1, 3)

GLRaV 1,3	DAS ELISA	
KS 1107 - 0500R	Reagent Set	500 test
KS 1107 - 1000R	Reagent Set	1000 test
KS 1107 - 5000R	Reagent Set	5000 test
KS 1107 - 0500S	Full Kit	500 test
KS 1107 - 1000S	Full Kit	1000 test
KS 1107 - 5000S	Full Kit	5000 test
KS 1107 - P	Positive Control	
KS 1107 - N	Negative Control	

Grapevine Leaf Roll associated Virus 2

GLRaV 2	DAS ELISA	
KS 1105 - 0500R	Reagent Set	500 test
KS 1105 - 1000R	Reagent Set	1000 test
KS 1105 - 5000R	Reagent Set	5000 test
KS 1105 - P	Positive Control	
KS 1105 - N	Negative Control	

Grapevine Leaf Roll associated Virus 2

GLRaV 2	TAS ELISA	
KS 1435 - 0500R	Reagent Set	500 test
KS 1435 - 1000R	Reagent Set	1000 test
KS 1435 - 5000R	Reagent Set	5000 test
KS 1435 - 0500S	Full Kit	500 test
KS 1435 - 1000S	Full Kit	1000 test
KS 1435 - 5000S	Full Kit	5000 test
KS 1435 - P	Positive Control	

KS 1435 - N	Negative Control
-------------	------------------

Grapevine Leaf Roll associated Virus 3

GLRaV 3	DAS ELISA	
KS 1106 - 0500R	Reagent Set	500 test
KS 1106 - 1000R	Reagent Set	1000 test
KS 1106 - 5000R	Reagent Set	5000 test
KS 1106 - 0500S	Full Kit	500 test
KS 1106 - 1000S	Full Kit	1000 test
KS 1106 - 5000S	Full Kit	5000 test
KS 1106 - P	Positive Control	
KS 1106 - N	Negative Control	

Grapevine Leaf Roll associated Virus 5

GLRaV 5	DAS ELISA	
KS 1108 - 0500R	Reagent Set	500 test
KS 1108 - 1000R	Reagent Set	1000 test
KS 1108 - 5000R	Reagent Set	5000 test
KS 1108 - 0500S	Full Kit	500 test
KS 1108 - 1000S	Full Kit	1000 test
KS 1108 - 5000S	Full Kit	5000 test
KS 1108 - P	Positive Control	
KS 1108 - N	Negative Control	

Grapevine Leaf Roll associated Virus 6

GLRaV 6	DAS ELISA	
KS 1109 - 0500R	Reagent Set	500 test
KS 1109 - 1000R	Reagent Set	1000 test
KS 1109 - 5000R	Reagent Set	5000 test
KS 1109 - P	Positive Control	
KS 1109 - N	Negative Control	

Grapevine Leaf Roll associated Virus 7

GLRaV 7	DAS ELISA	
KS 1110 - 0500R	Reagent Set	500 test
KS 1110 - 1000R	Reagent Set	1000 test
KS 1110 - 5000R	Reagent Set	5000 test
KS 1110 - 0500S	Full Kit	500 test
KS 1110 - 1000S	Full Kit	1000 test
KS 1110 - 5000S	Full Kit	5000 test
KS 1110 - P	Positive Control	
KS 1110 - N	Negative Control	

Grapevine Virus A

GVA	DAS ELISA	
KS 1111 - 0500R	Reagent Set	500 test
KS 1111 - 1000R	Reagent Set	1000 test
KS 1111 - 5000R	Reagent Set	5000 test
KS 1111 - 0500S	Full Kit	500 test
KS 1111 - 1000S	Full Kit	1000 test

KS 1111 - 5000S	Full Kit	5000 test
KS 1111 - P	Positive Control	
KS 1111 - N	Negative Control	

Grapevine Vitivirus B

GVB	DAS ELISA	
KS 1411 - 0500R	Reagent Set	500 test
KS 1411 - 1000R	Reagent Set	1000 test
KS 1411 - 5000R	Reagent Set	5000 test
KS 1411 - 0500S	Full Kit	500 test
KS 1411 - 1000S	Full Kit	1000 test
KS 1411 - 5000S	Full Kit	5000 test
KS 1411 - P	Positive Control	
KS 1411 - N	Negative Control	

Groundnut Bud Necrosis Virus (See Watermelon Silver Mottle Virus)

Groundnut Ringspot Virus

GRSV	DAS ELISA	
KS 1113 - 0500R	Reagent Set	500 test
KS 1113 - 1000R	Reagent Set	1000 test
KS 1113 - 5000R	Reagent Set	5000 test
KS 1113 - P	Positive Control	
KS 1113 - N	Negative Control	

Groundnut Ringspot Virus and Tomato Chlorotic Spot Virus

GRSV+TCSV	Compound Direct ELISA	Mixture
CAB 61000/0500	Coating Ab	500 test
CAB 61000/1000	Coating Ab	1000 test
CAB 61000/5000	Coating Ab	5000 test
ECA 61000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 61000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 61000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 61000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 61000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 61000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 61000/0500	Reagent Set,Alk.Phos.	500 test
SRA 61000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 61000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 61000	Negative Control	

Helenium Virus S

HVS	DAS ELISA	
KS 1342 - 0500R	Reagent Set	500 test
KS 1342 - 1000R	Reagent Set	1000 test
KS 1342 - 5000R	Reagent Set	5000 test
KS 1342 - P	Positive Control	
KS 1342 - N	Negative Control	

식물병 Viruses

Hibiscus Chlorotic Ringspot Virus

HCRSV	DAS ELISA	Poly/Poly
CAB 12300/0500	Coating Ab	500 test
CAB 12300/1000	Coating Ab	1000 test
CAB 12300/5000	Coating Ab	5000 test
ECA 12300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 12300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 12300/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 12300/0500	Conjugated Ab,Peroxidase	500 test
ECP 12300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 12300/5000	Conjugated Ab,Peroxidase	5000 test
SRA 12300/0500	Reagent Set,Alk.Phos.	500 test
SRA 12300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 12300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 12300/0500	Reagent Set,Peroxidase	500 test
SRP 12300/1000	Reagent Set,Peroxidase	1000 test
SRP 12300/5000	Reagent Set,Peroxidase	5000 test
LPC 12300	Positive Control	

Hop American Latent Virus

HALV	DAS ELISA	
KS 1116 - 0500R	Reagent Set	500 test
KS 1116 - 1000R	Reagent Set	1000 test
KS 1116 - 5000R	Reagent Set	5000 test
KS 1116 - P	Positive Control	
KS 1116 - N	Negative Control	

Hop Latent Virus

HpLV	DAS ELISA	
KS 1412 - 0500R	Reagent Set	500 test
KS 1412 - 1000R	Reagent Set	1000 test
KS 1412 - 5000R	Reagent Set	5000 test
KS 1412 - P	Positive Control	
KS 1412 - N	Negative Control	

Hop Mosaic Virus

HpMV	DAS ELISA	
KS 1117 - 0500R	Reagent Set	500 test
KS 1117 - 1000R	Reagent Set	1000 test
KS 1117 - 5000R	Reagent Set	5000 test
KS 1117 - P	Positive Control	
KS 1117 - N	Negative Control	

Hosta Virus X

HsVX	DAS ELISA	Poly/Poly
CAB 16600/0500	Coating Ab	500 test
CAB 16600/1000	Coating Ab	1000 test
CAB 16600/5000	Coating Ab	5000 test
ECA 16600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 16600/1000	Conjugated Ab,Alk.Phos.	1000 test

ECA 16600/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 16600/0500	Reagent Set,Alk.Phos.	500 test
SRA 16600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 16600/5000	Reagent Set,Alk.Phos.	5000 test
LNC 16600	Negative Control	
LPC 16600	Positive Control	

Hydrangea Ringspot Virus

HRSV	DAS ELISA	
KS 1119 - 0500R	Reagent Set	500 test
KS 1119 - 1000R	Reagent Set	1000 test
KS 1119 - 5000R	Reagent Set	5000 test
KS 1119 - P	Positive Control	
KS 1119 - N	Negative Control	

Impatiens Necrotic Spot Virus

INSV	DAS ELISA	Poly/Poly
CAB 20500/0500	Coating Ab	500 test
CAB 20500/1000	Coating Ab	1000 test
CAB 20500/5000	Coating Ab	5000 test
ECA 20500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 20500/0500	Conjugated Ab,Peroxidase	500 test
ECP 20500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 20500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 20500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 20500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 20500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 20500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 20500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 20500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 20500/0500	Reagent Set,Alk.Phos.	500 test
SRA 20500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 20500/0500	Reagent Set,Peroxidase	500 test
SRP 20500/1000	Reagent Set,Peroxidase	1000 test
SRP 20500/5000	Reagent Set,Peroxidase	5000 test
*ISK 20500/0005	ImmunoStrip Test Kit	5 strips
*ISK 20500/0025	ImmunoStrip Test Kit	25 strips
*ISK 20500/0250	ImmunoStrip Test Kit	250 strips
*STX 20500/0005	ImmunoStrip Test	5 strips
*STX 20500/0025	ImmunoStrip Test	25 strips
*STX 20500/0250	ImmunoStrip Test	250 strips
LNC 20500	Negative Control	
LPC 20500	Positive Control	

* ImmunoStrip Test kit

Including strips and buffer filled mesh bags

*** INSV ImmunoStrip Buffer**

INSV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

CMV, INSV, TSWV, TMV ImmunoStrip Comb

*ISK 37800/0024	ImmunoStrip Comb Kit, 24 combs and 25 buffer filled mesh bags	
*STX 37800/0024	ImmunoStrip Comb	24 combs

* ImmunoStrip Comb test system for CMV, INSV, TSWV, TMV are packaged in sets of 24. ImmunoStrip combs must be used with SEB1 buffer. Please see accessories for a listing of SEB1 buffer products. ImmunoStrip kits include sample extract bags.

Impatiens Necrotic Spot Virus

INSV	TAS ELISA	
KS 1413 - 0500R	Reagent Set	500 test
KS 1413 - 1000R	Reagent Set	1000 test
KS 1413 - 5000R	Reagent Set	5000 test
KS 1413 - P	Positive Control	
KS 1413 - N	Negative Control	

Indian Cassava Mosaic Virus

ICMV	DAS ELISA	
KS 1385 - 0500R	Reagent Set	500 test
KS 1385 - 1000R	Reagent Set	1000 test
KS 1385 - 5000R	Reagent Set	5000 test
KS 1385 - P	Positive Control	
KS 1385 - N	Negative Control	

Iris Fulva Mosaic Virus

IFMV	DAS ELISA	
KS 1121 - 0500R	Reagent Set	500 test
KS 1121 - 1000R	Reagent Set	1000 test
KS 1121 - 5000R	Reagent Set	5000 test
KS 1121 - P	Positive Control	
KS 1121 - N	Negative Control	

Iris Mild Mosaic Virus

IMMV	DAS ELISA	
KS 1122 - 0500R	Reagent Set	500 test
KS 1122 - 1000R	Reagent Set	1000 test
KS 1122 - 5000R	Reagent Set	5000 test
KS 1122 - P	Positive Control	
KS 1122 - N	Negative Control	

Iris Severe Mosaic Virus

ISMV	DAS ELISA	
KS 1123 - 0500R	Reagent Set	500 test

KS 1123 - 1000R	Reagent Set	1000 test
KS 1123 - 5000R	Reagent Set	5000 test
KS 1123 - P	Positive Control	
KS 1123 - N	Negative Control	

Iris Yellow Spot Virus

IYSV	DAS ELISA	Poly / Poly
CAB 60500/0500	Coating Ab	500 test
CAB 60500/1000	Coating Ab	1000 test
CAB 60500/5000	Coating Ab	5000 test
ECA 60500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 60500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 60500/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 60500/0500	Reagent Set,Alk.Phos.	500 test
SRA 60500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 60500/5000	Reagent Set,Alk.Phos.	5000 test
LNC 60500	Negative Control	

Jhonson Grass Chlorotic Strip Mosaic Virus

JGCSMV	DAS ELISA	
KS 1386 - 0500R	Reagent Set	500 test
KS 1386 - 1000R	Reagent Set	1000 test
KS 1386 - 5000R	Reagent Set	5000 test
KS 1386 - P	Positive Control	
KS 1386 - N	Negative Control	

Jhonson Grass Mosaic Virus

JGMV	DAS ELISA	
KS 1125 - 0500R	Reagent Set	500 test
KS 1125 - 1000R	Reagent Set	1000 test
KS 1125 - 5000R	Reagent Set	5000 test
KS 1125 - P	Positive Control	
KS 1125 - N	Negative Control	

Jhonsongrass Mosaic Virus

JgMV	Compound Direct ELISA	Poly / Poly
CAB 18300/0500	Coating Ab	500 test
CAB 18300/1000	Coating Ab	1000 test
CAB 18300/5000	Coating Ab	5000 test
ECA 18300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18300/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 18300/0500	Reagent Set,Alk.Phos.	500 test
SRA 18300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18300/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18300	Negative Control	
LPC 18300	Positive Control	

Kalanchoe Latent Virus

KLV	DAS ELISA	Poly/Poly
CAB 62000/0500	Coating Ab	500 test
CAB 62000/1000	Coating Ab	1000 test
CAB 62000/5000	Coating Ab	5000 test
ECA 62000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 62000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 62000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 62000/0500	Conjugated Ab,Peroxidase	500 test
ECP 62000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 62000/5000	Conjugated Ab,Peroxidase	5000 test
SRA 62000/0500	Reagent Set,Alk.Phos.	500 test
SRA 62000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 62000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 62000/0500	Reagent Set,Peroxidase	500 test
SRP 62000/1000	Reagent Set,Peroxidase	1000 test
SRP 62000/5000	Reagent Set,Peroxidase	5000 test
LNC 62000	Negative Control	
LPC 62000	Positive Control	

Kyuri Green Mottle Mosaic Virus

KGMMV	DAS ELISA	Poly/Poly
CAB 65000/0500	Coating Ab	500 test
CAB 65000/1000	Coating Ab	1000 test
CAB 65000/5000	Coating Ab	5000 test
ECA 65000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 65000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 65000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 65000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 65000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 65000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 65000/0500	Reagent Set,Alk.Phos.	500 test
SRA 65000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 65000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 65000	Negative Control	
LPC 65000	Positive Control	

Leek Yellow Stripe Virus

LYSV	DAS ELISA	Poly/Poly
CAB 20301/0500	Coating Ab	500 test
CAB 20301/1000	Coating Ab	1000 test
CAB 20301/5000	Coating Ab	5000 test
ECA 20301/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20301/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20301/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 20301/0500	Reagent Set,Alk.Phos.	500 test
SRA 20301/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20301/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20301	Negative Control	

Lettuce Mosaic Virus

LMV	DAS ELISA	Mono/Mono
CAB 45800/0500	Coating Ab	500 test
CAB 45800/1000	Coating Ab	1000 test
CAB 45800/5000	Coating Ab	5000 test
ECA 45800/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 45800/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 45800/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 45800/0500	Conjugated Ab,Peroxidase	500 test
ECP 45800/1000	Conjugated Ab,Peroxidase	1000 test
ECP 45800/5000	Conjugated Ab,Peroxidase	5000 test
PSA 45800/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 45800/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 45800/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 45800/0096	PathoScreen Kit,Peroxidase	96 test
PSP 45800/0288	PathoScreen Kit,Peroxidase	288 test
PSP 45800/0480	PathoScreen Kit,Peroxidase	480 test
SRA 45800/0500	Reagent Set,Alk.Phos.	500 test
SRA 45800/1000	Reagent Set,Alk.Phos.	1000 test
SRA 45800/5000	Reagent Set,Alk.Phos.	5000 test
SRP 45800/0500	Reagent Set,Peroxidase	500 test
SRP 45800/1000	Reagent Set,Peroxidase	1000 test
SRP 45800/5000	Reagent Set,Peroxidase	5000 test
LNC 45800	Negative Control	
LPC 45800	Positive Control	

Lilac Chlorotic Leafspot Virus

LCLV	DAS ELISA	
KS 1343 - 0500R	Reagent Set	500 test
KS 1343 - 1000R	Reagent Set	1000 test
KS 1343 - 5000R	Reagent Set	5000 test
KS 1343 - P	Positive Control	
KS 1343 - N	Negative Control	

Lily Mottle Virus

LMoV	DAS ELISA	
KS 1130 - 0500R	Reagent Set	500 test
KS 1130 - 1000R	Reagent Set	1000 test
KS 1130 - 5000R	Reagent Set	5000 test
KS 1130 - P	Positive Control	
KS 1130 - N	Negative Control	

Lily Symptomless Virus

LSV	DAS ELISA	Poly/Poly
CAB 62300/0500	Coating Ab	500 test
CAB 62300/1000	Coating Ab	1000 test
CAB 62300/5000	Coating Ab	5000 test
ECA 62300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 62300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 62300/5000	Conjugated Ab,Alk.Phos.	5000 test

ECP 62300/0500	Conjugated Ab,Peroxidase	500 test
ECP 62300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 62300/5000	Conjugated Ab,Peroxidase	5000 test
PSA 62300/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 62300/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 62300/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 62300/0096	PathoScreen Kit,Peroxidase	96 test
PSP 62300/0288	PathoScreen Kit,Peroxidase	288 test
PSP 62300/0480	PathoScreen Kit,Peroxidase	480 test
SRA 62300/0500	Reagent Set,Alk.Phos.	500 test
SRA 62300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 62300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 62300/0500	Reagent Set,Peroxidase	500 test
SRP 62300/1000	Reagent Set,Peroxidase	1000 test
SRP 62300/5000	Reagent Set,Peroxidase	5000 test
LNC 62300	Negative Control	
LPC 62300	Positive Control	

Lily Virus X

LVX	DAS ELISA	
KS 1132 - 0500R	Reagent Set	500 test
KS 1132 - 1000R	Reagent Set	1000 test
KS 1132 - 5000R	Reagent Set	5000 test
KS 1132 - P	Positive Control	
KS 1132 - N	Negative Control	

Line Pattern Virus

LPV	DAS ELISA	
KS 1133 - 0500R	Reagent Set	500 test
KS 1133 - 1000R	Reagent Set	1000 test
KS 1133 - 5000R	Reagent Set	5000 test
KS 1133 - P	Positive Control	
KS 1133 - N	Negative Control	

Maize Chlorotic Dwarf Virus

MCDV	DAS ELISA	
KS 1387 - 0500R	Reagent Set	500 test
KS 1387 - 1000R	Reagent Set	1000 test
KS 1387 - 5000R	Reagent Set	5000 test
KS 1387 - P	Positive Control	
KS 1387 - N	Negative Control	

Maize Chlorotic Mottle Virus

MCMV	DAS ELISA	Poly/Poly
CAB 17001/0500	Coating Ab	500 test
CAB 17001/1000	Coating Ab	1000 test
CAB 17001/5000	Coating Ab	5000 test
ECA 17001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 17001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 17001/5000	Conjugated Ab,Alk.Phos.	5000 test

PSA 17001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 17001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 17001/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 17001/0500	Reagent Set,Alk.Phos.	500 test
SRA 17001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 17001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 17001	Negative Control	
LPC 17001	Positive Control	

Maize Dwarf Mosaic Virus

MDMV	DAS ELISA	Poly/Poly
CAB 18000/0500	Coating Ab	500 test
CAB 18000/1000	Coating Ab	1000 test
CAB 18000/5000	Coating Ab	5000 test
ECA 18000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18000/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 18000/0500	Reagent Set,Alk.Phos.	500 test
SRA 18000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18000	Negative Control	
LPC 18000	Positive Control	

**Maize Dwarf Mosaic Virus - strain A
(= Maize Dwarf Mosaic Virus)****Maize Dwarf Mosaic Virus - strain B
(= Sugarcane Mosaic Virus)****Maize Dwarf Mosaic Virus - strain O
(= Johnsongrass Mosaic Virus)****Maize Mosaic Virus**

MMV	DAS ELISA	Poly/Poly
CAB 18900/0500	Coating Ab	500 test
CAB 18900/1000	Coating Ab	1000 test
CAB 18900/5000	Coating Ab	5000 test
ECA 18900/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18900/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18900/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 18900/0500	Reagent Set,Alk.Phos.	500 test
SRA 18900/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18900/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18900	Negative Control	
LPC 18900	Positive Control	

Maize Rayado Fino Virus

MRFV	DAS ELISA	
KS 1388 - 0500R	Reagent Set	500 test
KS 1388 - 1000R	Reagent Set	1000 test

식물병 Viruses

KS 1388 - 5000R	Reagent Set	5000 test
KS 1388 - P	Positive Control	
KS 1388 - N	Negative Control	

Maize Rough Dwarf Virus

MRDV	DAS ELISA	
KS 1142 - 0500R	Reagent Set	500 test
KS 1142 - 1000R	Reagent Set	1000 test
KS 1142 - 5000R	Reagent Set	5000 test
KS 1142 - P	Positive Control	
KS 1142 - N	Negative Control	

Maize Streak Virus

MSV	DAS ELISA	
KS 1143 - 0500R	Reagent Set	500 test
KS 1143 - 1000R	Reagent Set	1000 test
KS 1143 - 5000R	Reagent Set	5000 test
KS 1143 - P	Positive Control	
KS 1143 - N	Negative Control	

Maize Streak Virus

MSV	Indirect ELISA	
KS 1414 - 0500R	Reagent Set	500 test
KS 1414 - 1000R	Reagent Set	1000 test
KS 1414 - 5000R	Reagent Set	5000 test
KS 1414 - P	Positive Control	
KS 1414 - N	Negative Control	

Maize Stripe Virus

MSPV	DAS ELISA	Poly/Poly
CAB 17300/0500	Coating Ab	500 test
CAB 17300/1000	Coating Ab	1000 test
CAB 17300/5000	Coating Ab	5000 test
ECA 17300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 17300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 17300/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 17300/0500	Reagent Set,Alk.Phos.	500 test
SRA 17300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 17300/5000	Reagent Set,Alk.Phos.	5000 test
LNC 17300	Negative Control	
LPC 17300	Positive Control	

Maize White Line Mosaic Virus

MWLMV	DAS ELISA	Poly/Poly
CAB 18200/0500	Coating Ab	500 test
CAB 18200/1000	Coating Ab	1000 test
CAB 18200/5000	Coating Ab	5000 test
ECA 18200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18200/5000	Conjugated Ab,Alk.Phos.	5000 test

SRA 18200/0500	Reagent Set,Alk.Phos.	500 test
SRA 18200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18200/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18200	Negative Control	
LPC 18200	Positive Control	

Maracuja Mosaic Virus

MMV	DAS ELISA	
KS 1344 - 0500R	Reagent Set	500 test
KS 1344 - 1000R	Reagent Set	1000 test
KS 1344 - 5000R	Reagent Set	5000 test
KS 1344 - P	Positive Control	
KS 1344- N	Negative Control	

Moroccan Watermelon Mosaic Virus

MWMV	DAS ELISA	
KS 1345 - 0500R	Reagent Set	500 test
KS 1345 - 1000R	Reagent Set	1000 test
KS 1345 - 5000R	Reagent Set	5000 test
KS 1345 - P	Positive Control	
KS 1345- N	Negative Control	

Melon Necrotic Spot Virus

MNSV	DAS ELISA	Poly/Poly
CAB 12401/0500	Coating Ab	500 test
CAB 12401/1000	Coating Ab	1000 test
CAB 12401/5000	Coating Ab	5000 test
ECA 12401/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 12401/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 12401/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 12401/0500	Reagent Set,Alk.Phos.	500 test
SRA 12401/1000	Reagent Set,Alk.Phos.	1000 test
SRA 12401/5000	Reagent Set,Alk.Phos.	5000 test
LNC 12401	Negative Control	
LPC 12401	Positive Control	

Mirafiori Lettuce Big-Vein Virus

MiLBVV	DAS ELISA	
KS 1346 - 0500R	Reagent Set	500 test
KS 1346 - 1000R	Reagent Set	1000 test
KS 1346 - 5000R	Reagent Set	5000 test
KS 1346 - P	Positive Control	
KS 1346 - N	Negative Control	

Mirafiori Lettuce Mosaic Virus

MiLMV	DAS ELISA	
KS 1389 - 0500R	Reagent Set	500 test
KS 1389 - 1000R	Reagent Set	1000 test
KS 1389 - 5000R	Reagent Set	5000 test
KS 1389 - P	Positive Control	

KS 1389 - N	Negative Control
-------------	------------------

Mirafiori Lettuce Virus

MiLV	DAS ELISA	
KS 1347 - 0500R	Reagent Set	500 test
KS 1347 - 1000R	Reagent Set	1000 test
KS 1347 - 5000R	Reagent Set	5000 test
KS 1347 - P	Positive Control	
KS 1347 - N	Negative Control	

Mushroom Virus 1

MV1	DAS ELISA	
KS 1147 - 0500R	Reagent Set	500 test
KS 1147 - 1000R	Reagent Set	1000 test
KS 1147 - 5000R	Reagent Set	5000 test
KS 1147 - P	Positive Control	
KS 1147 - N	Negative Control	

Mushroom Virus 1

MV1	Direct ELISA	
KS 1397 - 0500R	Reagent Set	500 test
KS 1397 - 1000R	Reagent Set	1000 test
KS 1397 - 5000R	Reagent Set	5000 test
KS 1397 - P	Positive Control	
KS 1397 - N	Negative Control	

Nandina Mosaic Virus

NaMV	DAS ELISA	Poly/Poly
CAB 72900/0500	Coating Ab	500 test
CAB 72900/1000	Coating Ab	1000 test
CAB 72900/5000	Coating Ab	5000 test
ECA 72900/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 72900/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 72900/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 72900/0500	Reagent Set,Alk.Phos.	500 test
SRA 72900/1000	Reagent Set,Alk.Phos.	1000 test
SRA 72900/5000	Reagent Set,Alk.Phos.	5000 test

Narcissus Latent Virus

NLV	DAS ELISA	
KS 1150 - 0500R	Reagent Set	500 test
KS 1150 - 1000R	Reagent Set	1000 test
KS 1150 - 5000R	Reagent Set	5000 test
KS 1150 - P	Positive Control	
KS 1150 - N	Negative Control	

Narcissus Latent Virus

NLV	Indirect ELISA	
KS 1415 - 0500R	Reagent Set	500 test
KS 1415 - 1000R	Reagent Set	1000 test

KS 1415 - 5000R	Reagent Set	5000 test
KS 1415 - P	Positive Control	
KS 1415 - N	Negative Control	

Narcissus Late Season Yellows Virus

NLSYV	Indirect ELISA	
KS 1151 - 0500R	Reagent Set	500 test
KS 1151 - 1000R	Reagent Set	1000 test
KS 1151 - 5000R	Reagent Set	5000 test
KS 1151 - P	Positive Control	
KS 1151 - N	Negative Control	

Narcissus Mosaic Virus

NMV	DAS ELISA	
KS 1152 - 0500R	Reagent Set	500 test
KS 1152 - 1000R	Reagent Set	1000 test
KS 1152 - 5000R	Reagent Set	5000 test
KS 1152 - P	Positive Control	
KS 1152 - N	Negative Control	

Narcissus Mosaic Virus

NMV	Indirect ELISA	
KS 1418 - 0500R	Reagent Set	500 test
KS 1418 - 1000R	Reagent Set	1000 test
KS 1418 - 5000R	Reagent Set	5000 test
KS 1418 - P	Positive Control	
KS 1418 - N	Negative Control	

Narcissus Q Virus

NQV	Indirect ELISA	
KS 1153 - 0500R	Reagent Set	500 test
KS 1153 - 1000R	Reagent Set	1000 test
KS 1153 - 5000R	Reagent Set	5000 test
KS 1153 - P	Positive Control	
KS 1153 - N	Negative Control	

Narcissus Tip Necrosis Virus

NTNV	Indirect ELISA	
KS 1154 - 0500R	Reagent Set	500 test
KS 1154 - 1000R	Reagent Set	1000 test
KS 1154 - 5000R	Reagent Set	5000 test
KS 1154 - P	Positive Control	
KS 1154 - N	Negative Control	

Narcissus Yellow Stripe Virus

NYSV	Indirect ELISA	
KS 1155 - 0500R	Reagent Set	500 test
KS 1155 - 1000R	Reagent Set	1000 test
KS 1155 - 5000R	Reagent Set	5000 test
KS 1155 - P	Positive Control	

KS 1155 - N	Negative Control
-------------	------------------

**Narcissus Virus Set
(NLV+NLSYV+NMV+NQV+NTNV+NYSV)**

NV Set	Indirect ELISA	
KS 1422 - 0500R	Reagent Set	500 test
KS 1422 - 1000R	Reagent Set	1000 test
KS 1422 - 5000R	Reagent Set	5000 test
KS 1422 - P	Positive Control	
KS 1422 - N	Negative Control	

Neckar River Virus

DAS ELISA

KS 1348 - 0500R	Reagent Set	500 test
KS 1348 - 1000R	Reagent Set	1000 test
KS 1348 - 5000R	Reagent Set	5000 test
KS 1348 - P	Positive Control	
KS 1348 - N	Negative Control	

Nerine Latent Virus

NeLV	DAS ELISA	
KS 1156 - 0500R	Reagent Set	500 test
KS 1156 - 1000R	Reagent Set	1000 test
KS 1156 - 5000R	Reagent Set	5000 test
KS 1156 - P	Positive Control	
KS 1156 - N	Negative Control	

Nerine Virus X

NeVX	DAS ELISA	
KS 1157 - 0500R	Reagent Set	500 test
KS 1157 - 1000R	Reagent Set	1000 test
KS 1157 - 5000R	Reagent Set	5000 test
KS 1157 - P	Positive Control	
KS 1157 - N	Negative Control	

Nerine Virus Y

NeVY	DAS ELISA	
KS 1158 - 0500R	Reagent Set	500 test
KS 1158 - 1000R	Reagent Set	1000 test
KS 1158 - 5000R	Reagent Set	5000 test
KS 1158 - P	Positive Control	
KS 1158 - N	Negative Control	

Nerine Yellow Stripe Virus

NeYSV	DAS ELISA	
KS 1159 - 0500R	Reagent Set	500 test
KS 1159 - 1000R	Reagent Set	1000 test
KS 1159 - 5000R	Reagent Set	5000 test
KS 1159 - P	Positive Control	
KS 1159 - N	Negative Control	

Odontoglossum Ringspot Virus

ORSV	DAS ELISA	Poly/Poly
CAB 54300/0500	Coating Ab	500 test
CAB 54300/1000	Coating Ab	1000 test
CAB 54300/5000	Coating Ab	5000 test
ECA 54300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 54300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 54300/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 54300/0500	Conjugated Ab,Peroxidase	500 test
ECP 54300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 54300/5000	Conjugated Ab,Peroxidase	5000 test
PSA 54300/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 54300/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 54300/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 54300/0096	PathoScreen Kit,Peroxidase	96 test
PSP 54300/0288	PathoScreen Kit,Peroxidase	288 test
PSP 54300/0480	PathoScreen Kit,Peroxidase	480 test
SRA 54300/0500	Reagent Set,Alk.Phos.	500 test
SRA 54300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 54300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 54300/0500	Reagent Set,Peroxidase	500 test
SRP 54300/1000	Reagent Set,Peroxidase	1000 test
SRP 54300/5000	Reagent Set,Peroxidase	5000 test
LNC 54300	Negative Control	
LPC 54300	Positive Control	

Okra Mosaic Virus

DAS ELISA

KS 1349 - 0500R	Reagent Set	500 test
KS 1349 - 1000R	Reagent Set	1000 test
KS 1349 - 5000R	Reagent Set	5000 test
KS 1349 - P	Positive Control	
KS 1349 - N	Negative Control	

Onion Yellow Dwarf Virus

OYDV	DAS ELISA	Poly/Poly
CAB 20101/0500	Coating Ab	500 test
CAB 20101/1000	Coating Ab	1000 test
CAB 20101/5000	Coating Ab	5000 test
ECA 20101/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20101/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20101/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 20101/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 20101/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 20101/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 20101/0500	Reagent Set,Alk.Phos.	500 test
SRA 20101/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20101/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20101	Negative Control	
LPC 20101	Positive Control	

Papaya Mosaic Virus

PapMV	DAS ELISA	Poly/Poly
CAB 88500/0500	Coating Ab	500 test
CAB 88500/1000	Coating Ab	1000 test
CAB 88500/5000	Coating Ab	5000 test
ECA 88500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 88500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 88500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 88500/0500	Conjugated Ab,Peroxidase	500 test
ECP 88500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 88500/5000	Conjugated Ab,Peroxidase	5000 test
SRA 88500/0500	Reagent Set,Alk.Phos.	500 test
SRA 88500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 88500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 88500/0500	Reagent Set,Peroxidase	500 test
SRP 88500/1000	Reagent Set,Peroxidase	1000 test
SRP 88500/5000	Reagent Set,Peroxidase	5000 test
LPC 88500	Positive Control	
LNC 88500	Negative Control	

Papaya Ringspot Virus

PRSV	DAS ELISA	Poly/Poly
CAB 53500/0500	Coating Ab	500 test
CAB 53500/1000	Coating Ab	1000 test
CAB 53500/5000	Coating Ab	5000 test
ECA 53500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 53500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 53500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 53500/0500	Conjugated Ab,Peroxidase	500 test
ECP 53500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 53500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 53500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 53500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 53500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 53500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 53500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 53500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 53500/0500	Reagent Set,Alk.Phos.	500 test
SRA 53500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 53500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 53500/0500	Reagent Set,Peroxidase	500 test
SRP 53500/1000	Reagent Set,Peroxidase	1000 test
SRP 53500/5000	Reagent Set,Peroxidase	5000 test
LNC 53500	Negative Control	
LPC 53500	Positive Control	

Paprika Mild Mottle Virus

PaMMV	DAS ELISA	Poly/Poly
KS 1390 - 0500R	Reagent Set	500 test
KS 1390 - 1000R	Reagent Set	1000 test

KS 1390 - 5000R	Reagent Set	5000 test
KS 1390 - P	Positive Control	
KS 1390 - N	Negative Control	

Passiflora Latent Virus

PLV	DAS ELISA	
KS 1350 - 0500R	Reagent Set	500 test
KS 1350 - 1000R	Reagent Set	1000 test
KS 1350 - 5000R	Reagent Set	5000 test
KS 1350 - P	Positive Control	
KS 1350 - N	Negative Control	

Pea Early-Browning Virus

PEBV	DAS ELISA	
KS 1164 - 0500R	Reagent Set	500 test
KS 1164 - 1000R	Reagent Set	1000 test
KS 1164 - 5000R	Reagent Set	5000 test
KS 1164 - P	Positive Control	
KS 1164 - N	Negative Control	

Pea Enation Mosaic Virus

PEMV	DAS ELISA	
KS 1165 - 0500R	Reagent Set	500 test
KS 1165 - 1000R	Reagent Set	1000 test
KS 1165 - 5000R	Reagent Set	5000 test
KS 1165 - P	Positive Control	
KS 1165 - N	Negative Control	

Pea Mosaic Virus

PMV	DAS ELISA	
KS 1166 - 0500R	Reagent Set	500 test
KS 1166 - 1000R	Reagent Set	1000 test
KS 1166 - 5000R	Reagent Set	5000 test
KS 1166 - P	Positive Control	
KS 1166 - N	Negative Control	

Pea Seed-Borne Mosaic Virus

PSbMV	DAS ELISA	
KS 1167 - 0500R	Reagent Set	500 test
KS 1167 - 1000R	Reagent Set	1000 test
KS 1167 - 5000R	Reagent Set	5000 test
KS 1167 - P	Positive Control	
KS 1167 - N	Negative Control	

Peach Rosette Mosaic Virus

PRMV	DAS ELISA	Poly/Poly
CAB 87700/0500	Coating Ab	500 test
CAB 87700/1000	Coating Ab	1000 test
CAB 87700/5000	Coating Ab	5000 test

식물병 Viruses

ECA 87700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 87700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 87700/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 87700/0500	Reagent Set,Alk.Phos.	500 test
SRA 87700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 87700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 87700	Negative Control	
LPC 87700	Positive Control	

Peanut Mottle Virus

PMoV	DAS ELISA	
KS 1351 - 0500R	Reagent Set	500 test
KS 1351 - 1000R	Reagent Set	1000 test
KS 1351 - 5000R	Reagent Set	5000 test
KS 1351 - P	Positive Control	
KS 1351 - N	Negative Control	

Peanut Stunt Virus

PSV	DAS ELISA	Poly/Poly
CAB 44600/0500	Coating Ab	500 test
CAB 44600/1000	Coating Ab	1000 test
CAB 44600/5000	Coating Ab	5000 test
ECA 44600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44600/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 44600/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 44600/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 44600/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 44600/0500	Reagent Set,Alk.Phos.	500 test
SRA 44600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 44600/5000	Reagent Set,Alk.Phos.	5000 test
LNC 44600	Negative Control	
LPC 44600	Positive Control	

Pelargonium Flower Break Virus

PFBV	DAS ELISA	Poly/Poly
CAB 90000/0500	Coating Ab	500 test
CAB 90000/1000	Coating Ab	1000 test
CAB 90000/5000	Coating Ab	5000 test
ECA 90000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 90000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 90000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 90000/0500	Conjugated Ab,Peroxidase	500 test
ECP 90000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 90000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 90000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 90000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 90000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 90000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 90000/0288	PathoScreen Kit,Peroxidase	288 test

PSP 90000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 90000/0500	Reagent Set,Alk.Phos.	500 test
SRA 90000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 90000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 90000/0500	Reagent Set,Peroxidase	500 test
SRP 90000/1000	Reagent Set,Peroxidase	1000 test
SRP 90000/5000	Reagent Set,Peroxidase	5000 test
LNC 90000	Negative Control	
LPC 90000	Positive Control	

Pelargonium Leaf Curl Virus

PLCV	DAS ELISA	Poly/Poly
CAB 92700/0500	Coating Ab	500 test
CAB 92700/1000	Coating Ab	1000 test
CAB 92700/5000	Coating Ab	5000 test
ECA 92700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 92700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 92700/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 92700/0500	Reagent Set,Alk.Phos.	500 test
SRA 92700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 92700/5000	Reagent Set,Alk.Phos.	5000 test

Pelargonium Line Pattern Virus

PLPV	DAS ELISA	
KS 1172 - 0500R	Reagent Set	500 test
KS 1172 - 1000R	Reagent Set	1000 test
KS 1172 - 5000R	Reagent Set	5000 test
KS 1172 - P	Positive Control	
KS 1172 - N	Negative Control	

Pelargonium Ringspot Virus

(See Elderberry Latent Virus)

Pelargonium Zonate Spot Virus

PZSV	DAS ELISA	Poly/Poly
CAB 90100/0500	Coating Ab	500 test
CAB 90100/1000	Coating Ab	1000 test
CAB 90100/5000	Coating Ab	5000 test
ECA 90100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 90100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 90100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 90100/0500	Reagent Set,Alk.Phos.	500 test
SRA 90100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 90100/5000	Reagent Set,Alk.Phos.	5000 test
LPC 90100	Positive Control	
LNC 90100	Negative Control	

Pepino Mosaic Virus

PepMV	DAS ELISA	poly/Poly
CAB 13000/0500	Coating Ab	500 test
CAB 13000/1000	Coating Ab	1000 test
CAB 13000/5000	Coating Ab	5000 test
ECA 13000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 13000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 13000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 13000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 13000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 13000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 13000/0500	Reagent Set,Alk.Phos.	500 test
SRA 13000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 13000/5000	Reagent Set,Alk.Phos.	5000 test
*ISK 13000/0005	ImmunoStrip Test Kit	5 strips
*ISK 13000/0025	ImmunoStrip Test Kit	25 strips
*ISK 13000/0250	ImmunoStrip Test Kit	250 strips
*STX 13000/0005	ImmunoStrip Test	5 strips
*STX 13000/0025	ImmunoStrip Test	25 strips
*STX 13000/0250	ImmunoStrip Test	250 strips
LPC 13000	Positive Control	
LNC 13000	Negative Control	

* ImmunoStrip Test kit

Including strips and buffer filled mesh bags

* PepMV ImmunoStrip Buffer

PepMV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

Pepino Mosaic Virus

PepMV	TAS ELISA	
KS 1424 - 0500R	Reagent Set	500 test
KS 1424 - 1000R	Reagent Set	1000 test
KS 1424 - 5000R	Reagent Set	5000 test
KS 1424 - P	Positive Control	
KS 1424 - N	Negative Control	

Pepper Mild Mottle Virus

PMMoV	DAS ELISA	
KS 1175 - 0500R	Reagent Set	500 test
KS 1175 - 1000R	Reagent Set	1000 test
KS 1175 - 5000R	Reagent Set	5000 test
KS 1175 - P	Positive Control	
KS 1175 - N	Negative Control	

Pepper Mild Mottle Virus

PMMoV	Compound Direct ELISA	Poly/Mono
CAB 37500/0500	Coating Ab	500 test
CAB 37500/1000	Coating Ab	1000 test
CAB 37500/5000	Coating Ab	5000 test

ECA 37500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 37500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 37500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 37500/0500	Conjugated Ab,Peroxidase	500 test
ECP 37500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 37500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 37500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 37500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 37500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 37500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 37500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 37500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 37500/0500	Reagent Set,Alk.Phos.	500 test
SRA 37500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 37500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 37500/0500	Reagent Set,Peroxidase	500 test
SRP 37500/1000	Reagent Set,Peroxidase	1000 test
SRP 37500/5000	Reagent Set,Peroxidase	5000 test
LNC 37500	Negative Control	
LPC 37500	Positive Control	

Pepper Mottle Virus

PepMoV	DAS ELISA	Poly/Poly
CAB 38000/0500	Coating Ab	500 test
CAB 38000/1000	Coating Ab	1000 test
CAB 38000/5000	Coating Ab	5000 test
ECA 38000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 38000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 38000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 38000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 38000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 38000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 38000/0500	Reagent Set,Alk.Phos.	500 test
SRA 38000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 38000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 38000	Negative Control	
LPC 38000	Positive Control	

Pepper Veinal Mottle Virus

PVMoV	DAS ELISA	
KS 1178 - 0500R	Reagent Set	500 test
KS 1178 - 1000R	Reagent Set	1000 test
KS 1178 - 5000R	Reagent Set	5000 test
KS 1178 - P	Positive Control	
KS 1178 - N	Negative Control	

Peru Tomato Mosaic Virus

PToMV	DAS ELISA	
KS 1352 - 0500R	Reagent Set	500 test
KS 1352 - 1000R	Reagent Set	1000 test

식물병 Viruses

KS 1352 - 5000R	Reagent Set	5000 test
KS 1352 - P	Positive Control	
KS 1352 - N	Negative Control	

Petunia Asteroid Mosaic Virus

PeAMV	DAS ELISA	
KS 1179 - 0500R	Reagent Set	500 test
KS 1179 - 1000R	Reagent Set	1000 test
KS 1179 - 5000R	Reagent Set	5000 test
KS 1179 - P	Positive Control	
KS 1179 - N	Negative Control	

Pineapple Mealybug Wilt associated Virus - 1

PMWaV-1	Tissue Blot	
CCS 12700	Monoclonal Ab supernatant, 10 ml	

Pineapple Mealybug Wilt associated Virus - 2

PMWaV-2	Tissue Blot	
CCS 12800	Monoclonal Ab supernatant, 10 ml	

Plum Pox Virus

PPV	DAS ELISA	
KS 1180 - 0500R	Reagent Set	500 test
KS 1180 - 1000R	Reagent Set	1000 test
KS 1180 - 5000R	Reagent Set	5000 test
KS 1180 - P	Positive Control	
KS 1180 - N	Negative Control	

Plum Pox Virus

PPV	Compound Direct ELISA	Mono/Poly
DRA 31503/1000	Durviz Reagent Set	1000 test
DRA 31503/5000	Durviz Reagent Set	5000 test
SPC 31503/0020	Positive control strip	20 strips
SPC 31503/0050	Positive control strip	50 strips

Plum Pox Virus

PPV	TAS ELISA	
KS 1181 - 0500R	Reagent Set	500 test
KS 1181 - 1000R	Reagent Set	1000 test
KS 1181 - 5000R	Reagent Set	5000 test
KS 1181 - P	Positive Control	
KS 1181 - N	Negative Control	

Poinsettia Cryptic Virus

PnCV	DAS ELISA	
KS 1353 - 0500R	Reagent Set	500 test
KS 1353 - 1000R	Reagent Set	1000 test
KS 1353 - 5000R	Reagent Set	5000 test
KS 1353 - P	Positive Control	
KS 1353 - N	Negative Control	

Poinsettia Mosaic Virus

PnMV	DAS ELISA	Poly/Poly
CAB 90700/0500	Coating Ab	500 test
CAB 90700/1000	Coating Ab	1000 test
CAB 90700/5000	Coating Ab	5000 test
ECA 90700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 90700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 90700/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 90700/0500	Conjugated Ab,Peroxidase	500 test
ECP 90700/1000	Conjugated Ab,Peroxidase	1000 test
ECP 90700/5000	Conjugated Ab,Peroxidase	5000 test
SRA 90700/0500	Reagent Set, Alk.Phos.	500 test
SRA 90700/1000	Reagent Set, Alk.Phos.	1000 test
SRA 90700/5000	Reagent Set, Alk.Phos.	5000 test
SRP 90700/0500	Reagent Set, Peroxidase	500 test
SRP 90700/1000	Reagent Set, Peroxidase	1000 test
SRP 90700/5000	Reagent Set, Peroxidase	5000 test
LPC 90700	Positive Control	

Poplar Mosaic Virus

PopMV	DAS ELISA	
KS 1184 - 0500R	Reagent Set	500 test
KS 1184 - 1000R	Reagent Set	1000 test
KS 1184 - 5000R	Reagent Set	5000 test
KS 1184 - P	Positive Control	
KS 1184 - N	Negative Control	

Potato Andean Mottle Virus (See Andean Potato Mottle Virus)

Potato Aucuba Mosaic Virus

PAMV	DAS ELISA	Poly/Poly
CAB 21000/0500	Coating Ab	500 test
CAB 21000/1000	Coating Ab	1000 test
CAB 21000/5000	Coating Ab	5000 test
ECA 21000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 21000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 21000/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 21000/0500	Reagent Set,Alk.Phos.	500 test
SRA 21000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 21000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 21000	Negative Control	
LPC 21000	Positive Control	

Potato Black Ringspot Virus

PBRV	DAS ELISA	
KS 1186 - 0500R	Reagent Set	500 test
KS 1186 - 1000R	Reagent Set	1000 test
KS 1186 - 5000R	Reagent Set	5000 test
KS 1186 - P	Positive Control	

KS 1186 - N	Negative Control
-------------	------------------

Potato Latent Virus

PotLV	Compound Direct ELISA	Poly/Mono
CAB 40501/0500	Coating Ab	500 test
CAB 40501/1000	Coating Ab	1000 test
CAB 40501/5000	Coating Ab	5000 test
ECA 40501/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 40501/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 40501/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 40501/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 40501/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 40501/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 40501/0500	Reagent Set,Alk.Phos.	500 test
SRA 40501/1000	Reagent Set,Alk.Phos.	1000 test
SRA 40501/5000	Reagent Set,Alk.Phos.	5000 test
LNC 40501	Negative Control	
LPC 40501	Positive Control	

Potato Leaf Roll Virus

PLRV	DAS ELISA	Poly/Poly
CAB 30002/0500	Coating Ab	500 test
CAB 30002/1000	Coating Ab	1000 test
CAB 30002/5000	Coating Ab	5000 test
ECA 30002/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 30002/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 30002/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 30002/0500	Conjugated Ab,Peroxidase	500 test
ECP 30002/1000	Conjugated Ab,Peroxidase	1000 test
ECP 30002/5000	Conjugated Ab,Peroxidase	5000 test
PSA 30002/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 30002/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 30002/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 30002/0096	PathoScreen Kit,Alk.Peroxidas	96 test
PSP 30002/0288	PathoScreen Kit,Alk.Peroxidas	288 test
PSP 30002/0480	PathoScreen Kit,Alk.Peroxidas	480 test
SRA 30002/0500	Reagent Set,Alk.Phos.	500 test
SRA 30002/1000	Reagent Set,Alk.Phos.	1000 test
SRA 30002/5000	Reagent Set,Alk.Phos.	5000 test
SRP 30002/0500	Reagent Set,Peroxidase	500 test
SRP 30002/1000	Reagent Set,Peroxidase	1000 test
SRP 30002/5000	Reagent Set,Peroxidase	5000 test
LNC 30002	Negative Control	
LPC 30002	Positive Control	

Potato Leaf Roll Virus

PLRV	TAS ELISA	
KS 1189 - 0500R	Reagent Set	500 test
KS 1189 - 1000R	Reagent Set	1000 test
KS 1189 - 5000R	Reagent Set	5000 test

KS 1189 - P	Positive Control
-------------	------------------

KS 1189 - N	Negative Control
-------------	------------------

Potato Mop Top Virus

PMTV	DAS ELISA	
KS 1354 - 0500R	Reagent Set	500 test
KS 1354 - 1000R	Reagent Set	1000 test
KS 1354 - 5000R	Reagent Set	5000 test
KS 1354 - P	Positive Control	
KS 1354 - N	Negative Control	

Potato Mop Top Virus

PMTV	TAS ELISA	
KS 1190 - 0500R	Reagent Set	500 test
KS 1190 - 1000R	Reagent Set	1000 test
KS 1190 - 5000R	Reagent Set	5000 test
KS 1190 - P	Positive Control	
KS 1190 - N	Negative Control	

Potato Virus A

PVA	DAS ELISA	
KS 1191 - 0500R	Reagent Set	500 test
KS 1191 - 1000R	Reagent Set	1000 test
KS 1191 - 5000R	Reagent Set	5000 test
KS 1191 - P	Positive Control	
KS 1191 - N	Negative Control	

Potato Virus A

PVA	Compound Direct ELISA	Poly/Mono
CAB 60000/0500	Coating Ab	500 test
CAB 60000/1000	Coating Ab	1000 test
CAB 60000/5000	Coating Ab	5000 test
ECA 60000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 60000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 60000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 60000/0500	Conjugated Ab,Peroxidase	500 test
ECP 60000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 60000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 60000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 60000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 60000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 60000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 60000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 60000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 60000/0500	Reagent Set,Alk.Phos.	500 test
SRA 60000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 60000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 60000/0500	Reagent Set,Peroxidase	500 test
SRP 60000/1000	Reagent Set,Peroxidase	1000 test
SRP 60000/5000	Reagent Set,Peroxidase	5000 test

식물병 Viruses

LNC 60000	Negative Control
LPC 60000	Positive Control

Potato Virus B - OCA strain

PVB-OCA	DAS ELISA	
KS 1355 - 0500R	Reagent Set	500 test
KS 1355 - 1000R	Reagent Set	1000 test
KS 1355 - 5000R	Reagent Set	5000 test
KS 1355 - P	Positive Control	
KS 1355 - N	Negative Control	

Potato Virus M

PVM	DAS ELISA	Poly/Poly
CAB 50000/0500	Coating Ab	500 test
CAB 50000/1000	Coating Ab	1000 test
CAB 50000/5000	Coating Ab	5000 test
ECA 50000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 50000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 50000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 50000/0500	Conjugated Ab,Peroxidase	500 test
ECP 50000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 50000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 50000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 50000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 50000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 50000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 50000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 50000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 50000/0500	Reagent Set,Alk.Phos.	500 test
SRA 50000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 50000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 50000/0500	Reagent Set,Peroxidase	500 test
SRP 50000/1000	Reagent Set,Peroxidase	1000 test
SRP 50000/5000	Reagent Set,Peroxidase	5000 test
LNC 50000	Negative Control	
LPC 50000	Positive Control	

Potato Virus S

PVS	DAS ELISA	Poly/Poly
CAB 40000/0500	Coating Ab	500 test
CAB 40000/1000	Coating Ab	1000 test
CAB 40000/5000	Coating Ab	5000 test
ECA 40000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 40000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 40000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 40000/0500	Conjugated Ab,Peroxidase	500 test
ECP 40000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 40000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 40000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 40000/0288	PathoScreen Kit,Alk.Phos.	288 test

PSA 40000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 40000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 40000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 40000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 40000/0500	Reagent Set,Alk.Phos.	500 test
SRA 40000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 40000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 40000/0500	Reagent Set,Peroxidase	500 test
SRP 40000/1000	Reagent Set,Peroxidase	1000 test
SRP 40000/5000	Reagent Set,Peroxidase	5000 test
LNC 40000	Negative Control	
LPC 40000	Positive Control	

Potato Virus T

PVT	DAS ELISA	
KS 1195 - 0500R	Reagent Set	500 test
KS 1195 - 1000R	Reagent Set	1000 test
KS 1195 - 5000R	Reagent Set	5000 test
KS 1195 - P	Positive Control	
KS 1195 - N	Negative Control	

Potato Virus T

PVT	Compound Direct ELISA	Poly / Poly
CAB 40100/0500	Coating Ab	500 test
CAB 40100/1000	Coating Ab	1000 test
CAB 40100/5000	Coating Ab	5000 test
ECA 40100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 40100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 40100/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 40100/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 40100/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 40100/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 40100/0500	Reagent Set,Alk.Phos.	500 test
SRA 40100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 40100/5000	Reagent Set,Alk.Phos.	5000 test
LNC 40100	Negative Control	
LPC 40100	Positive Control	

Potato Virus V

PVV	DAS ELISA	Poly / Poly
CAB 10100/0500	Coating Ab	500 test
CAB 10100/1000	Coating Ab	1000 test
CAB 10100/5000	Coating Ab	5000 test
ECA 10100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 10100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 10100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 10100/0500	Reagent Set,Alk.Phos.	500 test
SRA 10100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 10100/5000	Reagent Set,Alk.Phos.	5000 test
LNC 10100	Negative Control	

LPC 10100	Positive Control
-----------	------------------

Potato Virus X

PVX	DAS ELISA	Poly/Poly
CAB 10000/0500	Coating Ab	500 test
CAB 10000/1000	Coating Ab	1000 test
CAB 10000/5000	Coating Ab	5000 test
ECA 10000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 10000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 10000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 10000/0500	Conjugated Ab,Peroxidase	500 test
ECP 10000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 10000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 10000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 10000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 10000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 10000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 10000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 10000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 10000/0500	Reagent Set,Alk.Phos.	500 test
SRA 10000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 10000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 10000/0500	Reagent Set,Peroxidase	500 test
SRP 10000/1000	Reagent Set,Peroxidase	1000 test
SRP 10000/5000	Reagent Set,Peroxidase	5000 test
LNC 10000	Negative Control	
LPC 10000	Positive Control	

Potato Virus X

PVX	TAS ELISA	
KS 1198 - 0500R	Reagent Set	500 test
KS 1198 - 1000R	Reagent Set	1000 test
KS 1198 - 5000R	Reagent Set	5000 test
KS 1198 - P	Positive Control	
KS 1198 - N	Negative Control	

Potato Virus X - HB strain

PVX-HB	TAS ELISA	
KS 1199 - 0500R	Reagent Set	500 test
KS 1199 - 1000R	Reagent Set	1000 test
KS 1199 - 5000R	Reagent Set	5000 test
KS 1199 - P	Positive Control	
KS 1199 - N	Negative Control	

Potato Virus X+S

PVX+PVS	DAS ELISA	
KS 1200 - 0500R	Reagent Set	500 test
KS 1200 - 1000R	Reagent Set	1000 test
KS 1200 - 5000R	Reagent Set	5000 test
KS 1200 - P	Positive Control	

KS 1200 - N	Negative Control
-------------	------------------

Potato Virus Y

PVY	DAS ELISA	Poly/Poly
CAB 20000/0500	Coating Ab	500 test
CAB 20000/1000	Coating Ab	1000 test
CAB 20000/5000	Coating Ab	5000 test
ECP 20000/0500	Conjugated Ab,Peroxidase	500 test
ECP 20000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 20000/5000	Conjugated Ab,Peroxidase	5000 test
PSP 20000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 20000/0288	PathoScreen Kit, .Peroxidase	288 test
PSP 20000/0480	PathoScreen Kit, .Peroxidase	480 test
SRP 20000/0500	Reagent Set,Peroxidase	500 test
SRP 20000/1000	Reagent Set,Peroxidase	1000 test
SRP 20000/5000	Reagent Set,Peroxidase	5000 test
LNC 20000	Negative Control	
LPC 20000	Positive Control	

Potato Virus Y

PVY	Compound Direct ELISA	Poly/Mono
CAB 20001/0500	Coating Ab	500 test
CAB 20001/1000	Coating Ab	1000 test
CAB 20001/5000	Coating Ab	5000 test
ECA 20001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20001/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 20001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 20001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 20001/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 20001/0500	Reagent Set,Alk.Phos.	500 test
SRA 20001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20001	Negative Control	
LPC 20001	Positive Control	

Potato Virus Y - c

PVY-c	DAS ELISA	
KS 1356 - 0500R	Reagent Set	500 test
KS 1356 - 1000R	Reagent Set	1000 test
KS 1356 - 5000R	Reagent Set	5000 test
KS 1356 - P	Positive Control	
KS 1356 - N	Negative Control	

Potato Virus Y - strain c

PVY-c	Compound Direct ELISA	Poly/Mono
CAB 20700/0500	Coating Ab	500 test
CAB 20700/1000	Coating Ab	1000 test
CAB 20700/5000	Coating Ab	5000 test
ECA 20700/0500	Conjugated Ab,Alk.Phos.	500 test

식물병 Viruses

ECA 20700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 20700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 20700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 20700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 20700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 20700/0500	Reagent Set,Alk.Phos.	500 test
SRA 20700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20700	Negative Control	
LPC 20700	Positive Control	

Potato Virus Y - necrotic strain

PVY-n	DAS ELISA	
KS 1205 - 0500R	Reagent Set	500 test
KS 1205 - 1000R	Reagent Set	1000 test
KS 1205 - 5000R	Reagent Set	5000 test
KS 1205 - P	Positive Control	
KS 1205 - N	Negative Control	

Potato Virus Y - necrotic strain

PVY-n	TAS ELISA	
KS 1206 - 0500R	Reagent Set	500 test
KS 1206 - 1000R	Reagent Set	1000 test
KS 1206 - 5000R	Reagent Set	5000 test
KS 1206 - P	Positive Control	
KS 1206 - N	Negative Control	

Potato Virus Y - strain o

PVY-o	TAS ELISA	
KS 1208 - 0500R	Reagent Set	500 test
KS 1208 - 1000R	Reagent Set	1000 test
KS 1208 - 5000R	Reagent Set	5000 test
KS 1208 - P	Positive Control	
KS 1208 - N	Negative Control	

Potato Virus Y - strain o, c

PVY-o+c	DAS ELISA	
KS 1209 - 0500R	Reagent Set	500 test
KS 1209 - 1000R	Reagent Set	1000 test
KS 1209 - 5000R	Reagent Set	5000 test
KS 1209 - P	Positive Control	
KS 120 - N	Negative Control	

Potato Virus Y - strain o, c

PVY-o+c	Compound Direct ELISA	Poly/Mono
CAB 20600/0500	Coating Ab	500 test
CAB 20600/1000	Coating Ab	1000 test
CAB 20600/5000	Coating Ab	5000 test
ECA 20600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 20600/1000	Conjugated Ab,Alk.Phos.	1000 test

ECA 20600/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 20600/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 20600/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 20600/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 20600/0500	Reagent Set,Alk.Phos.	500 test
SRA 20600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 20600/5000	Reagent Set,Alk.Phos.	5000 test
LNC 20600	Negative Control	
LPC 20600	Positive Control	

Potato Virus Y - strain o, c, n

PVY-o+c+n	DAS ELISA	
KS 1211 - 0500R	Reagent Set	500 test
KS 1211 - 1000R	Reagent Set	1000 test
KS 1211 - 5000R	Reagent Set	5000 test
KS 1211 - P	Positive Control	
KS 1211 - N	Negative Control	

Potato Virus Y - strain o, c, n

PVY-o+c+n	TAS ELISA	
KS 1212 - 0500R	Reagent Set	500 test
KS 1212 - 1000R	Reagent Set	1000 test
KS 1212 - 5000R	Reagent Set	5000 test
KS 1212 - P	Positive Control	
KS 1212 - N	Negative Control	

Potato Yellow Dwarf Virus

PYDV	DAS ELISA	
KS 1213 - 0500R	Reagent Set	500 test
KS 1213 - 1000R	Reagent Set	1000 test
KS 1213 - 5000R	Reagent Set	5000 test
KS 1213 - P	Positive Control	
KS 1213 - N	Negative Control	

Potato Yellowing Virus

PYV	DAS ELISA	
KS 1391 - 0500R	Reagent Set	500 test
KS 1391 - 1000R	Reagent Set	1000 test
KS 1391 - 5000R	Reagent Set	5000 test
KS 1391 - P	Positive Control	
KS 1391 - N	Negative Control	

Potyvirus Group (Aphid Transmitted)

POTV	DAS ELISA	
KS 1215 - 0500R	Reagent Set	500 test
KS 1215 - 1000R	Reagent Set	1000 test
KS 1215 - 5000R	Reagent Set	5000 test
KS 1215 - P	Positive Control	
KS 1215 - N	Negative Control	

Potyvirus Group

POTY	Indirect ELISA	-- / Mono
CAB 27200/0500	Coating Ab	500 test
CAB 27200/1000	Coating Ab	1000 test
CAB 27200/5000	Coating Ab	5000 test
ECA 27200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 27200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 27200/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 27200/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 27200/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 27200/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 27200/0500	Reagent Set,Alk.Phos.	500 test
SRA 27200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 27200/5000	Reagent Set,Alk.Phos.	5000 test
LNC 27200	Negative Control	
LPC 27200	Positive Control	

Prune Dwarf Virus

PDV	DAS ELISA	Poly/Poly
CAB 98700/0500	Coating Ab	500 test
CAB 98700/1000	Coating Ab	1000 test
CAB 98700/5000	Coating Ab	5000 test
ECA 98700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 98700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 98700/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 98700/0500	Conjugated Ab,Peroxidase	500 test
ECP 98700/1000	Conjugated Ab,Peroxidase	1000 test
ECP 98700/5000	Conjugated Ab,Peroxidase	5000 test
PSA 98700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 98700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 98700/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 98700/0096	PathoScreen Kit,Peroxidase	96 test
PSP 98700/0288	PathoScreen Kit,Peroxidase	288 test
PSP 98700/0480	PathoScreen Kit,Peroxidase	480 test
SRA 98700/0500	Reagent Set,Alk.Phos.	500 test
SRA 98700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 98700/5000	Reagent Set,Alk.Phos.	5000 test
SRP 98700/0500	Reagent Set,Peroxidase	500 test
SRP 98700/1000	Reagent Set,Peroxidase	1000 test
SRP 98700/5000	Reagent Set,Peroxidase	5000 test
LNC 98700	Negative Control	
LPC 98700	Positive Control	

Prunus Necrotic Ringspot Virus

PNRSV	DAS ELISA	
KS 1217 - 0500R	Reagent Set	500 test
KS 1217 - 1000R	Reagent Set	1000 test
KS 1217 - 5000R	Reagent Set	5000 test
KS 1217 - P	Positive Control	
KS 1217 - N	Negative Control	

Prunus Necrotic Ringspot Virus

PNRSV	Compound Direct ELISA	Poly/Mono
CAB 30500/0500	Coating Ab	500 test
CAB 30500/1000	Coating Ab	1000 test
CAB 30500/5000	Coating Ab	5000 test
ECA 30500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 30500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 30500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 30500/0500	Conjugated Ab,Peroxidase	500 test
ECP 30500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 30500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 30500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 30500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 30500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 30500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 30500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 30500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 30500/0500	Reagent Set,Alk.Phos.	500 test
SRA 30500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 30500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 30500/0500	Reagent Set,Peroxidase	500 test
SRP 30500/1000	Reagent Set,Peroxidase	1000 test
SRP 30500/5000	Reagent Set,Peroxidase	5000 test
LNC 30500	Negative Control	
LPC 30500	Positive Control	

Radish Mosaic Virus

RaMV	DAS ELISA	
KS 1219 - 0500R	Reagent Set	500 test
KS 1219 - 1000R	Reagent Set	1000 test
KS 1219 - 5000R	Reagent Set	5000 test
KS 1219 - P	Positive Control	
KS 1219 - N	Negative Control	

Raspberry Bushy Dwarf Virus

RBDV	DAS ELISA	
KS 1220 - 0500R	Reagent Set	500 test
KS 1220 - 1000R	Reagent Set	1000 test
KS 1220 - 5000R	Reagent Set	5000 test
KS 1220 - P	Positive Control	
KS 1220 - N	Negative Control	

Raspberry Bushy Dwarf Virus

RBDV	Compound Direct ELISA	Poly/Mono
CAB 58200/0500	Coating Ab	500 test
CAB 58200/1000	Coating Ab	1000 test
CAB 58200/5000	Coating Ab	5000 test
ECA 58200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 58200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 58200/5000	Conjugated Ab,Alk.Phos.	5000 test

식물병 Viruses

PSA 58200/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 58200/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 58200/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 58200/0500	Reagent Set,Alk.Phos.	500 test
SRA 58200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 58200/5000	Reagent Set,Alk.Phos.	5000 test
LPC 58200	Positive Control	
LNC 58200	Negative Control	

Raspberry Ringspot Virus

RpRSV	DAS ELISA	Poly/Poly
CAB 59401/0500	Coating Ab	500 test
CAB 59401/1000	Coating Ab	1000 test
CAB 59401/5000	Coating Ab	5000 test
ECA 59401/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 59401/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 59401/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 59401/0500	Reagent Set,Alk.Phos.	500 test
SRA 59401/1000	Reagent Set,Alk.Phos.	1000 test
SRA 59401/5000	Reagent Set,Alk.Phos.	5000 test
LNC 59401	Negative Control	
LPC 59401	Positive Control	

Raspberry Ringspot Virus - cherry isolate

RpRSV-cherry	DAS ELISA	
KS 1224 - 0500R	Reagent Set	500 test
KS 1224 - 1000R	Reagent Set	1000 test
KS 1224 - 5000R	Reagent Set	5000 test
KS 1224 - P	Positive Control	
KS 1224 - N	Negative Control	

Raspberry Ringspot Virus - grapevine isolate

RpRSV-grapevine	DAS ELISA	
KS 1225 - 0500R	Reagent Set	500 test
KS 1225 - 1000R	Reagent Set	1000 test
KS 1225 - 5000R	Reagent Set	5000 test
KS 1225 - P	Positive Control	
KS 1225 - N	Negative Control	

Raspberry Ringspot Virus - raspberry

RpRSV-raspberry	DAS ELISA	
KS 1226 - 0500R	Reagent Set	500 test
KS 1226 - 1000R	Reagent Set	1000 test
KS 1226 - 5000R	Reagent Set	5000 test
KS 1226 - P	Positive Control	
KS 1226 - N	Negative Control	

Red Clover Cryptic Virus 2

RCCV2	DAS ELISA	
KS 1357 - 0500R	Reagent Set	500 test
KS 1357 - 1000R	Reagent Set	1000 test
KS 1357 - 5000R	Reagent Set	5000 test
KS 1357 - P	Positive Control	
KS 1357 - N	Negative Control	

Red Clover Mottle Virus

RCMV	DAS ELISA	
KS 1227 - 0500R	Reagent Set	500 test
KS 1227 - 1000R	Reagent Set	1000 test
KS 1227 - 5000R	Reagent Set	5000 test
KS 1227 - P	Positive Control	
KS 1227 - N	Negative Control	

Red Clover Vein Mosaic Virus

RCVMV	DAS ELISA	
KS 1228 - 0500R	Reagent Set	500 test
KS 1228 - 1000R	Reagent Set	1000 test
KS 1228 - 5000R	Reagent Set	5000 test
KS 1228 - P	Positive Control	
KS 1228 - N	Negative Control	

Red LaSoda Virus

RLSV	DAS ELISA	
KS 1229 - 0500R	Reagent Set	500 test
KS 1229 - 1000R	Reagent Set	1000 test
KS 1229 - 5000R	Reagent Set	5000 test
KS 1229 - P	Positive Control	
KS 1229 - N	Negative Control	

Red LaSoda Virus(See Potato Latent Virus)

Ribgrass Mosaic Virus

RMV	DAS ELISA	Poly/Poly
CAB 57600/0500	Coating Ab	500 test
CAB 57600/1000	Coating Ab	1000 test
CAB 57600/5000	Coating Ab	5000 test
ECA 57600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 57600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 57600/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 57600/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 57600/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 57600/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 57600/0500	Reagent Set,Alk.Phos.	500 test
SRA 57600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 57600/5000	Reagent Set,Alk.Phos.	5000 test
LPC 57600	Positive Control	
LNC 57600	Negative Control	

Rice Dwarf Virus

RDV	DAS ELISA	
KS 1232 - 0500R	Reagent Set	500 test
KS 1232 - 1000R	Reagent Set	1000 test
KS 1232 - 5000R	Reagent Set	5000 test
KS 1232 - P	Positive Control	
KS 1232 - N	Negative Control	

Rice Ragged Stunt Virus

RRSV	DAS ELISA	
KS 1358 - 0500R	Reagent Set	500 test
KS 1358 - 1000R	Reagent Set	1000 test
KS 1358 - 5000R	Reagent Set	5000 test
KS 1358 - P	Positive Control	
KS 1358 - N	Negative Control	

Rice Stripe Virus

RSV	DAS ELISA	
KS 1233 - 0500R	Reagent Set	500 test
KS 1233 - 1000R	Reagent Set	1000 test
KS 1233 - 5000R	Reagent Set	5000 test
KS 1233 - P	Positive Control	
KS 1233 - N	Negative Control	

Rice Stripe Virus

RSV	Indirect ELISA	
KS 1425 - 0500R	Reagent Set	500 test
KS 1425 - 1000R	Reagent Set	1000 test
KS 1425 - 5000R	Reagent Set	5000 test
KS 1425 - P	Positive Control	
KS 1425 - N	Negative Control	

Rice Yellow Mottle Virus

RYMV	DAS ELISA	
KS 1236 - 0500R	Reagent Set	500 test
KS 1236 - 1000R	Reagent Set	1000 test
KS 1236 - 5000R	Reagent Set	5000 test
KS 1236 - P	Positive Control	
KS 1236 - N	Negative Control	

Ryegrass Cryptic Virus

RCV	DAS ELISA	
KS 1359 - 0500R	Reagent Set	500 test
KS 1359 - 1000R	Reagent Set	1000 test
KS 1359 - 5000R	Reagent Set	5000 test
KS 1359 - P	Positive Control	
KS 1359 - N	Negative Control	

Ryegrass Mosaic Virus

RMV	DAS ELISA	
KS 1360 - 0500R	Reagent Set	500 test
KS 1360 - 1000R	Reagent Set	1000 test
KS 1360 - 5000R	Reagent Set	5000 test
KS 1360 - P	Positive Control	
KS 1360 - N	Negative Control	

Saguaro Cactus Virus

SCV	DAS ELISA	
KS 1392 - 0500R	Reagent Set	500 test
KS 1392 - 1000R	Reagent Set	1000 test
KS 1392 - 5000R	Reagent Set	5000 test
KS 1392 - P	Positive Control	
KS 1392 - N	Negative Control	

Scrophularia Mottle Virus

SMoV	DAS ELISA	
KS 1393 - 0500R	Reagent Set	500 test
KS 1393 - 1000R	Reagent Set	1000 test
KS 1393 - 5000R	Reagent Set	5000 test
KS 1393 - P	Positive Control	
KS 1393 - N	Negative Control	

Satsuma Dwarf Virus

SDV	DAS ELISA	
KS 1237 - 0500R	Reagent Set	500 test
KS 1237 - 1000R	Reagent Set	1000 test
KS 1237 - 5000R	Reagent Set	5000 test
KS 1237 - P	Positive Control	
KS 1237 - N	Negative Control	

Shallot Latent Virus

SLV	DAS ELISA	Poly/Poly
CAB 34200/0500	Coating Ab	500 test
CAB 34200/0501	Coating Ab	1000 test
CAB 34200/0502	Coating Ab	5000 test
ECP 34200/0500	Conjugated Ab, Peroxidase	500 test
ECP 34200/1000	Conjugated Ab, Peroxidase	1000 test
ECP 34200/5000	Conjugated Ab, Peroxidase	5000 test
SRP 34200/0500	Reagent Set, Peroxidase	500 test
SRP 34200/1000	Reagent Set, Peroxidase	1000 test
SRP 34200/5000	Reagent Set, Peroxidase	5000 test
LNC 34200	Negative Control	
LPC 34200	Positive Control	

Shallot Latent Virus

ShVX	TAS ELISA	
KS 1426 - 0500R	Reagent Set	500 test
KS 1426 - 1000R	Reagent Set	1000 test

식물병 Viruses

KS 1426 - 5000R	Reagent Set	5000 test
KS 1426 - P	Positive Control	
KS 1426 - N	Negative Control	

Shallot Virus X

ShvX	DAS ELISA	
KS 1239 - 0500R	Reagent Set	500 test
KS 1239 - 1000R	Reagent Set	1000 test
KS 1239 - 5000R	Reagent Set	5000 test
KS 1239 - P	Positive Control	
KS 1239 - N	Negative Control	

Shallot Yellow Stripe Virus

SYSV	DAS ELISA	
KS 1240 - 0500R	Reagent Set	500 test
KS 1240 - 1000R	Reagent Set	1000 test
KS 1240 - 5000R	Reagent Set	5000 test
KS 1240 - P	Positive Control	
KS 1240 - N	Negative Control	

Shallot Yellow Stripe Virus

SYSV	TAS ELISA	
KS 1427 - 0500R	Reagent Set	500 test
KS 1427 - 1000R	Reagent Set	1000 test
KS 1427 - 5000R	Reagent Set	5000 test
KS 1427 - P	Positive Control	
KS 1427 - N	Negative Control	

Shallot Yellow Stripe Virus

SYSV	Compound Direct ELISA	Poly/Mono
CAB 34600/0500	Coating Ab	500 test
CAB 34600/1000	Coating Ab	1000 test
CAB 34600/5000	Coating Ab	5000 test
ECA 34600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 34600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 34600/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 34600/0500	Reagent Set,Alk.Phos.	500 test
SRA 34600/1000	Reagent Set,Alk.Phos.	1000 test
SRA 34600/5000	Reagent Set,Alk.Phos.	5000 test
LNC 34600	Negative Control	

Sieg River Virus

SRV	DAS ELISA	
KS 1376 - 0500R	Reagent Set	500 test
KS 1376 - 1000R	Reagent Set	1000 test
KS 1376 - 5000R	Reagent Set	5000 test
KS 1376 - P	Positive Control	
KS 1376 - N	Negative Control	

Soil-borne Cereal Mosaic Virus

SBCMV	DAS ELISA	
KS 1361 - 0500R	Reagent Set	500 test
KS 1361 - 1000R	Reagent Set	1000 test
KS 1361 - 5000R	Reagent Set	5000 test
KS 1361 - P	Positive Control	
KS 1361 - N	Negative Control	

Soil-borne Wheat Mosaic Virus

SBWMV	DAS ELISA	
KS 1362 - 0500R	Reagent Set	500 test
KS 1362 - 1000R	Reagent Set	1000 test
KS 1362 - 5000R	Reagent Set	5000 test
KS 1362 - P	Positive Control	
KS 1362 - N	Negative Control	

Soil-borne Wheat Mosaic Virus

SBWMV	TAS ELISA	
KS 1241 - 0500R	Reagent Set	500 test
KS 1241 - 1000R	Reagent Set	1000 test
KS 1241 - 5000R	Reagent Set	5000 test
KS 1241 - P	Positive Control	
KS 1241 - N	Negative Control	

Soil-borne Wheat Mosaic Virus

SBWMV	Compound Direct ELISA	Poly/Mono
CAB 42001/0500	Coating Ab	500 test
CAB 42001/1000	Coating Ab	1000 test
CAB 42001/5000	Coating Ab	5000 test
ECA 42001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 42001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 42001/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 42001/0500	Reagent Set,Alk.Phos.	500 test
SRA 42001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 42001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 42001	Negative Control	
LPC 42001	Positive Control	

Sonchus Yellow Net Virus

SYNV	DAS ELISA	
KS 1377 - 0500R	Reagent Set	500 test
KS 1377 - 1000R	Reagent Set	1000 test
KS 1377 - 5000R	Reagent Set	5000 test
KS 1377 - P	Positive Control	
KS 1377 - N	Negative Control	

Sorghum Mosaic Virus

SrMV	DAS ELISA	
KS 1243 - 0500R	Reagent Set	500 test

KS 1243 - 1000R	Reagent Set	1000 test
KS 1243 - 5000R	Reagent Set	5000 test
KS 1243 - P	Positive Control	
KS 1243 - N	Negative Control	

Southern Bean Mosaic Virus

SBMV	DAS ELISA	Poly/Poly
CAB 24100/0500	Coating Ab	500 test
CAB 24100/1000	Coating Ab	1000 test
CAB 24100/5000	Coating Ab	5000 test
ECA 24100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 24100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 24100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 24100/0500	Reagent Set,Alk.Phos.	500 test
SRA 24100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 24100/5000	Reagent Set,Alk.Phos.	5000 test
LNC 24100	Negative Control	
LPC 24100	Positive Control	

Soybean Mosaic Virus

SMV	DAS ELISA	Poly/Poly
CAB 33300/0500	Coating Ab	500 test
CAB 33300/1000	Coating Ab	1000 test
CAB 33300/5000	Coating Ab	5000 test
ECA 33300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 33300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 33300/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 33300/0500	Conjugated Ab,Peroxidase	500 test
ECP 33300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 33300/5000	Conjugated Ab,Peroxidase	5000 test
PSA 33300/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 33300/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 33300/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 33300/0096	PathoScreen Kit,Peroxidase	96 test
PSP 33300/0288	PathoScreen Kit,Peroxidase	288 test
PSP 33300/0480	PathoScreen Kit,Peroxidase	480 test
SRA 33300/0500	Reagent Set,Alk.Phos.	500 test
SRA 33300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 33300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 33300/0500	Reagent Set,Peroxidase	500 test
SRP 33300/1000	Reagent Set,Peroxidase	1000 test
SRP 33300/5000	Reagent Set,Peroxidase	5000 test
LNC 33300	Negative Control	
LPC 33300	Positive Control	

Squash Leaf Curl Virus

SLCV	DAS ELISA	
KS 1363 - 0500R	Reagent Set	500 test
KS 1363 - 1000R	Reagent Set	1000 test
KS 1363 - 5000R	Reagent Set	5000 test

KS 1363 - P	Positive Control
KS 1363 - N	Negative Control

Squash Mosaic Virus

SqMV	DAS ELISA	Poly/Poly
CAB 26400/0500	Coating Ab	500 test
CAB 26400/1000	Coating Ab	1000 test
CAB 26400/5000	Coating Ab	5000 test
ECA 26400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 26400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 26400/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 26400/0500	Conjugated Ab,Peroxidase	500 test
ECP 26400/1000	Conjugated Ab,Peroxidase	1000 test
ECP 26400/5000	Conjugated Ab,Peroxidase	5000 test
PSA 26400/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 26400/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 26400/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 26400/0096	PathoScreen Kit, Peroxidase	96 test
PSP 26400/0288	PathoScreen Kit, Peroxidase	288 test
PSP 26400/0480	PathoScreen Kit, Peroxidase	480 test
SRA 26400/0500	Reagent Set,Alk.Phos.	500 test
SRA 26400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 26400/5000	Reagent Set,Alk.Phos.	5000 test
SRP 26400/0500	Reagent Set,Peroxidase	500 test
SRP 26400/1000	Reagent Set,Peroxidase	1000 test
SRP 26400/5000	Reagent Set,Peroxidase	5000 test
*ISK 26400/0005	ImmunoStrip Test Kit	5 strips
*ISK 26400/0025	ImmunoStrip Test Kit	25 strips
*ISK 26400/0250	ImmunoStrip Test Kit	250 strips
*STX 26400/0005	ImmunoStrip Test	5 strips
*STX 26400/0025	ImmunoStrip Test	25 strips
*STX 26400/0250	ImmunoStrip Test	250 strips
LNC 26400	Negative Control	
LPC 26400	Positive Control	

*** ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** SqMV ImmunoStrip Buffer**

SqMV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

Strawberry Latent Ringspot Virus

SLRSV	DAS ELISA	Poly/Poly
CAB 14000/0500	Coating Ab	500 test
CAB 14000/1000	Coating Ab	1000 test
CAB 14000/5000	Coating Ab	5000 test
ECA 14000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 14000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 14000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 14000/0500	Conjugated Ab,Peroxidase	500 test

ECP 14000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 14000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 14000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 14000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 14000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 14000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 14000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 14000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 14000/0500	Reagent Set,Alk.Phos.	500 test
SRA 14000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 14000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 14000/0500	Reagent Set,Peroxidase	500 test
SRP 14000/1000	Reagent Set,Peroxidase	1000 test
SRP 14000/5000	Reagent Set,Peroxidase	5000 test
LNC 14000	Negative Control	
LPC 14000	Positive Control	

Strawberry Latent Ringspot Virus - chestnut isolate

SLRSV-chestnut DAS ELISA

KS 1428 - 0500R	Reagent Set	500 test
KS 1428 - 1000R	Reagent Set	1000 test
KS 1428 - 5000R	Reagent Set	5000 test
KS 1428 - P	Positive Control	
KS 1428 - N	Negative Control	

Strawberry Mild Yellow Edge Virus

SMYEV DAS ELISA

KS 1248 - 0500R	Reagent Set	500 test
KS 1248 - 1000R	Reagent Set	1000 test
KS 1248 - 5000R	Reagent Set	5000 test
KS 1248 - P	Positive Control	
KS 1248 - N	Negative Control	

Strawberry Mild Yellow Edge Virus

SMYEV Compound Direct ELISA Poly/Mono

CAB 14500/0500	Coating Ab	500 test
CAB 14500/1000	Coating Ab	1000 test
CAB 14500/5000	Coating Ab	5000 test
ECA 14500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 14500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 14500/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 14500/0500	Reagent Set,Alk.Phos.	500 test
SRA 14500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 14500/5000	Reagent Set,Alk.Phos.	5000 test
LNC 14500	Negative Control	
LPC 14500	Positive Control	

Strawberry Mild Yellow Edge Potex Virus

SMYEPV DAS ELISA

KS 1394 - 0500R	Reagent Set	500 test
KS 1394 - 1000R	Reagent Set	1000 test
KS 1394 - 5000R	Reagent Set	5000 test
KS 1394 - P	Positive Control	
KS 1394 - N	Negative Control	

Streptocarpus Flower-Break Virus

SFBV DAS ELISA

KS 1364 - 0500R	Reagent Set	500 test
KS 1364 - 1000R	Reagent Set	1000 test
KS 1364 - 5000R	Reagent Set	5000 test
KS 1364 - P	Positive Control	
KS 1364 - N	Negative Control	

Sugarcane Bacilliform Virus

SCBV DAS ELISA Poly/Poly

CAB 72200/0500	Coating Ab	500 test
CAB 72200/1000	Coating Ab	1000 test
CAB 72200/5000	Coating Ab	5000 test
ECA 72200/500	Conjugated Ab,Alk.Phos.	500 test
ECA 72200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 72200/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 72200/0500	Reagent Set,Alk.Phos.	500 test
SRA 72200/1000	Reagent Set,Alk.Phos.	1000 test
SRA 72200/5000	Reagent Set,Alk.Phos.	5000 test
LNC 72200	Negative Control	
LPC 72200	Positive Control	

Sugarcane Mosaic Virus

SCMV DAS ELISA Poly/Poly

CAB 18100/0500	Coating Ab	500 test
CAB 18100/1000	Coating Ab	1000 test
CAB 18100/5000	Coating Ab	5000 test
ECA 18100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 18100/0500	Reagent Set,Alk.Phos.	500 test
SRA 18100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18100/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18100	Negative Control	
LPC 18100	Positive Control	

Sweetpotato NCM- ELISA Kits

KS 1369 - 0200R	Reagent Set	200 samples for 8 viruses
-----------------	-------------	---------------------------

- 1) Sweetpotato Chlorotic Fleck Virus(SPCFV)
- 2) Sweetpotato Chlorotic Stunt Virus(SPCSV)
- 3) Sweetpotato Feathery Mottle Virus(SPFMV)

- 4) Sweetpotato Latent Virus(SwPLV)
- 5) Sweetpotato Mild Mottle Virus(SPMMV)
- 6) Sweetpotato Mild Speckling Virus(SPMSV)
- 7) Sweetpotato Caulimovirus(SPCaLV)
- 8) C - 6 Virus

Tobacco Etch Virus

TEV	DAS ELISA	Poly/Poly
CAB 49500/0500	Coating Ab	500 test
CAB 49500/1000	Coating Ab	1000 test
CAB 49500/5000	Coating Ab	5000 test
ECA 49500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 49500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 49500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 49500/0500	Conjugated Ab,Peroxidase	500 test
ECP 49500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 49500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 49500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 49500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 49500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 49500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 49500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 49500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 49500/0500	Reagent Set,Alk.Phos.	500 test
SRA 49500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 49500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 49500/0500	Reagent Set,Peroxidase	500 test
SRP 49500/1000	Reagent Set,Peroxidase	1000 test
SRP 49500/5000	Reagent Set,Peroxidase	5000 test
LNC 49500	Negative Control	
LPC 49500	Positive Control	

Tobacco Mild Green Mosaic Virus

TMGMV	DAS ELISA	
KS 1365 - 0500R	Reagent Set	500 test
KS 1365 - 1000R	Reagent Set	1000 test
KS 1365 - 5000R	Reagent Set	5000 test
KS 1365 - P	Positive Control	
KS 1365 - N	Negative Control	

Tobacco Mosaic Virus

TMV	DAS ELISA	Poly/Poly
CAB 57400/0500	Coating Ab	500 test
CAB 57400/1000	Coating Ab	1000 test
CAB 57400/5000	Coating Ab	5000 test
ECA 57400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 57400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 57400/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 57400/0500	Conjugated Ab,Peroxidase	500 test
ECP 57400/1000	Conjugated Ab,Peroxidase	1000 test

ECP 57400/5000	Conjugated Ab,Peroxidase	5000 test
PSA 57400/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 57400/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 57400/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 57400/0096	PathoScreen Kit, Peroxidase	96 test
PSP 57400/0288	PathoScreen Kit, Peroxidase	288 test
PSP 57400/0480	PathoScreen Kit, Peroxidase	480 test
SRA 57400/0500	Reagent Set,Alk.Phos.	500 test
SRA 57400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 57400/5000	Reagent Set,Alk.Phos.	5000 test
SRP 57400/0500	Reagent Set, Peroxidase	500 test
SRP 57400/1000	Reagent Set, Peroxidase	1000 test
SRP 57400/5000	Reagent Set, Peroxidase	5000 test
*ISK 57400/0005	ImmunoStrip Test Kit	5 strips
*ISK 57400/0025	ImmunoStrip Test Kit	25 strips
*ISK 57400/0250	ImmunoStrip Test Kit	250 strips
*STX 57400/0005	ImmunoStrip Test	5 strips
*STX 57400/0025	ImmunoStrip Test	25 strips
*STX 57400/0250	ImmunoStrip Test	250 strips
LNC 57400	Negative Control	
LPC 57400	Positive Control	

*** ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** TMV ImmunoStrip Buffer**

TMV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

CMV, INSV, TSWV, TMV ImmunoStrip Comb

*ISK 37800/0024	ImmunoStrip Comb Kit, 24 combs and 25 buffer filled mesh bags	
*STX 37800/0024	ImmunoStrip Comb	24 combs

* ImmunoStrip Comb test system for CMV, INSV, TSWV, TMV are packaged in sets of 24.

ImmunoStrip combs must be used with SEB1 buffer.

Please see accessories for a listing of SEB1 buffer products. ImmunoStrip kits include sample extract bags.

Tobacco Mosaic Virus

TMV	Indirect ELISA	
KS 1255 - 0500R	Reagent Set	500 test
KS 1255 - 1000R	Reagent Set	1000 test
KS 1255 - 5000R	Reagent Set	5000 test
KS 1255 - P	Positive Control	
KS 1255 - N	Negative Control	

Tobacco Mosaic Virus - c

TMV-c	DAS ELISA	Poly/Poly
CAB 57500/0500	Coating Ab	500 test
CAB 57500/1000	Coating Ab	1000 test

식물병 Viruses

CAB 57500/5000	Coating Ab	5000 test
ECA 57500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 57500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 57500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 57500/0500	Conjugated Ab,Peroxidase	500 test
ECP 57500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 57500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 57500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 57500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 57500/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 57500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 57500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 57500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 57500/0500	Reagent Set,Alk.Phos.	500 test
SRA 57500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 57500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 57500/0500	Reagent Set,Peroxidase	500 test
SRP 57500/1000	Reagent Set,Peroxidase	1000 test
SRP 57500/5000	Reagent Set,Peroxidase	5000 test
LNC 57500	Negative Control	
LPC 57500	Positive Control	

Tobacco Mosaic Virus - L

TMV-L DAS ELISA

KS 1257 - 0500R	Reagent Set	500 test
KS 1257 - 1000R	Reagent Set	1000 test
KS 1257 - 5000R	Reagent Set	5000 test
KS 1257 - P	Positive Control	
KS 1257 - N	Negative Control	

Tobacco Mosaic Virus - OM

TMV-OM DAS ELISA

KS 1258 - 0500R	Reagent Set	500 test
KS 1258 - 1000R	Reagent Set	1000 test
KS 1258 - 5000R	Reagent Set	5000 test
KS 1258 - P	Positive Control	
KS 1258 - N	Negative Control	

Tobacco Mosaic Virus - P

TMV-P DAS ELISA

KS 1259 - 0500R	Reagent Set	500 test
KS 1259 - 1000R	Reagent Set	1000 test
KS 1259 - 5000R	Reagent Set	5000 test
KS 1259 - P	Positive Control	
KS 1259 - N	Negative Control	

Tobacco Mosaic Virus - W

TMV-W DAS ELISA

KS 1260 - 0500R	Reagent Set	500 test
KS 1260 - 1000R	Reagent Set	1000 test

KS 1260 - 5000R	Reagent Set	5000 test
KS 1260 - P	Positive Control	
KS 1260 - N	Negative Control	

Tobacco Necrosis Virus

TNV DAS ELISA

KS 1261 - 0500R	Reagent Set	500 test
KS 1261 - 1000R	Reagent Set	1000 test
KS 1261 - 5000R	Reagent Set	5000 test
KS 1261 - P	Positive Control	
KS 1261 - N	Negative Control	

Tobacco Necrosis Virus-A type tomato isolate

TNV-A DAS ELISA

KS 1262 - 0500R	Reagent Set	500 test
KS 1262 - 1000R	Reagent Set	1000 test
KS 1262 - 5000R	Reagent Set	5000 test
KS 1262 - P	Positive Control	
KS 1262 - N	Negative Control	

Tobacco Necrosis Virus - cherry isolate

TNV-cherry DAS ELISA

KS 1263 - 0500R	Reagent Set	500 test
KS 1263 - 1000R	Reagent Set	1000 test
KS 1263 - 5000R	Reagent Set	5000 test
KS 1263 - P	Positive Control	
KS 1263 - N	Negative Control	

Tobacco Necrosis Virus D

TNV-D DAS ELISA

KS 1366 - 0500R	Reagent Set	500 test
KS 1366 - 1000R	Reagent Set	1000 test
KS 1366 - 5000R	Reagent Set	5000 test
KS 1366 - P	Positive Control	
KS 1366 - N	Negative Control	

Tobacco Necrosis Virus E

TNV-E DAS ELISA

KS 1264 - 0500R	Reagent Set	500 test
KS 1264 - 1000R	Reagent Set	1000 test
KS 1264 - 5000R	Reagent Set	5000 test
KS 1264 - P	Positive Control	
KS 1264 - N	Negative Control	

Tobacco Necrosis Virus Satellite

TNV-S DAS ELISA

KS 1367 - 0500R	Reagent Set	500 test
KS 1367 - 1000R	Reagent Set	1000 test
KS 1367 - 5000R	Reagent Set	5000 test
KS 1367 - P	Positive Control	

KS 1367 - N	Negative Control
-------------	------------------

Tobacco Rattle Virus

TRV	DAS ELISA	Poly/Poly
CAB 25000/0500	Coating Ab	500 test
CAB 25000/1000	Coating Ab	1000 test
CAB 25000/5000	Coating Ab	5000 test
ECA 25000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 25000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 25000/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 25000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 25000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 25000/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 25000/0500	Reagent Set,Alk.Phos.	500 test
SRA 25000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 25000/5000	Reagent Set,Alk.Phos.	5000 test
LNC 25000	Negative Control	
LPC 25000	Positive Control	

Tobacco Rattle Virus - Jeffersonia isolate

TRV-J	DAS ELISA	
KS 1266 - 0500R	Reagent Set	500 test
KS 1266 - 1000R	Reagent Set	1000 test
KS 1266 - 5000R	Reagent Set	5000 test
KS 1266 - P	Positive Control	
KS 1266 - N	Negative Control	

Tobacco Rattle Virus K - potato isolate

TRV-K	DAS ELISA	
KS 1267 - 0500R	Reagent Set	500 test
KS 1267 - 1000R	Reagent Set	1000 test
KS 1267 - 5000R	Reagent Set	5000 test
KS 1267 - P	Positive Control	
KS 1267 - N	Negative Control	

Tobacco Rattle Virus Mirow - potato isolate

TRV-M	DAS ELISA	
KS 1268 - 0500R	Reagent Set	500 test
KS 1268 - 1000R	Reagent Set	1000 test
KS 1268 - 5000R	Reagent Set	5000 test
KS 1268 - P	Positive Control	
KS 1268 - N	Negative Control	

Tobacco Ringspot Virus

TRSV	DAS ELISA	Poly/Poly
CAB 64000/0500	Coating Ab	500 test
CAB 64000/1000	Coating Ab	1000 test
CAB 64000/5000	Coating Ab	5000 test
ECA 64000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 64000/1000	Conjugated Ab,Alk.Phos.	1000 test

ECA 64000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 64000/0500	Conjugated Ab,Peroxidase	500 test
ECP 64000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 64000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 64000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 64000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 64000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 64000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 64000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 64000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 64000/0500	Reagent Set,Alk.Phos.	500 test
SRA 64000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 64000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 64000/0500	Reagent Set,Peroxidase	500 test
SRP 64000/1000	Reagent Set,Peroxidase	1000 test
SRP 64000/5000	Reagent Set,Peroxidase	5000 test
LNC 64000	Negative Control	
LPC 64000	Positive Control	

Tobacco Ringspot Virus -Calico strain

TRSV	DAS ELISA	
KS 1270 - 0500R	Reagent Set	500 test
KS 1270 - 1000R	Reagent Set	1000 test
KS 1270 - 5000R	Reagent Set	5000 test
KS 1270 - P	Positive Control	
KS 1270 - N	Negative Control	

Tobacco Streak Virus

TSV	DAS ELISA	
KS 1271 - 0500R	Reagent Set	500 test
KS 1271 - 1000R	Reagent Set	1000 test
KS 1271 - 5000R	Reagent Set	5000 test
KS 1271 - P	Positive Control	
KS 1271 - N	Negative Control	

Tobacco Streak Virus

TSV	Compound Direct ELISA	Poly/Mono
CAB 25500/0500	Coating Ab	500 test
CAB 25500/1000	Coating Ab	1000 test
CAB 25500/5000	Coating Ab	5000 test
ECA 25500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 25500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 25500/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 25500/0500	Conjugated Ab,Peroxidase	500 test
ECP 25500/1000	Conjugated Ab,Peroxidase	1000 test
ECP 25500/5000	Conjugated Ab,Peroxidase	5000 test
PSA 25500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 25500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 25500/0480	PathoScreen Kit,Alk.Phos.	480 test

식물병 Viruses

PSP 25500/0096	PathoScreen Kit,Peroxidase	96 test
PSP 25500/0288	PathoScreen Kit,Peroxidase	288 test
PSP 25500/0480	PathoScreen Kit,Peroxidase	480 test
SRA 25500/0500	Reagent Set,Alk.Phos.	500 test
SRA 25500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 25500/5000	Reagent Set,Alk.Phos.	5000 test
SRP 25500/0500	Reagent Set,Peroxidase	500 test
SRP 25500/1000	Reagent Set,Peroxidase	1000 test
SRP 25500/5000	Reagent Set,Peroxidase	5000 test
LNC 25500	Negative Control	
LPC 25500	Positive Control	

Tobacco Vein Banding Mosaic Virus

TVBMV	DAS ELISA	
KS 1395 - 0500R	Reagent Set	500 test
KS 1395 - 1000R	Reagent Set	1000 test
KS 1395 - 5000R	Reagent Set	5000 test
KS 1395 - P	Positive Control	
KS 1395 - N	Negative Control	

Tobacco Vein Mottling Virus

TVMV	DAS ELISA	Poly/Poly
CAB 46300/0500	Coating Ab	500 test
CAB 46300/1000	Coating Ab	1000 test
CAB 46300/5000	Coating Ab	5000 test
ECA 46300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 46300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 46300/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 46300/0500	Reagent Set,Alk.Phos.	500 test
SRA 46300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 46300/5000	Reagent Set,Alk.Phos.	5000 test
LNC 46300	Negative Control	
LPC 46300	Positive Control	

Tomato Aspermy Virus

TAV	DAS ELISA	Poly/Poly
CAB 37000/0500	Coating Ab	500 test
CAB 37000/1000	Coating Ab	1000 test
CAB 37000/5000	Coating Ab	5000 test
ECA 37000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 37000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 37000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 37000/0500	Conjugated Ab,Peroxidase	500 test
ECP 37000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 37000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 37000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 37000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 37000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 37000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 37000/0288	PathoScreen Kit,Peroxidase	288 test

PSP 37000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 37000/0500	Reagent Set,Alk.Phos.	500 test
SRA 37000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 37000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 37000/0500	Reagent Set,Peroxidase	500 test
SRP 37000/1000	Reagent Set,Peroxidase	1000 test
SRP 37000/5000	Reagent Set,Peroxidase	5000 test
LNC 37000	Negative Control	
LPC 37000	Positive Control	

Tomato Black Ring Virus

TBRV	DAS ELISA	
KS 1276 - 0500R	Reagent Set	500 test
KS 1276 - 1000R	Reagent Set	1000 test
KS 1276 - 5000R	Reagent Set	5000 test
KS 1276 - P	Positive Control	
KS 1276 - N	Negative Control	

Tomato Black Ring Virus - g strain

TBRV-g	DAS ELISA	
KS 1277 - 0500R	Reagent Set	500 test
KS 1277 - 1000R	Reagent Set	1000 test
KS 1277 - 5000R	Reagent Set	5000 test
KS 1277 - P	Positive Control	
KS 1277 - N	Negative Control	

Tomato Black Ring Virus - potato isolate

TBRV-potato	DAS ELISA	
KS 1278 - 0500R	Reagent Set	500 test
KS 1278 - 1000R	Reagent Set	1000 test
KS 1278 - 5000R	Reagent Set	5000 test
KS 1278 - P	Positive Control	
KS 1278 - N	Negative Control	

Tomato Black Ring Virus - S, G

TBRV-S,G	DAS ELISA	
KS 1368 - 0500R	Reagent Set	500 test
KS 1368 - 1000R	Reagent Set	1000 test
KS 1368 - 5000R	Reagent Set	5000 test
KS 1368 - P	Positive Control	
KS 1368 - N	Negative Control	

Tomato Bushy Stunt Virus

TBSV	DAS ELISA	Poly/Poly
CAB 45400/0500	Coating Ab	500 test
CAB 45400/1000	Coating Ab	1000 test
CAB 45400/5000	Coating Ab	5000 test
ECA 45400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 45400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 45400/5000	Conjugated Ab,Alk.Phos.	5000 test

PSA 45400/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 45400/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 45400/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 45400/0500	Reagent Set,Alk.Phos.	500 test
SRA 45400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 45400/5000	Reagent Set,Alk.Phos.	5000 test
LNC 45400	Negative Control	
LPC 45400	Positive Control	

Tomato Chlorotic Spot Virus

TCSV	DAS ELISA	
KS 1280 - 0500R	Reagent Set	500 test
KS 1280 - 1000R	Reagent Set	1000 test
KS 1280 - 5000R	Reagent Set	5000 test
KS 1280 - P	Positive Control	
KS 1280 - N	Negative Control	

Tomato Chlorotic Spot Virus (See Groundnut Ringspot Virus)

Tomato Mosaic Virus

ToMV	DAS ELISA	Poly/Poly
CAB 35400/0500	Coating Ab	500 test
CAB 35400/1000	Coating Ab	1000 test
CAB 35400/5000	Coating Ab	5000 test
ECA 35400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 35400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 35400/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 35400/0500	Conjugated Ab,Peroxidase	500 test
ECP 35400/1000	Conjugated Ab,Peroxidase	1000 test
ECP 35400/5000	Conjugated Ab,Peroxidase	5000 test
PSA 35400/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 35400/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 35400/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 35400/0096	PathoScreen Kit,Peroxidase	96 test
PSP 35400/0288	PathoScreen Kit,Peroxidase	288 test
PSP 35400/0480	PathoScreen Kit,Peroxidase	480 test
SRA 35400/0500	Reagent Set,Alk.Phos.	500 test
SRA 35400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 35400/5000	Reagent Set,Alk.Phos.	5000 test
SRP 35400/0500	Reagent Set,Peroxidase	500 test
SRP 35400/1000	Reagent Set,Peroxidase	1000 test
SRP 35400/5000	Reagent Set,Peroxidase	5000 test
LNC 35400	Negative Control	
LPC 35400	Positive Control	

Tomato Mosaic Virus - Dahlemense strain

ToMV	DAS ELISA	
KS 1429 - 0500R	Reagent Set	500 test
KS 1429 - 1000R	Reagent Set	1000 test

KS 1429 - 5000R	Reagent Set	5000 test
KS 1429 - P	Positive Control	
KS 1429 - N	Negative Control	

Tomato Mottle Virus (See Tomato Yellow Leaf Curl Virus)

Tomato Ringspot Virus

ToRSV	DAS ELISA	Poly/Poly
CAB 22000/0500	Coating Ab	500 test
CAB 22000/1000	Coating Ab	1000 test
CAB 22000/5000	Coating Ab	5000 test
ECA 22000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 22000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 22000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 22000/0500	Conjugated Ab,Peroxidase	500 test
ECP 22000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 22000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 22000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 22000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 22000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 22000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 22000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 22000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 22000/0500	Reagent Set,Alk.Phos.	500 test
SRA 22000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 22000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 22000/0500	Reagent Set,Peroxidase	500 test
SRP 22000/1000	Reagent Set,Peroxidase	1000 test
SRP 22000/5000	Reagent Set,Peroxidase	5000 test
LNC 22000	Negative Control	
LPC 22000	Positive Control	

Tomato Ringspot Virus - Chickadee strain

ToRSV-Ch	DAS ELISA	
KS 1285 - 0500R	Reagent Set	500 test
KS 1285 - 1000R	Reagent Set	1000 test
KS 1285 - 5000R	Reagent Set	5000 test
KS 1285 - P	Positive Control	
KS 1285 - N	Negative Control	

Tomato Ringspot Virus- Peach yellow bud mosaic strain

ToRSV-PYBM	DAS ELISA	
KS 1286 - 0500R	Reagent Set	500 test
KS 1286 - 1000R	Reagent Set	1000 test
KS 1286 - 5000R	Reagent Set	5000 test
KS 1286 - P	Positive Control	
KS 1286 - N	Negative Control	

Tomato Ringspot Virus - pelargonium isolate

ToRSV	DAS ELISA	
KS 1430 - 0500R	Reagent Set	500 test
KS 1430 - 1000R	Reagent Set	1000 test
KS 1430 - 5000R	Reagent Set	5000 test
KS 1430 - P	Positive Control	
KS 1430 - N	Negative Control	

Tomato Spotted Wilt Virus

TSWV	DAS ELISA	Mono/Mono
CAB 39300/0500	Coating Ab	500 test
CAB 39300/1000	Coating Ab	1000 test
CAB 39300/5000	Coating Ab	5000 test
ECA 39300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 39300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 39300/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 39300/0500	Conjugated Ab,Peroxidase	500 test
ECP 39300/1000	Conjugated Ab,Peroxidase	1000 test
ECP 39300/5000	Conjugated Ab,Peroxidase	5000 test
PSA 39300/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 39300/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 39300/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 39300/0096	PathoScreen Kit,Alk.Peroxidas	96 test
PSP 39300/0288	PathoScreen Kit,Alk.Peroxidas	288 test
PSP 39300/0480	PathoScreen Kit,Alk.Peroxidas	480 test
SRA 39300/0500	Reagent Set,Alk.Phos.	500 test
SRA 39300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 39300/5000	Reagent Set,Alk.Phos.	5000 test
SRP 39300/0500	Reagent Set,Peroxidase	500 test
SRP 39300/1000	Reagent Set,Peroxidase	1000 test
SRP 39300/5000	Reagent Set,Peroxidase	5000 test
*ISK 39300/0005	ImmunoStrip Test Kit	5 strips
*ISK 39300/0025	ImmunoStrip Test Kit	250 strips
*ISK 39300/0250	ImmunoStrip Test Kit	25 strips
*STX 39300/0005	ImmunoStrip Test	5 strips
*STX 39300/0025	ImmunoStrip Test	25 strips
*STX 39300/0250	ImmunoStrip Test	250 strips
LNC 39300	Negative Control	
LPC 39300	Positive Control	

*** ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** TSWV ImmunoStrip Buffer**

TSWV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

CMV, INSV, TSWV, TMV ImmunoStrip Comb

*ISK 37800/0024	ImmunoStrip Comb Kit, 24 combs and 25 buffer filled mesh bags	
*STX 37800/0024	ImmunoStrip Comb	24 combs

* ImmunoStrip Comb test system for CMV, INSV, TSWV, TMV are packaged in sets of 24. ImmunoStrip combs must be used with SEB1 buffer. Please see accessories for a listing of SEB1 buffer products. ImmunoStrip kits include sample extract bags.

Tomato Spotted Wilt Virus and Impatiens Necrotic Spot Virus

TSWV and INSV	DAS ELISA	Mixture
CAB 30400/0500	Coating Ab	500 test
CAB 30400/1000	Coating Ab	1000 test
CAB 30400/5000	Coating Ab	5000 test
ECA 30400/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 30400/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 30400/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 30400/0500	Reagent Set,Alk.Phos.	500 test
SRA 30400/1000	Reagent Set,Alk.Phos.	1000 test
SRA 30400/5000	Reagent Set,Alk.Phos.	5000 test
*ISK 30400/0005	ImmunoStrip Test Kit	5 strips
*ISK 30400/0025	ImmunoStrip Test Kit	25 strips
*ISK 30400/0250	ImmunoStrip Test Kit	250 strips
*STX 30400/0005	ImmunoStrip Test	5 strips
*STX 30400/0025	ImmunoStrip Test	25 strips
*STX 30400/0250	ImmunoStrip Test	250 strips

*** ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** TSWV/INSV ImmunoStrip Buffer**

TSWV/INSV ImmunoStrips must be used with SEB1 buffer which is offered separately. Please see accessories for a listing of SEB1 buffer products.

Tomato Spotted Wilt Virus

TSWV	TAS ELISA	
KS 1290 - 0500R	Reagent Set	500 test
KS 1290 - 1000R	Reagent Set	1000 test
KS 1290 - 5000R	Reagent Set	5000 test
KS 1290 - P	Positive Control	
KS 1290 - N	Negative Control	

Tomato Spotted Wilt Virus NSs

TSWV NSs	Indirect ELISA	
CAB 39500/0500	Coating Ab	500 test
CAB 39500/1000	Coating Ab	1000 test
CAB 39500/5000	Coating Ab	5000 test
ECA 39500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 39500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 39500/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 39500/0500	Reagent Set,Alk.Phos.	500 test
SRA 39500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 39500/5000	Reagent Set,Alk.Phos.	5000 test

LNC 39500	Negative Control
LPC 39500	Positive Control

Tomato Yellow Leaf Curl Virus

TYLCV	DAS ELISA	
KS 1431 - 0500R	Reagent Set	500 test
KS 1431 - 1000R	Reagent Set	1000 test
KS 1431 - 5000R	Reagent Set	5000 test
KS 1431 - P	Positive Control	
KS 1431 - N	Negative Control	

Tomato Yellow Leaf Curl Virus

TYLCV	TAS ELISA	
KS 1291 - 0500R	Reagent Set	500 test
KS 1291 - 1000R	Reagent Set	1000 test
KS 1291 - 5000R	Reagent Set	5000 test
KS 1291 - P	Positive Control	
KS 1291 - N	Negative Control	

Gemini Virus Screen**(Indian, European, African strains of TYLCV)**

TYLCV	TAS ELISA	
KS 1396 - 0100R	Reagent Set	100 test
KS 1396 - P	Positive Control	
KS 1396 - N	Negative Control	

Tospo Broad Range artificial mixture of GRSV, INSV, TCSV, TSWV

Tospo	DAS ELISA	
KS 1288 - 0500R	Reagent Set	500 test
KS 1288 - 1000R	Reagent Set	1000 test
KS 1288 - 5000R	Reagent Set	5000 test
KS 1288 - P	Positive Control	
KS 1288 - N	Negative Control	

Tospo Virus Group (See Tomato Spotted Wilt Virus and Impatiens Necrotic Spot Virus)**Tulip Breaking Virus (See Lily Mottle Virus)****Tulip Crinkle Virus**

TCV	DAS ELISA	
KS 1370 - 0500R	Reagent Set	500 test
KS 1370 - 1000R	Reagent Set	1000 test
KS 1370 - 5000R	Reagent Set	5000 test
KS 1370 - P	Positive Control	
KS 1370 - N	Negative Control	

Tulip Virus X

TVX	DAS ELISA	
KS 1294 - 0500R	Reagent Set	500 test
KS 1294 - 1000R	Reagent Set	1000 test
KS 1294 - 5000R	Reagent Set	5000 test
KS 1294 - P	Positive Control	
KS 1294 - N	Negative Control	

Tulip Virus X

TVX	Indirect ELISA	
KS 1432 - 0500R	Reagent Set	500 test
KS 1432 - 1000R	Reagent Set	1000 test
KS 1432 - 5000R	Reagent Set	5000 test
KS 1432 - P	Positive Control	
KS 1432 - N	Negative Control	

Turnip Mosaic Virus

TuMV	DAS ELISA	
KS 1295 - 0500R	Reagent Set	500 test
KS 1295 - 1000R	Reagent Set	1000 test
KS 1295 - 5000R	Reagent Set	5000 test
KS 1295 - P	Positive Control	
KS 1295 - N	Negative Control	

Turnip Mosaic Virus

TuMV	Compound Direct ELISA	Poly/Mono
CAB 18700/0500	Coating Ab	500 test
CAB 18700/1000	Coating Ab	1000 test
CAB 18700/5000	Coating Ab	5000 test
ECA 18700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 18700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 18700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 18700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 18700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 18700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 18700/0500	Reagent Set,Alk.Phos.	500 test
SRA 18700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 18700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 18700	Negative Control	
LPC 18700	Positive Control	

Turnip Yellow Mosaic Virus

TYMV	DAS ELISA	
KS 1297 - 0500R	Reagent Set	500 test
KS 1297 - 1000R	Reagent Set	1000 test
KS 1297 - 5000R	Reagent Set	5000 test
KS 1297 - P	Positive Control	
KS 1297 - N	Negative Control	

Turnip Yellows Virus
(See Beet Mild Yellowing Virus)

Viola Mottle Virus

VMV	DAS ELISA	
KS 1371 - 0500R	Reagent Set	500 test
KS 1371 - 1000R	Reagent Set	1000 test
KS 1371 - 5000R	Reagent Set	5000 test
KS 1371 - P	Positive Control	
KS 1371 - N	Negative Control	

Watermelon Fruit Blotch Virus
Spot Check

KS 1433 - 0010R	Spot Check	10 tests
KS 1433 - 0100R	Spot Check	100 tests

Watermelon Mosaic Virus 1
(See Papaya Ringspot Virus)

Watermelon Mosaic Virus 2

WMV2	DAS ELISA	Poly/Poly
CAB 54000/0500	Coating Ab	500 test
CAB 54000/1000	Coating Ab	1000 test
CAB 54000/5000	Coating Ab	5000 test
ECA 54000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 54000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 54000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 54000/0500	Conjugated Ab,Peroxidase	500 test
ECP 54000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 54000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 54000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 54000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 54000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 54000/0096	PathoScreen Peroxidase	96 test
PSP 54000/0288	PathoScreen Peroxidase	288 test
PSP 54000/0480	PathoScreen Peroxidase	480 test
SRA 54000/0500	Reagent Set,Alk.Phos.	500 test
SRA 54000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 54000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 54000/0500	Reagent Set,Peroxidase	500 test
SRP 54000/1000	Reagent Set,Peroxidase	1000 test
SRP 54000/5000	Reagent Set,Peroxidase	5000 test
LNC 54000	Negative Control	
LPC 54000	Positive Control	

Watermelon Sliver Mottle Virus

WSMoV	DAS ELISA	
KS 1372 - 0500R	Reagent Set	500 test
KS 1372 - 1000R	Reagent Set	1000 test
KS 1372 - 5000R	Reagent Set	5000 test

KS 1372 - P	Positive Control
KS 1372 - N	Negative Control

Watermelon Silver Mottle Virus
and Groundnut Bud Necrosis Virus

WSMoV/GBNV	Compound Direct ELISA	Poly/Mono
CAB 61500/0500	Coating Ab	500 test
CAB 61500/1000	Coating Ab	1000 test
CAB 61500/5000	Coating Ab	5000 test
ECA 61500/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 61500/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 61500/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 61500/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 61500/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 61500/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 61500/0500	Reagent Set,Alk.Phos.	500 test
SRA 61500/1000	Reagent Set,Alk.Phos.	1000 test
SRA 61500/5000	Reagent Set,Alk.Phos.	5000 test
LNC 61500	Negative Control	

Wheat Dwarf Virus

WDV	DAS ELISA	
KS 1301 - 0500R	Reagent Set	500 test
KS 1301 - 1000R	Reagent Set	1000 test
KS 1301 - 5000R	Reagent Set	5000 test
KS 1301 - P	Positive Control	
KS 1301 - N	Negative Control	

Wheat Soil-Borne Mosaic Virus
(See Soil-Borne Wheat Mosaic Virus)

Wheat Spindle Streak Mosaic Virus

WSSMV	DAS ELISA	Poly/Poly
CAB 43001/0500	Coating Ab	500 test
CAB 43001/1000	Coating Ab	1000 test
CAB 43001/5000	Coating Ab	5000 test
ECA 43001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 43001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 43001/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 43001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 43001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 43001/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 43001/0500	Reagent Set,Alk.Phos.	500 test
SRA 43001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 43001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 43001	Negative Control	
LPC 43001	Positive Control	

Wheat Streak Mosaic Virus

WSMV	DAS ELISA	Poly/Poly
CAB 47001/0500	Coating Ab	500 test
CAB 47001/1000	Coating Ab	1000 test
CAB 47001/5000	Coating Ab	5000 test
ECA 47001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 47001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 47001/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 47001/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 47001/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 47001/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 47001/0500	Reagent Set,Alk.Phos.	500 test
SRA 47001/1000	Reagent Set,Alk.Phos.	1000 test
SRA 47001/5000	Reagent Set,Alk.Phos.	5000 test
LNC 47001	Negative Control	
LPC 47001	Positive Control	

Wheat Yellow Mosaic Virus

WYMV	DAS ELISA	
KS 1305 - 0500R	Reagent Set	500 test
KS 1305 - 1000R	Reagent Set	1000 test
KS 1305 - 5000R	Reagent Set	5000 test
KS 1305 - P	Positive Control	
KS 1305 - N	Negative Control	

Whitefly-Transmitted (Begomoviruses) Geminiviruses

WTG	*Nucleic Acid Hybridization	
NAK 15000/0096	Nucleic acid assay, 1 membran	96 assays

White Clover Cryptic Virus 1

WCCV 1	DAS ELISA	
KS 1373 - 0500R	Reagent Set	500 test
KS 1373 - 1000R	Reagent Set	1000 test
KS 1373 - 5000R	Reagent Set	5000 test
KS 1373 - P	Positive Control	
KS 1373 - N	Negative Control	

White Lupine Mosaic Virus

WLMV	DAS ELISA	
KS 1306 - 0500R	Reagent Set	500 test
KS 1306 - 1000R	Reagent Set	1000 test
KS 1306 - 5000R	Reagent Set	5000 test
KS 1306 - P	Positive Control	
KS 1306 - N	Negative Control	

Wild Potato Mosaic Virus

WPMV	DAS ELISA	
KS 1374 - 0500R	Reagent Set	500 test

KS 1374 - 1000R	Reagent Set	1000 test
KS 1374 - 5000R	Reagent Set	5000 test
KS 1374 - P	Positive Control	
KS 1374 - N	Negative Control	

Yam Mosaic Virus

YMV	DAS ELISA	
KS 1375 - 0500R	Reagent Set	500 test
KS 1375 - 1000R	Reagent Set	1000 test
KS 1375 - 5000R	Reagent Set	5000 test
KS 1375 - P	Positive Control	
KS 1375 - N	Negative Control	

Zucchini Green Mottle Mosaic Virus

ZGMMV	DAS ELISA	
KS 1319 - 0500R	Reagent Set	500 test
KS 1319 - 1000R	Reagent Set	1000 test
KS 1319 - 5000R	Reagent Set	5000 test
KS 1319 - P	Positive Control	
KS 1319 - N	Negative Control	

Zucchini Yellow Fleck Virus

ZYFV	DAS ELISA	
KS 1307 - 0500R	Reagent Set	500 test
KS 1307 - 1000R	Reagent Set	1000 test
KS 1307 - 5000R	Reagent Set	5000 test
KS 1307 - P	Positive Control	
KS 1307 - N	Negative Control	

Zucchini Yellow Mosaic Virus

ZYMV	DAS ELISA	
KS 1308 - 0500R	Reagent Set	500 test
KS 1308 - 1000R	Reagent Set	1000 test
KS 1308 - 5000R	Reagent Set	5000 test
KS 1308 - P	Positive Control	
KS 1308 - N	Negative Control	

Zucchini Yellow Mosaic Virus

ZYMV	Indirect ELISA	-- /Mono
CAB 77700/0500	Coating Ab	500 test
CAB 77700/1000	Coating Ab	1000 test
CAB 77700/5000	Coating Ab	5000 test
ECA 77700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 77700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 77700/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 77700/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 77700/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 77700/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 77700/0500	Reagent Set,Alk.Phos.	500 test
SRA 77700/1000	Reagent Set,Alk.Phos.	1000 test

SRA 77700/5000	Reagent Set,Alk.Phos.	5000 test
LNC 77700	Negative Control	
LPC 77700	Positive Control	

Chrysanthemum Chlorotic Mottle Viroid

CChMVd *Nucleic acid hybridization

NAK 07000/0096	Nucleic acid assay	1membrane
----------------	--------------------	-----------

Chrysanthemum Stunt Viroid

CSVd *Nucleic acid hybridization

NAK 09000/0096	Nucleic acid assay	1membrane
----------------	--------------------	-----------

Potato Spindle Tuber Viroid

PSTVd *Nucleic acid hybridization

DLP 08000/0001	Probe, digoxigenin label. Add to 8ml prehybridization buffer	
----------------	--	--

NAK 08000/0096	Nucleic acid assay	1membrane
----------------	--------------------	-----------

*** Agdia nucleic acid hybridization kits include detailed protocols and all the reagents necessary to extract and enrich RNA from plant samples.**

Simply spot samples on the provided membrane and return it to Agdia. We complete the nucleic acid hybridization assay at Agdia lab and send you a report on your test results. Membrane processing is included in the assay price.

Acidovorax avenae subsp. citrulli

Aac	DAS ELISA	Poly/Poly
CAB 14800/0500	Coating Ab	500 test
CAB 14800/1000	Coating Ab	1000 test
CAB 14800/5000	Coating Ab	5000 test
ECA 14800/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 14800/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 14800/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 14800/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 14800/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 14800/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 14800/0500	Reagent Set,Alk.Phos.	500 test
SRA 14800/1000	Reagent Set,Alk.Phos.	1000 test
SRA 14800/5000	Reagent Set,Alk.Phos.	5000 test
LPC 14800	Positive Control	
LNC 14800	Negative Control	

Agrobacterium radiobacter subsp. tumefaciens

DAS ELISA		
KS 2001 - 0500R	Reagent Set	500 test
KS 2001 - 1000R	Reagent Set	1000 test
KS 2001 - 5000R	Reagent Set	5000 test
KS 2001 - P	Positive Control	
KS 2001 - N	Negative Control	

Agrobacterium tumefaciens (strain-specific)

At	DAS ELISA	
KS 2002 - 0500R	Reagent Set	500 test
KS 2002 - 1000R	Reagent Set	1000 test
KS 2002 - 5000R	Reagent Set	5000 test
KS 2002 - P	Positive Control	
KS 2002 - N	Negative Control	

Banana Blood Disease

(See *Ralstonia solanacearum*)

Burkholderia gladioli pv. allicola

Bga	Indirect ELISA/IF	
KS 2003 - 0500R	Reagent Set	500 test
KS 2003 - 1000R	Reagent Set	1000 test
KS 2003 - 5000R	Reagent Set	5000 test
KS 2003 - P	Positive Control	
KS 2003 - N	Negative Control	

Burkholderia glumae

Bg	DAS ELISA	
KS 2504 - 0500R	Reagent Set	500 test
KS 2504 - 1000R	Reagent Set	1000 test
KS 2504 - 5000R	Reagent Set	5000 test

KS 2504 - P	Positive Control
KS 2504 - N	Negative Control

Burkholderia cepacia

Bc	DAS ELISA	
KS 2505 - 0500R	Reagent Set	500 test
KS 2505 - 1000R	Reagent Set	1000 test
KS 2505 - 5000R	Reagent Set	5000 test
KS 2505 - P	Positive Control	
KS 2505 - N	Negative Control	

Cladosporium spp.

	Indirect ELISA	
KS 2104 - 0500R	Reagent Set	500 test
KS 2104 - 1000R	Reagent Set	1000 test
KS 2104 - 5000R	Reagent Set	5000 test
KS 2104 - P	Positive Control	
KS 2104 - N	Negative Control	

Clavibacter michiganensis subsp. insidiosum

Cmi	DAS ELISA	
KS 2004 - 0500R	Reagent Set	500 test
KS 2004 - 1000R	Reagent Set	1000 test
KS 2004 - 5000R	Reagent Set	5000 test
KS 2004 - P	Positive Control	
KS 2004 - N	Negative Control	

Clavibacter michiganensis subsp. insidiosum

Cmi	Indirect ELISA/IF	
KS 2098 - 0500R	Reagent Set	500 test
KS 2098 - 1000R	Reagent Set	1000 test
KS 2098 - 5000R	Reagent Set	5000 test
KS 2098 - P	Positive Control	
KS 2098 - N	Negative Control	

Clavibacter michiganensis subsp. michiganensis

Cmm	DAS ELISA	Poly/Poly
CAB 44000/0500	Coating Ab	500 test
CAB 44000/1000	Coating Ab	1000 test
CAB 44000/5000	Coating Ab	5000 test
ECA 44000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 44000/0500	Conjugated Ab,Peroxidase	500 test
ECP 44000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 44000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 44000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 44000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 44000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 44000/0096	PathoScreen Kit,Peroxidase	96 test

PSP 44000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 44000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 44000/0500	Reagent Set,Alk.Phos.	500 test
SRA 44000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 44000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 44000/0500	Reagent Set,Peroxidase	500 test
SRP 44000/1000	Reagent Set,Peroxidase	1000 test
SRP 44000/5000	Reagent Set,Peroxidase	5000 test
LNC 44000	Negative Control	
LPC 44000	Positive Control	

Clavibacter michiganensis subsp. michiganensis

Cmm	Indirect ELISA	-- /Mono
CAB 44001/0500	Coating Ab	500 test
CAB 44001/1000	Coating Ab	1000 test
CAB 44001/5000	Coating Ab	5000 test
ECA 44001/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44001/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44001/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 44001/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 44001/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 44001/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
BID 44001/0096	Bacterial ID	96 test
BID 44001/0288	Bacterial ID	288 test
BID 44001/0480	Bacterial ID	480 test
LNC 44001	Negative Control	
LPC 44001	Positive Control	

Clavibacter michiganensis subsp. sepedonicus

Cms	IF	-- /Mono
IFA 70001/0125	Immunofluorescence Assay	125 tests
IFA 70001/0250	Immunofluorescence Assay	250 tests
IFA 70001/0500	Immunofluorescence Assay	500 tests
IFA 70001/1000	Immunofluorescence Assay	1000 tests
LNC 70001	Negative Control	
LPC 70001	Positive Control	

Clavibacter michiganensis subsp. sepedonicus

Cms	DAS ELISA	Mono/Poly
CAB 70002/0500	Coating Ab	500 test
CAB 70002/1000	Coating Ab	1000 test
CAB 70002/5000	Coating Ab	5000 test
ECA 70002/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 70002/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 70002/5000	Conjugated Ab, Alk.Phos.	5000 test
PSA 70002/0096	PathoScreen Kit, Alk.Phos.	96 test
PSA 70002/0288	PathoScreen Kit, Alk.Phos.	288 test
PSA 70002/0480	PathoScreen Kit, Alk.Phos.	480 test
SRA 70002/0500	Reagent Set, Alk.Phos.	500 test
SRA 70002/1000	Reagent Set, Alk.Phos.	1000 test

식물병 Bacteria

SRA 70002/5000	Reagent Set, Alk.Phos.	5000 test
LNC 70002	Negative Control	
LPC 70002	Positive Control	

Clavibacter michiganensis subsp. *sepedonicus*

Cms	IMF, Canadian Protocol	-- /Mono
CAI 70004	Mab Clone 9A1, BRR Protocol	

Clavibacter michiganensis subsp. *sepedonicus*

Cms	ELISA, Canadian Protocol	Poly/Mono
CAE 70005/0500	Mab Clone 1H3, BRR Protocol	500 test

Clavibacter michiganensis subsp. *sepedonicus*

Cms	Indirect ELISA/IF	
KS 2008 - 0500R	Reagent Set	500 test
KS 2008 - 1000R	Reagent Set	1000 test
KS 2008 - 5000R	Reagent Set	5000 test
KS 2008 - P	Positive Control	
KS 2008 - N	Negative Control	

Clavibacter michiganensis subsp. *sepedonicus*-Cy3

Cms-Cy3	IF-kit	
KS 2507-10000R	Reagent Set	10000 test
KS 2507 - P	Positive Control	
KS 2507 - N	Negative Control	

Clavibacter michiganensis subsp. *sepedonicus*-FITC

Cms-FITC	IF-kit	
KS 2508-10000R	Reagent Set	10000 test
KS 2508 - P	Positive Control	
KS 2508 - N	Negative Control	

Clavibacter michiganensis subsp. *tessellarius*

Cmt	DAS ELISA	Poly/Poly
CAB 44100/0500	Coating Ab	500 test
CAB 44100/1000	Coating Ab	1000 test
CAB 44100/5000	Coating Ab	5000 test
ECA 44100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 44100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 44100/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 44100/0500	Reagent Set,Alk.Phos.	500 test
SRA 44100/1000	Reagent Set,Alk.Phos.	1000 test
SRA 44100/5000	Reagent Set,Alk.Phos.	5000 test
LPC 44100	Positive Control	
LNC 44100	Negative Control	

Curtobacterium flaccumfaciens pv. *flaccumfaciens*

Cff	DAS ELISA	
KS 2006 - 0500R	Reagent Set	500 test
KS 2006 - 1000R	Reagent Set	1000 test
KS 2006 - 5000R	Reagent Set	5000 test
KS 2006 - P	Positive Control	
KS 2006 - N	Negative Control	

Curtobacterium flaccumfaciens pv. *flaccumfaciens*

Cff	Indirect ELISA/IF	
KS 2532 - 0500R	Reagent Set	500 test
KS 2532 - 1000R	Reagent Set	1000 test
KS 2532 - 5000R	Reagent Set	5000 test
KS 2532 - P	Positive Control	
KS 2532 - N	Negative Control	

Curtobacterium flaccumfaciens pv. *oortii*

Cfo	DAS ELISA	
KS 2007 - 0500R	Reagent Set	500 test
KS 2007 - 1000R	Reagent Set	1000 test
KS 2007 - 5000R	Reagent Set	5000 test
KS 2007 - P	Positive Control	
KS 2007 - N	Negative Control	

Curtobacterium flaccumfaciens pv. *oortii*

Cfo	Indirect ELISA/IF	
KS 2533 - 0500R	Reagent Set	500 test
KS 2533 - 1000R	Reagent Set	1000 test
KS 2533 - 5000R	Reagent Set	5000 test
KS 2533 - P	Positive Control	
KS 2533 - N	Negative Control	

Erwinia amylovora

Ea	DAS ELISA	
KS 2009 - 0500R	Reagent Set	500 test
KS 2009 - 1000R	Reagent Set	1000 test
KS 2009 - 5000R	Reagent Set	5000 test
KS 2009 - P	Positive Control	
KS 2009 - N	Negative Control	

Erwinia amylovora

Ea	Indirect ELISA	
KS 2011 - 0500R	Reagent Set	500 test
KS 2011 - 1000R	Reagent Set	1000 test
KS 2011 - 5000R	Reagent Set	5000 test
KS 2011 - P	Positive Control	
KS 2011 - N	Negative Control	

Erwinia carotovora

Ec	Indirect ELISA	
KS 2525 - 0500R	Reagent Set	500 test
KS 2525 - 1000R	Reagent Set	1000 test
KS 2525 - 5000R	Reagent Set	5000 test
KS 2525 - P	Positive Control	
KS 2525 - N	Negative Control	

Erwinia carotovora

Ec	Immunoblot	-- /Poly
PBK 90500/0096	Immunoblot Kit, Alk. Phos.	l membrane
LNC 90500	Negative control	
LPC 90500	Positive Control	

Erwinia carotovora subsp. atroseptica

Eca	DAS ELISA	Poly/Mono
CAB 80003/0500	Coating Ab	500 test
CAB 80003/1000	Coating Ab	1000 test
CAB 80003/5000	Coating Ab	5000 test
ECA 80003/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 80003/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 80003/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 80003/0500	Conjugated Ab,Peroxidase	500 test
ECP 80003/1000	Conjugated Ab,Peroxidase	1000 test
ECP 80003/5000	Conjugated Ab,Peroxidase	5000 test
PSA 80003/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 80003/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 80003/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 80003/0096	PathoScreen Kit,Peroxidase	96 test
PSP 80003/0288	PathoScreen Kit,Peroxidase	288 test
PSP 80003/0480	PathoScreen Kit,Peroxidase	480 test
SRA 80003/0500	Reagent Set,Alk.Phos.	500 test
SRA 80003/1000	Reagent Set,Alk.Phos.	1000 test
SRA 80003/5000	Reagent Set,Alk.Phos.	5000 test
SRP 80003/0500	Reagent Set,Peroxidase	500 test
SRP 80003/1000	Reagent Set,Peroxidase	1000 test
SRP 80003/5000	Reagent Set,Peroxidase	5000 test
LNC 80003	Negative Control	
LPC 80003	Positive Control	

Erwinia carotovora subsp. atroseptica

Eca	TAS ELISA	
KS 2510 - 0500R	Reagent Set	500 test
KS 2510 - 1000R	Reagent Set	1000 test
KS 2510 - 5000R	Reagent Set	5000 test
KS 2510 - P	Positive Control	
KS 2510 - N	Negative Control	

Erwinia carotovora subsp. atroseptica

Eca	Indirect ELISA/IF	
KS 2501 - 0500R	Reagent Set	500 test
KS 2501 - 1000R	Reagent Set	1000 test
KS 2501 - 5000R	Reagent Set	5000 test
KS 2501 - P	Positive Control	
KS 2501 - N	Negative Control	

Erwinia carotovora subsp. carotovora

Ecc	DAS ELISA	
KS 2013 - 0500R	Reagent Set	500 test
KS 2013 - 1000R	Reagent Set	1000 test
KS 2013 - 5000R	Reagent Set	5000 test
KS 2013 - P	Positive Control	
KS 2013 - N	Negative Control	

Erwinia chrysanthemi

Echr	DAS ELISA	
KS 2014 - 0500R	Reagent Set	500 test
KS 2014 - 1000R	Reagent Set	1000 test
KS 2014 - 5000R	Reagent Set	5000 test
KS 2014 - P	Positive Control	
KS 2014 - N	Negative Control	

Erwinia chrysanthemi

Echr	Indirect ELISA	-- /Mono
CAB 73300/0500	Coating Ab	500 test
CAB 73300/1000	Coating Ab	1000 test
CAB 73300/5000	Coating Ab	5000 test
ECA 73300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 73300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 73300/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 73300/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 73300/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 73300/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LNC 73300	Negative Control	
LPC 73300	Positive Control	

Erwinia chrysanthemi pv. dianthicola

DAS ELISA		
KS 2017 - 0500R	Reagent Set	500 test
KS 2017 - 1000R	Reagent Set	1000 test
KS 2017 - 5000R	Reagent Set	5000 test
KS 2017 - P	Positive Control	
KS 2017 - N	Negative Control	

Erwinia herbicola

Indirect ELISA/IF		
KS 2018 - 0500R	Reagent Set	500 test
KS 2018 - 1000R	Reagent Set	1000 test

식물병 Bacteria

KS 2018 - 5000R	Reagent Set	5000 test
KS 2018 - P	Positive Control	
KS 2018 - N	Negative Control	

Erwinia salicis

Es	DAS ELISA	
KS 2019 - 0500R	Reagent Set	500 test
KS 2019 - 1000R	Reagent Set	1000 test
KS 2019 - 5000R	Reagent Set	5000 test
KS 2019 - P	Positive Control	
KS 2019 - N	Negative Control	

Erwinia salicis

Es	Indirect ELISA/IF	
KS 2534 - 0500R	Reagent Set	500 test
KS 2534 - 1000R	Reagent Set	1000 test
KS 2534 - 5000R	Reagent Set	5000 test
KS 2534 - P	Positive Control	
KS 2534 - N	Negative Control	

Erwinia stewartii

Es	DAS ELISA	
KS 2021 - 0500R	Reagent Set	500 test
KS 2021 - 1000R	Reagent Set	1000 test
KS 2021 - 5000R	Reagent Set	5000 test
KS 2021 - P	Positive Control	
KS 2021 - N	Negative Control	

Erwinia stewartii

Es	Compound Direct ELISA	Poly/Mono
CAB 52000/0500	Coating Ab	500 test
CAB 52000/1000	Coating Ab	1000 test
CAB 52000/5000	Coating Ab	5000 test
ECA 52000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 52000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 52000/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 52000/0500	Conjugated Ab,Peroxidase	500 test
ECP 52000/1000	Conjugated Ab,Peroxidase	1000 test
ECP 52000/5000	Conjugated Ab,Peroxidase	5000 test
PSA 52000/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 52000/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 52000/0480	PathoScreen Kit,Alk.Phos.	480 test
PSP 52000/0096	PathoScreen Kit,Peroxidase	96 test
PSP 52000/0288	PathoScreen Kit,Peroxidase	288 test
PSP 52000/0480	PathoScreen Kit,Peroxidase	480 test
SRA 52000/0500	Reagent Set,Alk.Phos.	500 test
SRA 52000/1000	Reagent Set,Alk.Phos.	1000 test
SRA 52000/5000	Reagent Set,Alk.Phos.	5000 test
SRP 52000/0500	Reagent Set,Peroxidase	500 test
SRP 52000/1000	Reagent Set,Peroxidase	1000 test

SRP 52000/5000	Reagent Set,Peroxidase	5000 test
LNC 52000	Negative Control	
LPC 52000	Positive Control	

Erwinia stewartii

Es	Indirect ELISA/IF	
KS 2020 - 0500R	Reagent Set	500 test
KS 2020 - 1000R	Reagent Set	1000 test
KS 2020 - 5000R	Reagent Set	5000 test
KS 2020 - P	Positive Control	
KS 2020 - N	Negative Control	

Pantoea stewartii (See *Erwinia stewartii*)

Pantoea stewartii subsp. *stewartii*

Pss	DAS ELISA	
KS 2506 - 0500R	Reagent Set	500 test
KS 2506 - 1000R	Reagent Set	1000 test
KS 2506 - 5000R	Reagent Set	5000 test
KS 2506 - P	Positive Control	
KS 2506 - N	Negative Control	

Pectobacterium chrysanthemi

	DAS ELISA/IF	
KS 2538 - 0500R	Reagent Set	500 test
KS 2538 - 1000R	Reagent Set	1000 test
KS 2538 - 5000R	Reagent Set	5000 test
KS 2538 - P	Positive Control	
KS 2538 - N	Negative Control	

Pierce's Disease (See *Xylella fastidiosa*)

Pseudomonas avenae

Pa	DAS ELISA	
KS 2526 - 0500R	Reagent Set	500 test
KS 2526 - 1000R	Reagent Set	1000 test
KS 2526 - 5000R	Reagent Set	5000 test
KS 2526 - P	Positive Control	
KS 2526 - N	Negative Control	

Pseudomonas avenae

Pa	Indirect ELISA	-- /Mono
CAB 41900/0500	Coating Ab	500 test
CAB 41900/1000	Coating Ab	1000 test
CAB 41900/5000	Coating Ab	5000 test
ECA 41900/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 41900/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 41900/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 41900/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 41900/1000	Bacterial Reagent Set,Alk.Phos.	1000 test

BRA 41900/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 41900	Positive Control	
LNC 41900	Negative Control	

Pseudomonas andropogonis

Pa	DAS ELISA	
KS 2023 - 0500R	Reagent Set	500 test
KS 2023 - 1000R	Reagent Set	1000 test
KS 2023 - 5000R	Reagent Set	5000 test
KS 2023 - P	Positive Control	
KS 2023 - N	Negative Control	

Pseudomonas cichorii

Pci	DAS ELISA	
KS 2024 - 0500R	Reagent Set	500 test
KS 2024 - 1000R	Reagent Set	1000 test
KS 2024 - 5000R	Reagent Set	5000 test
KS 2024 - P	Positive Control	
KS 2024 - N	Negative Control	

Pseudomonas cichorii

Pci	Indirect ELISA/IF	
KS 2502 - 0500R	Reagent Set	500 test
KS 2502 - 1000R	Reagent Set	1000 test
KS 2502 - 5000R	Reagent Set	5000 test
KS 2502 - P	Positive Control	
KS 2502 - N	Negative Control	

Pseudomonas fuscovaginae

Pf	DAS ELISA	
KS 2527 - 0500R	Reagent Set	500 test
KS 2527 - 1000R	Reagent Set	1000 test
KS 2527 - 5000R	Reagent Set	5000 test
KS 2527 - P	Positive Control	
KS 2527 - N	Negative Control	

Pseudomonas fuscovaginae

Pf	Indirect ELISA	-- /Mono
CAB 53200/0500	Coating Ab	500 test
CAB 53200/1000	Coating Ab	1000 test
CAB 53200/5000	Coating Ab	5000 test
ECA 53200/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 53200/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 53200/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 53200/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 53200/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 53200/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 53200	Positive Control	
LNC 53200	Negative Control	

Pseudomonas gladioli pv. allicola

Pga	Express Kit	
KS 2512 - 0050R	Express Kit	50 test
KS 2512 - 0100R	Express Kit	100 test

Pseudomonas glumae

Pg	DAS ELISA	
KS 2528 - 0500R	Reagent Set	500 test
KS 2528 - 1000R	Reagent Set	1000 test
KS 2528 - 5000R	Reagent Set	5000 test
KS 2528 - P	Positive Control	
KS 2528 - N	Negative Control	

Pseudomonas glumae

Pg	Indirect ELISA	-- /Mono
CAB 62100/0500	Coating Ab	500 test
CAB 62100/1000	Coating Ab	1000 test
CAB 62100/5000	Coating Ab	5000 test
ECA 62100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 62100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 62100/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 62100/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 62100/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 62100/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 62100	Positive Control	
LNC 62100	Negative Control	

Pseudomonas marginalis pv. marginalis

Pm	DAS ELISA	
KS 2097 - 0500R	Reagent Set	500 test
KS 2097 - 1000R	Reagent Set	1000 test
KS 2097 - 5000R	Reagent Set	5000 test
KS 2097 - P	Positive Control	
KS 2097 - N	Negative Control	

Pseudomonas savastanoi pv. phaseolicola

	DAS ELISA	
KS 2513 - 0500R	Reagent Set	500 test
KS 2513 - 1000R	Reagent Set	1000 test
KS 2513 - 5000R	Reagent Set	5000 test
KS 2513 - P	Positive Control	
KS 2513 - N	Negative Control	

Pseudomonas solanacearum - Ps1
(See *Ralstonia solanacearum*)

Pseudomonas solanacearum
(See *Ralstonia solanacearum*)

Pseudomonas syringae pv. *avellanae*

Psa	Indirect ELISA/IF	
KS 2029 - 0500R	Reagent Set	500 test
KS 2029 - 1000R	Reagent Set	1000 test
KS 2029 - 5000R	Reagent Set	5000 test
KS 2029 - P	Positive Control	
KS 2029 - N	Negative Control	

Pseudomonas syringae pv. *glycinea*

Psg	Indirect ELISA/IF	
KS 2031 - 0500R	Reagent Set	500 test
KS 2031 - 1000R	Reagent Set	1000 test
KS 2031 - 5000R	Reagent Set	5000 test
KS 2031 - P	Positive Control	
KS 2031 - N	Negative Control	

Pseudomonas syringae pv. *lachrymans*

Psl	DAS ELISA	
KS 2034 - 0500R	Reagent Set	500 test
KS 2034 - 1000R	Reagent Set	1000 test
KS 2034 - 5000R	Reagent Set	5000 test
KS 2034 - P	Positive Control	
KS 2034 - N	Negative Control	

Pseudomonas syringae pv. *lachrymans*

Psl	Indirect ELISA/IF	
KS 2033 - 0500R	Reagent Set	500 test
KS 2033 - 1000R	Reagent Set	1000 test
KS 2033 - 5000R	Reagent Set	5000 test
KS 2033 - P	Positive Control	
KS 2033 - N	Negative Control	

Pseudomonas syringae pv. *mors-prunorum*

Psm	DAS ELISA	
KS 2035 - 0500R	Reagent Set	500 test
KS 2035 - 1000R	Reagent Set	1000 test
KS 2035 - 5000R	Reagent Set	5000 test
KS 2035 - P	Positive Control	
KS 2035 - N	Negative Control	

Pseudomonas syringae pv. *mors-prunorum*

Psm	Indirect ELISA/IF	
KS 2535 - 0500R	Reagent Set	500 test
KS 2535 - 1000R	Reagent Set	1000 test
KS 2535 - 5000R	Reagent Set	5000 test
KS 2535 - P	Positive Control	
KS 2535 - N	Negative Control	

Pseudomonas syringae pv. *phaseolicola*

Psph	DAS ELISA		Poly/Poly
CAB 53300/0500	Coating Ab		500 test
CAB 53300/1000	Coating Ab		1000 test
CAB 53300/5000	Coating Ab		5000 test
ECA 53300/0500	Conjugated Ab,Alk.Phos.		500 test
ECA 53300/1000	Conjugated Ab,Alk.Phos.		1000 test
ECA 53300/5000	Conjugated Ab,Alk.Phos.		5000 test
PSA 53300/0096	PathoScreen Kit,Alk.Phos.		96 test
PSA 53300/0288	PathoScreen Kit,Alk.Phos.		288 test
PSA 53300/0480	PathoScreen Kit,Alk.Phos.		480 test
SRA 53300/0500	Reagent Set,Alk.Phos.		500 test
SRA 53300/1000	Reagent Set,Alk.Phos.		1000 test
SRA 53300/5000	Reagent Set,Alk.Phos.		5000 test
LNC 53300	Negative Control		
LPC 53300	Positive Control		

Pseudomonas syringae pv. *phaseolicola*

Psph	Indirect ELISA/IF	
KS 2036 - 0500R	Reagent Set	500 test
KS 2036 - 1000R	Reagent Set	1000 test
KS 2036 - 5000R	Reagent Set	5000 test
KS 2036 - P	Positive Control	
KS 2036 - N	Negative Control	

Pseudomonas syringae pv. *pisi*

Pspi	DAS ELISA	
KS 2038 - 0500R	Reagent Set	500 test
KS 2038 - 1000R	Reagent Set	1000 test
KS 2038 - 5000R	Reagent Set	5000 test
KS 2038 - P	Positive Control	
KS 2038 - N	Negative Control	

Pseudomonas syringae pv. *pisi*

Pspi	Indirect ELISA/IF	
KS 2039 - 0500R	Reagent Set	500 test
KS 2039 - 1000R	Reagent Set	1000 test
KS 2039 - 5000R	Reagent Set	5000 test
KS 2039 - P	Positive Control	
KS 2039 - N	Negative Control	

Pseudomonas syringae pv. *porri*

Pspo	DAS ELISA	
KS 2040 - 0500R	Reagent Set	500 test
KS 2040 - 1000R	Reagent Set	1000 test
KS 2040 - 5000R	Reagent Set	5000 test
KS 2040 - P	Positive Control	
KS 2040 - N	Negative Control	

Pseudomonas syringae* pv. *porri

Pspo	Indirect ELISA/IF	
KS 2536 - 0500R	Reagent Set	500 test
KS 2536 - 1000R	Reagent Set	1000 test
KS 2536 - 5000R	Reagent Set	5000 test
KS 2536 - P	Positive Control	
KS 2536 - N	Negative Control	

Pseudomonas syringae* pv. *syringae

Pss	DAS ELISA	
KS 2042 - 0500R	Reagent Set	500 test
KS 2042 - 1000R	Reagent Set	1000 test
KS 2042 - 5000R	Reagent Set	5000 test
KS 2042 - P	Positive Control	
KS 2042 - N	Negative Control	

Pseudomonas syringae* pv. *syringae

Pss	Indirect ELISA/IF	
KS 2041 - 0500R	Reagent Set	500 test
KS 2041 - 1000R	Reagent Set	1000 test
KS 2041 - 5000R	Reagent Set	5000 test
KS 2041 - P	Positive Control	
KS 2041 - N	Negative Control	

Pseudomonas syringae* pv. *tomato

Pst	DAS ELISA	
KS 2037 - 0500R	Reagent Set	500 test
KS 2037 - 1000R	Reagent Set	1000 test
KS 2037 - 5000R	Reagent Set	5000 test
KS 2037 - P	Positive Control	
KS 2037 - N	Negative Control	

Pseudomonas syringae* pv. *tomato

Pst	Indirect ELISA/IF	
KS 2516 - 0500R	Reagent Set	500 test
KS 2516 - 1000R	Reagent Set	1000 test
KS 2516 - 5000R	Reagent Set	5000 test
KS 2516 - P	Positive Control	
KS 2516 - N	Negative Control	

Pseudomonas syringae* subsp. *savastanoi

Pssa	DAS ELISA	
KS 2044 - 0500R	Reagent Set	500 test
KS 2044 - 1000R	Reagent Set	1000 test
KS 2044 - 5000R	Reagent Set	5000 test
KS 2044 - P	Positive Control	
KS 2044 - N	Negative Control	

Pseudomonas syringae* subsp. *savastanoi

Pssa	Indirect ELISA/IF	
KS 2537 - 0500R	Reagent Set	500 test
KS 2537 - 1000R	Reagent Set	1000 test
KS 2537 - 5000R	Reagent Set	5000 test
KS 2537 - P	Positive Control	
KS 2537 - N	Negative Control	

Pseudomonas viridiflava

Pv	DAS ELISA	
KS 2529 - 0500R	Reagent Set	500 test
KS 2529 - 1000R	Reagent Set	1000 test
KS 2529 - 5000R	Reagent Set	5000 test
KS 2529 - P	Positive Control	
KS 2529 - N	Negative Control	

Pseudomonas viridiflava

Pv	Indirect ELISA/IF	
KS 2043 - 0500R	Reagent Set	500 test
KS 2043 - 1000R	Reagent Set	1000 test
KS 2043 - 5000R	Reagent Set	5000 test
KS 2043 - P	Positive Control	
KS 2043 - N	Negative Control	

Ralstonia solanacearum

Rsol	DAS ELISA	
KS 2096 - 0500R	Reagent Set	500 test
KS 2096 - 1000R	Reagent Set	1000 test
KS 2096 - 5000R	Reagent Set	5000 test
KS 2096 - P	Positive Control	
KS 2096 - N	Negative Control	

Ralstonia solanacearum

Rsol	DAS ELISA	
KS 2045 - 0500R	Reagent Set	500 test
KS 2045 - 1000R	Reagent Set	1000 test
KS 2045 - 5000R	Reagent Set	5000 test
KS 2045 - P	Positive Control	
KS 2045 - N	Negative Control	

Ralstonia solanacearum

Rs	ELISA and ImmunoStrip -- /Mono	
CAB 33900/0500	Coating Ab	500 test
CAB 33900/1000	Coating Ab	1000 test
CAB 33900/5000	Coating Ab	5000 test
ECP 33900/0500	Conjugated Ab, Peroxidase	500 test
ECP 33900/1000	Conjugated Ab, Peroxidase	1000 test
ECP 33900/5000	Conjugated Ab, Peroxidase	5000 test
PSP 33900/0096	PathoScreen Kit, Peroxidase	96 test
PSP 33900/0288	PathoScreen Kit, Peroxidase	288 test

식물병 Bacteria

PSP 33900/0480	PathoScreen Kit, Peroxidase	480 test
SRP 33900/0500	Reagent Set, Peroxidase	500 test
SRP 33900/1000	Reagent Set, Peroxidase	1000 test
SRP 33900/5000	Reagent Set, Peroxidase	5000 test
*ISK 33900/0005	ImmunoStrip Test Kit	5 strips
*ISK 33900/0025	ImmunoStrip Test Kit	25 strips
*ISK 33900/0250	ImmunoStrip Test Kit	250 strips
LNC 33900	Negative Control	
LPC 33900	Positive Control	

* ImmunoStrip Test kit

Including strips and buffer filled mesh bags

Ralstonia solanacearum

Rs	Indirect ELISA	
KS 2046 - 0500R	Reagent Set	500 test
KS 2046 - 1000R	Reagent Set	1000 test
KS 2046 - 5000R	Reagent Set	5000 test
KS 2046 - P	Positive Control	
KS 2046 - N	Negative Control	

Ralstonia solanacearum

Rs	TAS ELISA/IF	
KS 2047 - 0500R	Reagent Set	500 test
KS 2047 - 1000R	Reagent Set	1000 test
KS 2047 - 5000R	Reagent Set	5000 test
KS 2047 - P	Positive Control	
KS 2047 - N	Negative Control	

Ralstonia solanacearum-Cy3

Rs-Cy3	IF-kit	
KS 2517 - 10000R	Reagent Set	10000 test
KS 2517 - P	Positive Control	
KS 2517 - N	Negative Control	

Ralstonia solanacearum-FITC

Rs-FITC	IF-kit	
KS 2518 - 10000R	Reagent Set	10000 test
KS 2518 - P	Positive Control	
KS 2518 - N	Negative Control	

Rhodococcus fascians

Rhf	DAS ELISA	
KS 2092 - 0500R	Reagent Set	500 test
KS 2092 - 1000R	Reagent Set	1000 test
KS 2092 - 5000R	Reagent Set	5000 test
KS 2092 - P	Positive Control	
KS 2092 - N	Negative Control	

Rhodococcus fascians

Rhf	Indirect ELISA/IF	
KS 2509 - 0500R	Reagent Set	500 test
KS 2509 - 1000R	Reagent Set	1000 test
KS 2509 - 5000R	Reagent Set	5000 test
KS 2509 - P	Positive Control	
KS 2509 - N	Negative Control	

Spongospora subterranea f. sp. subterranea

Sss	DAS ELISA	
KS 2531 - 0500R	Reagent Set	500 test
KS 2531 - 1000R	Reagent Set	1000 test
KS 2531 - 5000R	Reagent Set	5000 test
KS 2531 - P	Positive Control	
KS 2531 - N	Negative Control	

Spiroplasma citri

Sc	DAS ELISA	oly/Poly
CAB 14700/0500	Coating Ab	500 test
CAB 14700/1000	Coating Ab	1000 test
CAB 14700/5000	Coating Ab	5000 test
ECA 14700/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 14700/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 14700/5000	Conjugated Ab,Alk.Phos.	5000 test
ECP 14700/0500	Conjugated Ab,Peroxidase	500 test
ECP 14700/1000	Conjugated Ab,Peroxidase	1000 test
ECP 14700/5000	Conjugated Ab,Peroxidase	5000 test
SRA 14700/0500	Reagent Set,Alk.Phos.	500 test
SRA 14700/1000	Reagent Set,Alk.Phos.	1000 test
SRA 14700/5000	Reagent Set,Alk.Phos.	5000 test
SRP 14700/0500	Reagent Set,Peroxidase	500 test
SRP 14700/1000	Reagent Set,Peroxidase	1000 test
SRP 14700/5000	Reagent Set,Peroxidase	5000 test
LNC 14700	Negative Control	
LPC 14700	Positive Control	

Xanthomonas

Xan	Indirect ELISA	-- /Mono
CAB 14600/0500	Coating Ab	500 test
CAB 14600/1000	Coating Ab	1000 test
CAB 14600/5000	Coating Ab	5000 test
ECA 14600/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 14600/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 14600/5000	Conjugated Ab, Alk.Phos.	5000 test
BRA 14600/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 14600/1000	Bacterial Reagent Set, Alk.Phos.	1000 test
BRA 14600/5000	Bacterial Reagent Set, Alk.Phos.	5000 test

Xanthomonas albilineans

Xalb	DAS ELISA	
KS 2530 - 0500R	Reagent Set	500 test
KS 2530 - 1000R	Reagent Set	1000 test
KS 2530 - 5000R	Reagent Set	5000 test
KS 2530 - P	Positive Control	
KS 2530 - N	Negative Control	

Xanthomonas albilineans

Xalb	Indirect ELISA	-- /Mono
CAB 71100/0500	Coating Ab	500 test
CAB 71100/1000	Coating Ab	1000 test
CAB 71100/5000	Coating Ab	5000 test
ECA 71100/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 71100/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 71100/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 71100/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 71100/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 71100/5000	Bacterial Reagent Set,Alk.Phos.	5000 test

Xanthomonas arboricola pv. pruni

	DAS ELISA	
KS 2523 - 0500R	Reagent Set	500 test
KS 2523 - 1000R	Reagent Set	1000 test
KS 2523 - 5000R	Reagent Set	5000 test
KS 2523 - P	Positive Control	
KS 2523 - N	Negative Control	

Xanthomonas axonopodis pv. arboricola

	DAS ELISA	
KS 2524 - 0500R	Reagent Set	500 test
KS 2524 - 1000R	Reagent Set	1000 test
KS 2524 - 5000R	Reagent Set	5000 test
KS 2524 - P	Positive Control	
KS 2524 - N	Negative Control	

Xanthomonas axonopodis pv. begoniae

Xab	DAS ELISA	
KS 2093 - 0500R	Reagent Set	500 test
KS 2093 - 1000R	Reagent Set	1000 test
KS 2093 - 5000R	Reagent Set	5000 test
KS 2093 - P	Positive Control	
KS 2093 - N	Negative Control	

***Xanthomonas axonopodis pv. citri*
(See *Xanthomonas campestris pv. citri*)*****Xanthomonas campestris***

Xc	Immunoblot	-- /Poly
PBK 69000/0096	Immunoblot kit, Alk.Phos.	l membrane

LNC 69000	Negative Control
LPC 69000	Positive Control

Xanthomonas campestris pv. armoraciae

Xcarm	Indirect ELISA	-- /Mono
CAB 88800/0500	Coating Ab	500 test
CAB 88800/1000	Coating Ab	1000 test
CAB 88800/5000	Coating Ab	5000 test
ECA 88800/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 88800/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 88800/5000	Conjugated Ab, Alk.Phos.	5000 test
BRA 88800/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 88800/1000	Bacterial Reagent Set, Alk.Phos.	1000 test
BRA 88800/5000	Bacterial Reagent Set, Alk.Phos.	5000 test
LPC 88800	Positive Control	

Xanthomonas campestris - lobelia isolate

Xcl	DAS ELISA	
KS 2065 - 0500R	Reagent Set	500 test
KS 2065 - 1000R	Reagent Set	1000 test
KS 2065 - 5000R	Reagent Set	5000 test
KS 2065 - P	Positive Control	
KS 2065 - N	Negative Control	

Xanthomonas campestris - lobelia isolate

Xcl	Indirect ELISA/IF	
KS 2519 - 0500R	Reagent Set	500 test
KS 2519 - 1000R	Reagent Set	1000 test
KS 2519 - 5000R	Reagent Set	5000 test
KS 2519 - P	Positive Control	
KS 2519 - N	Negative Control	

Xanthomonas campestris pv. armoraciae

Xcarm	DAS ELISA	
KS 2520 - 0500R	Reagent Set	500 test
KS 2520 - 1000R	Reagent Set	1000 test
KS 2520 - 5000R	Reagent Set	5000 test
KS 2520 - P	Positive Control	
KS 2520 - N	Negative Control	

Xanthomonas campestris pv. armoraciae

Xcarm	Indirect ELISA	-- /Mono
CAB 88800/0500	Coating Ab	500 test
CAB 88800/1000	Coating Ab	1000 test
CAB 88800/5000	Coating Ab	5000 test
ECA 88800/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 88800/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 88800/5000	Conjugated Ab, Alk.Phos.	5000 test
BRA 88800/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 88800/1000	Bacterial Reagent Set, Alk.Phos.	1000 test

식물병 Bacteria

BRA 88800/5000	Bacterial Reagent Set, Alk.Phos.	5000 test
LNC 88800	Negative Control	
LPC 88800	Positive Control	

Xanthomonas campestris pv. begoniae

Xcb	DAS ELISA	
KS 2054 - 0500R	Reagent Set	500 test
KS 2054 - 1000R	Reagent Set	1000 test
KS 2054 - 5000R	Reagent Set	5000 test
KS 2054 - P	Positive Control	
KS 2054 - N	Negative Control	

Xanthomonas campestris pv. begoniae

Xcb	Indirect ELISA	-- /Mono
CAB 23700/0500	Coating Ab	500 test
CAB 23700/1000	Coating Ab	1000 test
CAB 23700/5000	Coating Ab	5000 test
ECA 23700/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 23700/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 23700/5000	Conjugated Ab, Alk.Phos.	5000 test
BRA 23700/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 23700/1000	Bacterial Reagent Set, Alk.Phos.	1000 test
BRA 23700/5000	Bacterial Reagent Set, Alk.Phos.	5000 test
LNC 23700	Negative Control	
LPC 23700	Positive Control	

Xanthomonas campestris pv. begoniae

Xcb	Indirect ELISA/IF	
KS 2521 - 0500R	Reagent Set	500 test
KS 2521 - 1000R	Reagent Set	1000 test
KS 2521 - 5000R	Reagent Set	5000 test
KS 2521 - P	Positive Control	
KS 2521 - N	Negative Control	

Xanthomonas campestris pv. campestris

Xcc	DAS ELISA	
KS 2056 - 0500R	Reagent Set	500 test
KS 2056 - 1000R	Reagent Set	1000 test
KS 2056 - 5000R	Reagent Set	5000 test
KS 2056 - P	Positive Control	
KS 2056 - N	Negative Control	

Xanthomonas campestris pv. campestris

Xcc	Indirect ELISA	-- /Mono
CAB 97000/0500	Coating Ab	500 test
CAB 97000/1000	Coating Ab	1000 test
CAB 97000/5000	Coating Ab	5000 test
ECA 97000/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 97000/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 97000/5000	Conjugated Ab, Alk.Phos.	5000 test

BRA 97000/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 97000/1000	Bacterial Reagent Set, Alk.Phos.	1000 test
BRA 97000/5000	Bacterial Reagent Set, Alk.Phos.	5000 test
LNC 97000	Negative Control	
LPC 97000	Positive Control	

Xanthomonas campestris pv. carotae

Xccar	DAS ELISA	
KS 2061 - 0500R	Reagent Set	500 test
KS 2061 - 1000R	Reagent Set	1000 test
KS 2061 - 5000R	Reagent Set	5000 test
KS 2061 - P	Positive Control	
KS 2061 - N	Negative Control	

Xanthomonas campestris pv. carotae

Xccar	Indirect ELISA/IF	
KS 2060 - 0500R	Reagent Set	500 test
KS 2060 - 1000R	Reagent Set	1000 test
KS 2060 - 5000R	Reagent Set	5000 test
KS 2060 - P	Positive Control	
KS 2060 - N	Negative Control	

Xanthomonas campestris pv. citri

Xccit	DAS ELISA	
KS 2522 - 0500R	Reagent Set	500 test
KS 2522 - 1000R	Reagent Set	1000 test
KS 2522 - 5000R	Reagent Set	5000 test
KS 2522 - P	Positive Control	
KS 2522 - N	Negative Control	

Xanthomonas campestris pv. citri

Xccit	Indirect ELISA	-- /Mono
CAB 92200/0500	Coating Ab	500 test
CAB 92200/1000	Coating Ab	1000 test
CAB 92200/5000	Coating Ab	5000 test
ECA 92200/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 92200/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 92200/5000	Conjugated Ab, Alk.Phos.	5000 test
BRA 92200/0500	Bacterial Reagent Set, Alk.Phos.	500 test
BRA 92200/1000	Bacterial Reagent Set, Alk.Phos.	1000 test
BRA 92200/5000	Bacterial Reagent Set, Alk.Phos.	5000 test
LPC 92200	Positive Control	
LNC 92200	Negative Control	

Xanthomonas campestris pv. dieffenbachiae

Xcd	DAS ELISA	
KS 2063 - 0500R	Reagent Set	500 test
KS 2063 - 1000R	Reagent Set	1000 test
KS 2063 - 5000R	Reagent Set	5000 test
KS 2063 - P	Positive Control	

KS 2063 - N	Negative Control
-------------	------------------

Xanthomonas campestris pv. dieffenbachiae

Xcd	Indirect ELISA	-- /Mono
CAB 47600/0500	Coating Ab	500 test
CAB 47600/1000	Coating Ab	1000 test
CAB 47600/5000	Coating Ab	5000 test
ECA 47600/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 47600/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 47600/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 47600/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 47600/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 47600/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 47600	Positive Control	
LNC 47600	Negative Control	

Xanthomonas campestris pv. hyacinthi

Xch	DAS ELISA	
KS 2094 - 0500R	Reagent Set	500 test
KS 2094 - 1000R	Reagent Set	1000 test
KS 2094 - 5000R	Reagent Set	5000 test
KS 2094 - P	Positive Control	
KS 2094 - N	Negative Control	

Xanthomonas campestris pv. jacinti

Xcj	Indirect ELISA/IF	
KS 2099 - 0500R	Reagent Set	500 test
KS 2099 - 1000R	Reagent Set	1000 test
KS 2099 - 5000R	Reagent Set	5000 test
KS 2099 - P	Positive Control	
KS 2099 - N	Negative Control	

Xanthomonas campestris pv. oryzae

Xco	DAS ELISA	
KS 2069 - 0500R	Reagent Set	500 test
KS 2069 - 1000R	Reagent Set	1000 test
KS 2069 - 5000R	Reagent Set	5000 test
KS 2069 - P	Positive Control	
KS 2069 - N	Negative Control	

Xanthomonas campestris pv. oryzae

Xco	Indirect ELISA	-- /Mono
CAB 85000/0500	Coating Ab	500 test
CAB 85000/1000	Coating Ab	1000 test
CAB 85000/5000	Coating Ab	5000 test
ECA 85000/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 85000/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 85000/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 85000/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 85000/1000	Bacterial Reagent Set,Alk.Phos.	1000 test

BRA 85000/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 85000	Positive Control	
LNC 85000	Negative Control	

Xanthomonas campestris pv. pelargonii

Xcp	DAS ELISA	ily/Mono
CAB 32503/0500	Coating Ab	500 test
CAB 32503/1000	Coating Ab	1000 test
CAB 32503/5000	Coating Ab	5000 test
ECA 32503/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 32503/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 32503/5000	Conjugated Ab,Alk.Phos.	5000 test
PSA 32503/0096	PathoScreen Kit,Alk.Phos.	96 test
PSA 32503/0288	PathoScreen Kit,Alk.Phos.	288 test
PSA 32503/0480	PathoScreen Kit,Alk.Phos.	480 test
SRA 32503/0500	Reagent Set,Alk.Phos.	500 test
SRA 32503/1000	Reagent Set,Alk.Phos.	1000 test
SRA 32503/5000	Reagent Set,Alk.Phos.	5000 test
*ISK 32503/0005	ImmunoStrip Test Kit	5 strips
*ISK 32503/0025	ImmunoStrip Test Kit	25 strips
*ISK 32503/0250	ImmunoStrip Test Kit	250 strips
*STX32503/0005	ImmunoStrip Test	5 strips
*STX32503/0025	ImmunoStrip Test	25 strips
*STX32503/0250	ImmunoStrip Test	250 strips
LNC 32503	Negative Control	
LPC 32503	Positive Control	

*** ImmunoStrip Test kit**

Including strips and buffer filled mesh bags

*** Xcp ImmunoStrip Buffer**

Xcp ImmunoStrips must be used with SEB1 buffer which is offered separately.

Please see accessories for a listing of SEB1 buffer products.

Xanthomonas campestris pv. pelargonii

Xcp	Indirect ELISA	-- /Mono
CAB 32502/0500	Coating Ab	500 test
CAB 32502/1000	Coating Ab	1000 test
CAB 32502/5000	Coating Ab	5000 test
ECA 32502/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 32502/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 32502/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 32502/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 32502/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 32502/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LPC 32502	Positive Control	

Xanthomonas campestris pv. phaseoli

Xcph	DAS ELISA	oly/Poly
CAB 36101/0500	Coating Ab	500 test
CAB 36101/1000	Coating Ab	1000 test

식물병 Bacteria

CAB 36101/5000	Coating Ab	5000 test
ECA 36101/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 36101/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 36101/5000	Conjugated Ab,Alk.Phos.	5000 test
SRA 36101/0500	Reagent Set,Alk.Phos.	500 test
SRA 36101/1000	Reagent Set,Alk.Phos.	1000 test
SRA 36101/5000	Reagent Set,Alk.Phos.	5000 test
LNC 36101	Negative Control	
LPC 36101	Positive Control	

Xanthomonas campestris pv. phaseoli

Xcph	Indirect ELISA	
KS 2077 - 0500R	Reagent Set	500 test
KS 2077 - 1000R	Reagent Set	1000 test
KS 2077 - 5000R	Reagent Set	5000 test
KS 2077 - P	Positive Control	
KS 2077 - N	Negative Control	

Xanthomonas campestris pv. phaseoli (fuscans-type)

Xcphf	DAS ELISA	
KS 2079 - 0500R	Reagent Set	500 test
KS 2079 - 1000R	Reagent Set	1000 test
KS 2079 - 5000R	Reagent Set	5000 test
KS 2079 - P	Positive Control	
KS 2079 - N	Negative Control	

Xanthomonas campestris pv. pruni

Xcpr	Indirect ELISA	
KS 2081 - 0500R	Reagent Set	500 test
KS 2081 - 1000R	Reagent Set	1000 test
KS 2081 - 5000R	Reagent Set	5000 test
KS 2081 - P	Positive Control	
KS 2081 - N	Negative Control	

Xanthomonas campestris pv. vesicatoria

Xcv	DAS ELISA	
KS 2084 - 0500R	Reagent Set	500 test
KS 2084 - 1000R	Reagent Set	1000 test
KS 2084 - 5000R	Reagent Set	5000 test
KS 2084 - P	Positive Control	
KS 2084 - N	Negative Control	

Xanthomonas campestris pv. vesicatoria

Xcv	Indirect ELISA	-- /Mono
CAB 42300/0500	Coating Ab	500 test
CAB 42300/1000	Coating Ab	1000 test
CAB 42300/5000	Coating Ab	5000 test
ECA 42300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 42300/1000	Conjugated Ab,Alk.Phos.	1000 test

ECA 42300/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 42300/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 42300/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 42300/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LNC 42300	Negative Control	
LPC 42300	Positive Control	

Xanthomonas fragariae

Xf	DAS ELISA	
KS 2085 - 0500R	Reagent Set	500 test
KS 2085 - 1000R	Reagent Set	1000 test
KS 2085 - 5000R	Reagent Set	5000 test
KS 2085 - P	Positive Control	
KS 2085 - N	Negative Control	

Xanthomonas fragariae

Xf	Indirect ELISA	
KS 2087 - 0500R	Reagent Set	500 test
KS 2087 - 1000R	Reagent Set	1000 test
KS 2087 - 5000R	Reagent Set	5000 test
KS 2087 - P	Positive Control	
KS 2087 - N	Negative Control	

Xanthomonas hortorum pv. pelargonii

Xhp	DAS ELISA	
KS 2095 - 0500R	Reagent Set	500 test
KS 2095 - 1000R	Reagent Set	1000 test
KS 2095 - 5000R	Reagent Set	5000 test
KS 2095 - P	Positive Control	
KS 2095 - N	Negative Control	

Xanthomonas hortorum pv. pelargonii

Xhp	Indirect ELISA/IF	
KS 2544 - 0500R	Reagent Set	500 test
KS 2544 - 1000R	Reagent Set	1000 test
KS 2544 - 5000R	Reagent Set	5000 test
KS 2544 - P	Positive Control	
KS 2544 - N	Negative Control	

Xanthomonas maltophilia

Xm	DAS ELISA	
KS 2538 - 0500R	Reagent Set	500 test
KS 2538 - 1000R	Reagent Set	1000 test
KS 2538 - 5000R	Reagent Set	5000 test
KS 2538 - P	Positive Control	
KS 2538 - N	Negative Control	

Xanthomonas maltophilia

Xm	Indirect ELISA	-- /Mono
CAB 31300/0500	Coating Ab	500 test

CAB 31300/1000	Coating Ab	1000 test
CAB 31300/5000	Coating Ab	5000 test
ECA 31300/0500	Conjugated Ab,Alk.Phos.	500 test
ECA 31300/1000	Conjugated Ab,Alk.Phos.	1000 test
ECA 31300/5000	Conjugated Ab,Alk.Phos.	5000 test
BRA 31300/0500	Bacterial Reagent Set,Alk.Phos.	500 test
BRA 31300/1000	Bacterial Reagent Set,Alk.Phos.	1000 test
BRA 31300/5000	Bacterial Reagent Set,Alk.Phos.	5000 test
LNC 31300	Negative Control	
LPC 31300	Positive Control	

Xylella fastidiosa

Xf	DAS ELISA	Poly/Poly
CAB 34501/0500	Coating Ab	500 test
CAB 34501/1000	Coating Ab	1000 test
CAB 34501/5000	Coating Ab	5000 test
ECP 34501/0500	Conjugated Ab,Peroxidase	500 test
ECP 34501/1000	Conjugated Ab,Peroxidase	1000 test
ECP 34501/5000	Conjugated Ab,Peroxidase	5000 test
PSP 34501/0096	PathoScreen Kit,Peroxidase	96 test
PSP 34501/0288	PathoScreen Kit,Peroxidase	288 test
PSP 34501/0480	PathoScreen Kit,Peroxidase	480 test
SRP 34501/0500	Reagent Set,Peroxidase	500 test
SRP 34501/1000	Reagent Set,Peroxidase	1000 test
SRP 34501/5000	Reagent Set,Peroxidase	5000 test
LNC 34501	Negative Control	
LPC 34501	Positive Control	

Xylophilus ampelinus (bacterial necrosis)

	DAS ELISA	
KS 2090 - 0500R	Reagent Set	500 test
KS 2090 - 1000R	Reagent Set	1000 test
KS 2090 - 5000R	Reagent Set	5000 test
KS 2090 - P	Positive Control	
KS 2090 - N	Negative Control	

Ascochyta pisi

	Indirect ELISA	
KS 2101 - 0500R	Reagent Set	500 test
KS 2101 - 1000R	Reagent Set	1000 test
KS 2101 - 5000R	Reagent Set	5000 test
KS 2101 - P	Positive Control	
KS 2101 - N	Negative Control	

Aspergillus & penicillium spp.

AsPe	Indirect ELISA	
KS 2102 - 0500R	Reagent Set	500 test
KS 2102 - 1000R	Reagent Set	1000 test
KS 2102 - 5000R	Reagent Set	5000 test
KS 2102 - 0500S	Full Kit	500 test
KS 2102 - 1000S	Full Kit	1000 test
KS 2102 - 5000S	Full Kit	5000 test
KS 2102 - P	Positive Control	
KS 2102 - N	Negative Control	

Botrytis allii

Boa	Indirect ELISA	
KS 2129 - 0500R	Reagent Set	500 test
KS 2129 - 1000R	Reagent Set	1000 test
KS 2129 - 5000R	Reagent Set	5000 test
KS 2129 - P	Positive Control	
KS 2129 - N	Negative Control	

Botrytis cinerea

Boc	DAS ELISA	
KS 2130 - 0500R	Reagent Set	500 test
KS 2130 - 1000R	Reagent Set	1000 test
KS 2130 - 5000R	Reagent Set	5000 test
KS 2130 - P	Positive Control	
KS 2130 - N	Negative Control	

Botrytis cinerea

Boc	Indirect ELISA	
KS 2103 - 0500R	Reagent Set	500 test
KS 2103 - 1000R	Reagent Set	1000 test
KS 2103 - 5000R	Reagent Set	5000 test
KS 2103 - P	Positive Control	
KS 2103 - N	Negative Control	

Cladosporium spp.

	Indirect ELISA	
KS 2104 - 0500R	Reagent Set	500 test
KS 2104 - 1000R	Reagent Set	1000 test
KS 2104 - 5000R	Reagent Set	5000 test
KS 2104 - P	Positive Control	
KS 2104 - N	Negative Control	

Colletotrichum acutatum

Coa	Indirect ELISA	
KS 2105 - 0500R	Reagent Set	500 test
KS 2105 - 1000R	Reagent Set	1000 test
KS 2105 - 5000R	Reagent Set	5000 test
KS 2105 - P	Positive Control	
KS 2105 - N	Negative Control	

Colletotrichum spp.

Cos	Indirect ELISA	
KS 2106 - 0500R	Reagent Set	500 test
KS 2106 - 1000R	Reagent Set	1000 test
KS 2106 - 5000R	Reagent Set	5000 test
KS 2106 - P	Positive Control	
KS 2106 - N	Negative Control	

Ergosterol (Total Fungus detection)

	Competitive ELISA	
KS 2107 - 0500R	Reagent Set	500 test
KS 2107 - 1000R	Reagent Set	1000 test
KS 2107 - 5000R	Reagent Set	5000 test
KS 2107 - P	Positive Control	
KS 2107 - N	Negative Control	

Humicola lanuginosa

	Indirect ELISA	
KS 2108 - 0500R	Reagent Set	500 test
KS 2108 - 1000R	Reagent Set	1000 test
KS 2108 - 5000R	Reagent Set	5000 test
KS 2108 - P	Positive Control	
KS 2108 - N	Negative Control	

Laccase

	Indirect ELISA	
KS 2109 - 0500R	Reagent Set	500 test
KS 2109- 1000R	Reagent Set	1000 test
KS 2109 - 5000R	Reagent Set	5000 test
KS 2109 - P	Positive Control	
KS 2109 - N	Negative Control	

Monilinia fructicola

Mof	Indirect ELISA	
KS 2110 - 0500R	Reagent Set	500 test
KS 2110 - 1000R	Reagent Set	1000 test
KS 2110 - 5000R	Reagent Set	5000 test
KS 2110 - P	Positive Control	
KS 2110 - N	Negative Control	

Monilinia spp.

Mos	Indirect ELISA	
KS 2111 - 0500R	Reagent Set	500 test
KS 2111 - 1000R	Reagent Set	1000 test
KS 2111 - 5000R	Reagent Set	5000 test
KS 2111 - P	Positive Control	
KS 2111 - N	Negative Control	

Mouldetect

	DAS ELISA	
KS 2112 - 0500R	Reagent Set	500 test
KS 2112 - 1000R	Reagent Set	1000 test
KS 2112 - 5000R	Reagent Set	5000 test
KS 2112 - P	Positive Control	
KS 2112 - N	Negative Control	

Mycosphasrella pinodes

	Indirect ELISA	
KS 2113 - 0500R	Reagent Set	500 test
KS 2113 - 1000R	Reagent Set	1000 test
KS 2113 - 5000R	Reagent Set	5000 test
KS 2113 - P	Positive Control	
KS 2113 - N	Negative Control	

Nectria galligena

	Indirect ELISA	
KS 2114 - 0500R	Reagent Set	500 test
KS 2114 - 1000R	Reagent Set	1000 test
KS 2114 - 5000R	Reagent Set	5000 test
KS 2114 - P	Positive Control	
KS 2114 - N	Negative Control	

Phoma medicaginis pinodella

Mos	Indirect ELISA	
KS 2115 - 0500R	Reagent Set	500 test
KS 2115 - 1000R	Reagent Set	1000 test
KS 2115 - 5000R	Reagent Set	5000 test
KS 2115 - P	Positive Control	
KS 2115 - N	Negative Control	

Phytophthora

Phyt	DAS ELISA	Poly/Mono
CAB 92600/0500	Coating ab	500 test
CAB 92600/1000	Coating ab	1000 test
CAB 92600/5000	Coating ab	5000 test
ECA 92600/0500	Conjugated ab,Alk.Phos.	500 test
ECA 92600/1000	Conjugated ab,Alk.Phos.	1000 test
ECA 92600/5000	Conjugated ab,Alk.Phos.	5000 test
SRA 92600/0500	Reagent Set,Alk.Phos.	500 test
SRA 92600/1000	Reagent Set,Alk.Phos.	1000 test

SRA 92600/5000	Reagent Set, Alk.Phos.	5000 test
LNC 92600	Negative Control	
LPC 92600	Positive Control	

Phytophthora infestans**DAS ELISA**

KS 2116 - 0500R	Reagent Set	500 test
KS 2116 - 1000R	Reagent Set	1000 test
KS 2116 - 5000R	Reagent Set	5000 test
KS 2116 - P	Positive Control	
KS 2116 - N	Negative Control	

Phythium species.**Express kit**

KS 2138 - 0090E	Agriscreen Express Kit	90 test
KS 2138 - 0008E	Alert Express Kit	8 test

Pseudocercospora (Eyespot)

Pse	Indirect ELISA	
KS 2117 - 0500R	Reagent Set	500 test
KS 2117 - 1000R	Reagent Set	1000 test
KS 2117 - 5000R	Reagent Set	5000 test
KS 2117 - 0096S	Full Kit precoated plate	96 test
KS 2117 - 0500S	Full Kit	500 test
KS 2117 - 1000S	Full Kit	1000 test
KS 2117 - 5000S	Full Kit	5000 test
KS 2117 - P	Positive Control	
KS 2117 - N	Negative Control	

Pythium ultimum

Pult	DAS ELISA	
KS 2141 - 0500R	Reagent Set	500 test
KS 2141 - 1000R	Reagent Set	1000 test
KS 2141 - 5000R	Reagent Set	5000 test

KS 2141 - P	Positive Control
KS 2141 - N	Negative Control

Pythium ultimum

Pult	Compound Direct ELISA	Poly/Mono
CAB 92500/0500	Coating ab	500 test
CAB 92500/1000	Coating ab	1000 test
CAB 92500/5000	Coating ab	5000 test
ECA 92500/0500	Conjugated ab, Alk.Phos.	500 test
ECA 92500/1000	Conjugated ab, Alk.Phos.	1000 test
ECA 92500/5000	Conjugated ab, Alk.Phos.	5000 test
SRA 92500/0500	Reagent Set, Alk.Phos.	500 test
SRA 92500/1000	Reagent Set, Alk.Phos.	1000 test
SRA 92500/5000	Reagent Set, Alk.Phos.	5000 test
LNC 92500	Negative Control	

Rhizoctonia solani I

Rh. Solani I	DAS ELISA	
KS 2126 - 0500R	Reagent Set	500 test
KS 2126 - 1000R	Reagent Set	1000 test
KS 2126 - 5000R	Reagent Set	5000 test
KS 2126 - P	Positive Control	
KS 2126 - N	Negative Control	

Rhizoctonia solani II

Rh. Solani II	DAS ELISA	
KS 2127 - 0500R	Reagent Set	500 test
KS 2127 - 1000R	Reagent Set	1000 test
KS 2127 - 5000R	Reagent Set	5000 test
KS 2127 - P	Positive Control	
KS 2127 - N	Negative Control	

Rhizoctonia spp.

Rhiz	DAS ELISA	
KS 2144 - 0500R	Reagent Set	500 test
KS 2144 - 1000R	Reagent Set	1000 test
KS 2144 - 5000R	Reagent Set	5000 test
KS 2144 - P	Positive Control	
KS 2144 - N	Negative Control	

Rhizoctonia species.

Rhiz	Express kit	
KS 2145 - 0090E	Agriscreen Express Kit	90 test
KS 2145 - 0008E	Alert Express Kit	8 test

Phytophthora +Pythium +Rhizoctonia species.

PPR	Express kit	
KS 2146 - 0012E	Alert Express Kit	3 x 4 test

Sclerotinia sclerotiorum

Scs	Indirect ELISA	
KS 2118 - 0500R	Reagent Set	500 test
KS 2118 - 1000R	Reagent Set	1000 test
KS 2118 - 5000R	Reagent Set	5000 test
KS 2118 - 0096S	Full Kit precoated plate	96 test
KS 2118 - 0500S	Full Kit	500 test
KS 2118 - 1000S	Full Kit	1000 test
KS 2118 - 5000S	Full Kit	5000 test
KS 2118 - P	Positive Control	
KS 2118 - N	Negative Control	

Septoria nodorum + Septoria tritici

Sent	DAS ELISA	
KS 2119 - 0096S	Full Kit precoated plate	2 x 96 testls

Septoria nodorum

Sen	DAS ELISA	
KS 2120 - 0500R	Reagent Set	500 test
KS 2120 - 1000R	Reagent Set	1000 test
KS 2120 - 5000R	Reagent Set	5000 test
KS 2120 - P	Positive Control	
KS 2120 - N	Negative Control	

Septoria tritici

Set	DAS ELISA	
KS 2128 - 0500R	Reagent Set	500 test
KS 2128 - 1000R	Reagent Set	1000 test
KS 2128 - 5000R	Reagent Set	5000 test
KS 2128 - P	Positive Control	
KS 2128 - N	Negative Control	

Tapesia spp. (Eyespot)

	DAS ELISA	
KS 2149 - 0500R	Reagent Set	500 test
KS 2149 - 1000R	Reagent Set	1000 test
KS 2149 - 5000R	Reagent Set	5000 test
KS 2149 - P	Positive Control	
KS 2149 - N	Negative Control	

Verticillium dahliae

Vdah	DAS ELISA	
KS 2123 - 0500R	Reagent Set	500 test
KS 2123 - 1000R	Reagent Set	1000 test
KS 2123 - 5000R	Reagent Set	5000 test
KS 2123 - P	Positive Control	
KS 2123 - N	Negative Control	

Verticillium dahliae

Vdah	Compound Direct	
CAB 92300/0500	Coating ab	500 test
CAB 92300/1000	Coating ab	1000 test
CAB 92300/5000	Coating ab	5000 test
ECA 92300/0500	Conjugated ab,Alk.Phos.	500 test
ECA 92300/1000	Conjugated ab,Alk.Phos.	1000 test
ECA 92300/5000	Conjugated ab,Alk.Phos.	5000 test
SRA 92300/0500	Reagent Set,Alk.Phos.	500 test
SRA 92300/1000	Reagent Set,Alk.Phos.	1000 test
SRA 92300/5000	Reagent Set,Alk.Phos.	5000 test
LNC 92300	Negative Control	
LPC 92300	Positive Control	

Verticillium spec. rape isolate

	DAS ELISA	
KS 2124 - 0500R	Reagent Set	500 test
KS 2124 - 1000R	Reagent Set	1000 test
KS 2124 - 5000R	Reagent Set	5000 test
KS 2124 - P	Positive Control	
KS 2124 - N	Negative Control	

Verticillium spec. strawberry isolate

	Indirect ELISA	
KS 2125 - 0500R	Reagent Set	500 test
KS 2125 - 1000R	Reagent Set	1000 test
KS 2125 - 5000R	Reagent Set	5000 test
KS 2125 - P	Positive Control	
KS 2125 - N	Negative Control	

Verticillium spp.

	DAS ELISA	
KS 2148 - 0500R	Reagent Set	500 test
KS 2148 - 1000R	Reagent Set	1000 test
KS 2148 - 5000R	Reagent Set	5000 test
KS 2148 - P	Positive Control	
KS 2148 - N	Negative Control	

Apple Proliferation Phytoplasma

ApP	DAS ELISA	
KS 2303 - 0500R	Reagent Set	500 test
KS 2303 - 1000R	Reagent Set	1000 test
KS 2303 - 5000R	Reagent Set	5000 test
KS 2303 - P	Positive Control	
KS 2303 - N	Negative Control	

Aster Yellows Phytoplasma**(Set with primers, dNTP and controls)**

KS 2301 - 0100R	Reagent Set	100 test
KS 2301 - P	Positive Control	
KS 2301 - N	Negative Control	

Corn Stunt Spiroplasma

CSS	DAS ELISA	Poly/Poly
CAB 16000/0500	Coating Ab	500 test
CAB 16000/1000	Coating Ab	1000 test
CAB 16000/5000	Coating Ab	5000 test
ECA 16000/0500	Conjugated Ab, Alk.Phos.	500 test
ECA 16000/1000	Conjugated Ab, Alk.Phos.	1000 test
ECA 16000/5000	Conjugated Ab, Alk.Phos.	5000 test
SRA 16000/0500	Reagent Set, Alk.Phos.	500 test
SRA 16000/1000	Reagent Set, Alk.Phos.	1000 test
SRA 16000/5000	Reagent Set, Alk.Phos.	5000 test
LNC 16000	Negative control	
LPC 16000	Positive Control	
LNC 16000	Negative Control	

Elm Yellows Phytoplasma**(Set with primers, dNTP and controls)**

KS 2307 - 0100R	Reagent Set	100 test
KS 2307 - P	Positive Control	
KS 2307 - N	Negative Control	

Flavescence doree Phytoplasma

FD	DAS ELISA	
KS 2312 - 0500R	Reagent Set, Biotin	500 test
KS 2312 - 1000R	Reagent Set, Biotin	1000 test
KS 2312 - 5000R	Reagent Set, Biotin	5000 test
KS 2312 - 0500S	Full Kit, Biotin	500 test
KS 2312 - 1000S	Full Kit, Biotin	1000 test
KS 2312 - 5000S	Full Kit, Biotin	5000 test
KS 2312 - P	Positive Control	
KS 2312 - N	Negative Control	

* A special buffer is required for grapevine extraction.

Stolbur Phytoplasma (Bosi noir)

BN	IF	
KS 2308 - 0100R	Reagent Set	100 test

KS 2308 - P	Positive Control
KS 2308 - N	Negative Control

Stolbur Phytoplasma**(including bois noir, Grapevine Yellows)**

KS 2309 - 0100R	Reagent Set	100 test
KS 2309 - P	Positive Control	
KS 2309 - N	Negative Control	

Universal Phytoplasma**DAS ELISA**

KS 2310 - 0100R	Reagent Set	100 test
KS 2310 - P	Positive Control	
KS 2310 - N	Negative Control	

Universal Phytoplasma nested**DAS ELISA**

KS 2311 - 0100R	Reagent Set	100 test
KS 2311 - P	Positive Control	
KS 2311 - N	Negative Control	

GMO 진단 Kit

Transgenic crops에서 발현되는 특정 단백질을 detection할 수 있는 kit 입니다.

Proteins Detected	Trade Name/Crop	Format	Cat.No
Bt-Cry1Ab/Ac	Bollgard & Bollgard II Cotton YieldGard Corn Borer Corn NatureGard(E176), KnockOut (E176) Corn Bt-Xtra Corn	ELISA & ImmunoStrip	06200
Bt-Cry2A	Bollgard II Cotton	ELISA & ImmunoStrip	05800
Bt-Cry2A and Bt-Cry1Ab/Ac	Bollgard II Cotton	ImmunoStrip	06800
Bt-Cry3Bb1	YieldGard Rootworm Corn	ELISA & ImmunoStrip	06100
Bt-Cry9C	Starlink Corn	ELISA & ImmunoStrip	05600
CP4 EPSPS	Roundup Ready, Corn, Cotton, Soybean, Canola, Alfalfa, and Sugarbeet	ELISA & ImmunoStrip	74000
Bt-Cry1Ab/Ac and CP4 EPSPS	Bollgard Cotton with Roundup Ready Yieldgard Corn with Roundup Ready	ImmunoStrip	74500
Bt-Cry2A, Bt-Cry1Ab/1Ac and CP4 EPSPS	Bollgard II Cotton with Roundup Ready	ImmunoStrip	51000

Bt-Cry1Ab/1Ac

Bt-Cry1Ab/1Ac	DAS ELISA/ImmunoStrip
PSP 06200/0096	PathoScreen Kit, Peroxidase 96test
PSP 06200/0288	PathoScreen Kit, Peroxidase 288test
PSP	
PSP 06200/4800	PathoScreen Kit, Peroxidase 4800test
*STX 06200/0012	Comb Format ImmunoStrip Test 12 combs
*STX 06200/0050	ImmunoStrip Test 50 strips
*STX 06200/0250	ImmunoStrip Test 250 strips
LNC 06200	Negative Control

LPC 06200 Positive Control

*** Bt-Cry1Ab/1Ac ImmunoStrip Buffer**

Bt-Cry1Ab/1Ac ImmunoStrips must be used with SEB4 buffer which is offered separately.

Please see accessories for a listing of SEB4 products and sizes.

Bt-Cry1Ab/1Ac and RoundupReady (CP4 EPSPS) ImmunoStrip

*STX 74500/0012	Comb Format ImmunoStrip Test	12 combs
*STX 74500/0050	ImmunoStrip Test	50 strips
*STX 74500/0250	ImmunoStrip Test	250 strips

LPC 05700	Positive Control
-----------	------------------

Bt-Cry1F protein

Bt-Cry1F	ImmunoStrip	
*STX 10300/0012	Comb Format ImmunoStrip Test	12 combs
*STX 10300/0050	ImmunoStrip Test	50 strips
*STX 10300/0250	ImmunoStrip Test	250 strips

*Sample Extract Buffer sold separately

Bt-Cry1F and Bt-Cry 1Ac protein

ImmunoStrip		
*STX 12900/0012	Comb Format ImmunoStrip Test	12 combs
*STX 12900/0050	ImmunoStrip Test	50 strips
*STX 12900/0250	ImmunoStrip Test	250 strips

*Sample Extract Buffer sold separately

Bt-Cry2A protein

Bt-Cry2A	DAS ELISA/ImmunoStrip Poly/Poly	
PSB 05800/0096	PathoScreen Kit,Alk.Phos, BP substrate	96test
PSB 05800/0288	PathoScreen Kit,Alk.Phos, BP substrate	288test
PSB 05800/0480	PathoScreen Kit,Alk.Phos, BP substrate	480test
PSB 05800/4800	PathoScreen Kit,Alk.Phos, BP substrate	4800test
*STX 05800/0012	Comb Format ImmunoStrip Test	12 combs
*STX 05800/0050	ImmunoStrip Test	50 strips
*STX 05800/0250	ImmunoStrip Test	250 strips
LNC 05800	Negative Control	
LPC 05800	Positive Control	

* Bt-Cry2A ImmunoStrip Buffer

Bt-Cry2A ImmunoStrips must be used with SEB4 buffer which is offered separately. Please see accessories for a listing of SEB4 buffer products and sizes.

Bt-Cry2A and Bt-Cry 1Ac

ImmunoStrip		
*STX 06800/0012	Comb Format ImmunoStrip Test	12 combs
*STX 06800/0050	ImmunoStrip Test	50 strips
*STX 06800/0250	ImmunoStrip Test	250 strips

*Sample Extract Buffer sold separately

Bt-Cry2A, Bt-Cry 1Ac and Roundup Ready

ImmunoStrip		
*STX 51000/0012	Comb Format ImmunoStrip Test	12 combs
*STX 51000/0050	ImmunoStrip Test	50 strips
*STX 51000/0250	ImmunoStrip Test	250 strips

*Sample Extract Buffer sold separately

Bt-Cry3A protein

Bt-Cry3A	DAS ELISA/ImmunoStrip Poly/Poly	
PSA 05900/0096	PathoScreen Kit, Alk.Phos	96test
PSA 05900/0288	PathoScreen Kit, Alk.Phos	288test
PSA 05900/0480	PathoScreen Kit, Alk.Phos	480test
PSA 05900/4800	PathoScreen Kit, Alk.Phos	4800test
PSB 05900/0096	PathoScreen Kit,Alk.Phos,	96test
PSB 05900/0288	PathoScreen Kit,Alk.Phos,	288test
PSB 05900/0480	PathoScreen Kit,Alk.Phos,	480test
PSB 05900/4800	PathoScreen Kit,Alk.Phos,	4800test
*STX 05900/0050	ImmunoStrip Test	50 strips
*STX 05900/0250	ImmunoStrip Test	250 strips
LNC 05900	Negative Control	
LPC 05900	Positive Control	

*Bt-Cry3A ImmunoStrip Buffer

Bt-Cry3A ImmunoStrips must be used with SEB4 buffer which is offered separately. Please use one of the SEB4 products with the Bt-Cry3A ImmunoStrips.

Bt-Cry3A protein

Bt-Cry3A	Eeay-to-use kit	
NMA 05901/0096	Bt-Cry3A kit, Alk.Phos	96test
LPC 05901	Positive Control	

Bt-Cry3Bb1 protein

Bt-Cry3Bb1	DAS ELISA/ImmunoStrip	
PSP 06100/0096	PathoScreen Kit, Peroxidase	96test
PSP 06100/0288	PathoScreen Kit, Peroxidase	288test
PSP 06100/0480	PathoScreen Kit, Peroxidase	480test
PSP 06100/4800	PathoScreen Kit, Peroxidase	4800test
*STX 06100/0012	Comb Format ImmunoStrip Test	12 combs
*STX 06100/0050	ImmunoStrip Test	50 strips
*STX 06100/0250	ImmunoStrip Test	250 strips
LNC 06100	Negative Control	
LPC 06100	Positive Control	

*Bt-Cry3Bb1 ImmunoStrip Buffer

Bt-Cry3Bb1 ImmunoStrips must be used with SEB4 buffer which is offered separately. Please use one of the SEB4 products with the Bt-Cry3Bb1 ImmunoStrips.

Bt-Cry3Bb1 and Bt-Cry1Ab

ImmunoStrip		
*STX 11200/0012	Comb Format ImmunoStrip Test	12 combs
*STX 11200/0050	ImmunoStrip Test	50 strips
*STX 11200/0250	ImmunoStrip Test	250 strips

*Sample Extract Buffer sold separately

Bt-Cry9C protein

Bt-Cry9C	DAS ELISA/ImmunoStrip Poly/Poly		
PSP 05600/0096	PathoScreen Kit, Peroxidase		96test
PSP 05600/0288	PathoScreen Kit, Peroxidase		288test
PSP 05600/0480	PathoScreen Kit, Peroxidase		480test
PSP 05600/4800	PathoScreen Kit, Peroxidase		4800test
*STX 05600/0012	Comb Format ImmunoStrip Test		12 combs
*STX 05600/0050	ImmunoStrip Test		50 strips
*STX 05600/0250	ImmunoStrip Test		250 strips
LNC 05600	Negative Control		
LPC 05600	Positive Control		

Bacillus thuringiensis subsp. kurstaki

Btk	DAS ELISA		
KS 2401 - 0500R	Reagent Set		500 test
KS 2401 - 1000R	Reagent Set		1000 test
KS 2401 - 5000R	Reagent Set		5000 test
KS 2401 - P	Positive Control		
KS 2401 - N	Negative Control		

Neomycin Phosphotransferase II

NPT II	Compound Direct ELISA		Poly/Mono
PSP 73000/0096	PathoScreen Kit, Peroxidase		96test
PSP 73000/0288	PathoScreen Kit, Peroxidase		288test
PSP 73000/0480	PathoScreen Kit, Peroxidase		480test
PSP 73000/4800	PathoScreen Kit, Peroxidase		4800test
LNC 73000	Negative Control		
LPC 73000	Positive Control		

Roundup Ready (CP4 EPSPS)

CP4-EPSPS	DAS ELISA/ImmunoStrip		
PSP 74000/0096	PathoScreen Kit, Peroxidase		96test
PSP 74000/0288	PathoScreen Kit, Peroxidase		288test
PSP 74000/0480	PathoScreen Kit, Peroxidase		480test
PSP 74000/4800	PathoScreen Kit, Peroxidase		4800test
*STX 74000/0012	Comb Format ImmunoStrip Test		12 combs
*STX 74000/0050	ImmunoStrip Test		50 strips
*STX 74000/0250	ImmunoStrip Test		250 strips
LNC 74000	Negative Control		
LPC 74000	Positive Control		

*Roundup Ready ImmunoStrip Buffer

Roundup Ready ImmunoStrips must be used with SEB4 buffer which is offered separately. Please use one of the SEB4 products with the Roundup Ready ImmunoStrips. Please see accessories for SEB4 buffer products and sizes.

INSECT ID

Helicoverpa/Heliothis - ID

Hel-ID	Indirect ELISA	-- /Mono
AGD 72500/0088	Insect egg ID Kit	

NEMATODES

Globodera rostochiensis/pallida

Indirect ELISA		
KS 2201 - 0500R	Reagent Set	500 test
KS 2201 - 1000R	Reagent Set	1000 test
KS 2201 - P	Positive Control	
KS 2201 - N	Negative Control	


Plant Growth Hormone



식물에 **Plant Growth Hormone** 이 어느 정도 있는지 **Phytodetek Kit**를 이용하여 **ELISA** 방식으로 정량하는 제품입니다.

- Stress, mutation events, tissue differentiation and other aspects of plant development에 의해 분비되는 PGH levels을 monitoring 할 수 있습니다.
- Accuracy – 비슷한 compounds에 대한 cross reactivity가 거의 없습니다.
- Sensitivity - 0.2-5.0 picomoles 정도까지 PGH 농도를 잴 수 있습니다.

<제품구성>

- 1) phytodetek kit – 실험에 필요한 모든 제품이 들어있습니다.
Antibody coated testwells, Tracer, Substrate, Wash solution, Stop reagent
- 2) phytodetek, Tracer (conjugated antibody)
- 3) phytodetek, monoclonal Antibody (coating antibody)



- 1) Antibody coated well 
- 2) Put plant extract 
- 3) Put tracer (Antibody-conjugate)
- 4) Read

Abscisic Acid

ABA	Competitive ELISA	1ono/Mono
PDM 09347/0096	Phytodetek, Monoclonal Ab	96test
PDA 09347/0032	Phytodetek, Tracer	32test
PDK 09347/0096	Phytodetek Kit, Alk. Phos.	96test

Dihydrozeatin Riboside

DHZR	Competitive ELISA	1ono/Mono
PDM 09349/0096	Phytodetek, Monoclonal Ab	96test
PDA 09349/0032	Phytodetek, Tracer	32test
PDK 09349/0096	Phytodetek Kit, Alk. Phos.	96test

Indole-3-acetic Acid

IAA	Competitive ELISA	1ono/Mono
PDM 09346/0096	Phytodetek, Monoclonal Ab	96test
PDA 09346/0032	Phytodetek, Tracer	32test
PDK 09346/0096	Phytodetek Kit, Alk. Phos.	96test

Isopentenyladenosine

IPA	Competitive ELISA	1ono/Mono
PDM 09351/0096	Phytodetek, Monoclonal Ab	96test
PDA 09351/0032	Phytodetek, Tracer	32test
PDK 09351/0096	Phytodetek Kit, Alk. Phos.	96test

trans-Zeatin Riboside

t-ZR	Competitive ELISA	1ono/Mono
PDM 09348/0096	Phytodetek, Monoclonal Ab	96test
PDA 09348/0032	Phytodetek, Tracer	32test
PDK 09348/0096	Phytodetek Kit, Alk. Phos.	96test

CHEMICALS AND BUFFERS

3M sodium hydroxide

Reaction stopping solution for Alk. Phos.

Conjugated tests

ACC 00409/0015	stop solution, 15ml
ACC 00409/0030	stop solution, 30ml
ACC 00409/0125	stop solution, 125ml

3M sulfuric acid

Reaction stopping solution for Peroxidase

conjugated tests

ACC 00402/0015	stop solution, 15ml
ACC 00402/0030	stop solution, 30ml
ACC 00402/0125	stop solution, 125ml

Ames buffer

Buffer for CSVd and PSTVd tests

ACC 00910	1X, ready for use, 60ml
-----------	-------------------------

APLPV EC buffer

ACC 00320	powder, 8.5g, makes 200ml
-----------	---------------------------

BCIP substrate

Substrate used in Xcp Immunoblot Kit

ACC 00425	1X, ready for use, 40ml
-----------	-------------------------

Bovine serum albumin

ACC 00300/0005	BSA, 5g
ACC 00300/0010	BSA, 10g
ACC 00300/0050	BSA, 50g
ACC 00300/0100	BSA, 100g
ACC 00300/0500	BSA, 500g
ACC 00300/1000	BSA, 1000g

BP substrate solution, alkaline phosphatase

BP substrate solution - Alk. Phos. substrate

ACC 00290/0015	1X, bottles A and B, 15ml each
ACC 00290/0033	1X, bottles A and B, 33ml each
ACC 00290/0055	1X, bottles A and B, 55ml each
ACC 00290/0500	1X, bottles A and B, 500ml each

Buffer pack, Bacterial

For use with Bacterial Reagent sets. Includes coating buffer, substrate, extraction buffer, PBST wash buffer, conjugate diluent, and stop solution.

ACC 00112/0500	Alk. Phos., 500 tests
ACC 00112/1000	Alk. Phos., 1000 tests
ACC 00112/5000	Alk. Phos., 5000 tests

Buffer pack, DAS/Compound, Alk. Phos.

For use with alk. phos. DAS and compound direct Reagent Sets. Includes coating buffer, substrate, extraction buffer, PBST wash buffer, conjugate diluent, and stop solution.

ACC 00111/0500	Alk. Phos., 500 tests
ACC 00111/1000	Alk. Phos., 1000 tests
ACC 00111/5000	Alk. Phos., 5000 tests

Buffer pack, DAS/Compound

For use with peroxidase DAS and compound direct Reagent Sets. Includes coating buffer, substrate, extraction buffer, PBST wash buffer, conjugate diluent, and stop solution.

ACC 00222/0500	Peroxidase, 500 tests
ACC 00222/1000	Peroxidase, 1000 tests
ACC 00222/5000	Peroxidase, 5000 tests

Buffer pack, Durviz PPV

Buffer for use with Durviz reagent and typing sets

ACC 00432/1000	Alk. Phos., 1000 tests
ACC 00432/1500	Alk. Phos., 1500 tests
ACC 00432/5000	Alk. Phos., 5000 tests

Buffer pack, Grape

ACC 00667/0500	Alk. Phos., 500 tests
ACC 00667/1000	Alk. Phos., 1000 tests
ACC 00667/5000	Alk. Phos., 5000 tests

Buffer pack, Indirect

ACC 00444/0500	Alk. Phos., 500 tests
ACC 00444/1000	Alk. Phos., 1000 tests
ACC 00444/5000	Alk. Phos., 5000 tests

Carbonate coating buffer

ACC 00413/0030	10X concentrate, 30ml
ACC 00413/0075	10X concentrate, 75ml
ACC 00413/0750	10X concentrate, 750ml

Chicken egg albumin

ACC 00310/0050	50g
ACC 00310/0100	100g
ACC 00310/0500	500g

Diethyldithiocarbamate

ACC 00644	DIECA, 10g
-----------	------------

ECI buffer

Conjugate diluent for Alk.Phos. tests

ACC 00120	powder, 10g, makes 300ml
-----------	--------------------------

ECI buffer

ACC 00523/0003	5X concentrate, 3ml
ACC 00523/006.5	5X concentrate, 6.5ml
ACC 00523/0011	5X concentrate, 11ml
ACC 00523/0015	5X concentrate, 15ml
ACC 00523/0030	5X concentrate, 30ml

ECG buffer

ACC 00623/0003	5X concentrate, 3ml
ACC 00623/006.5	5X concentrate, 6.5ml
ACC 00623/0011	5X concentrate, 11ml
ACC 00623/0021	5X concentrate, 21ml
ACC 00623/0102	5X concentrate, 102ml

Extraction buffer, BEB1

ACC 00066/025.4	powder, 25.4g, makes 1 liter
-----------------	------------------------------

Extraction buffer, blueberry

ACC 00880/022.2	powder, 22.2g, makes 500ml
ACC 00880/044.4	powder, 44.4g, makes 1000ml

Extraction buffer, general

ACC 00955/016.5	powder, 16.5g, makes 500 ml
ACC 00955/0033	powder, 33g, makes 1000 ml
ACC 00955/0121	powder, 121g, makes 3667 ml
ACC 00953/016.5	powder, with Tween, 16.5g, makes 500ml
ACC 00953/0033	powder, with Tween, 33g, makes 1000ml
ACC 00953/0121	powder, with Tween, 121 g

Extraction buffer, grape

ACC 00970/049.5	powder, 49.5g, makes 500ml
ACC 00970/098.8	powder, 98.8g, makes 1000ml
ACC 00970/0150	powder, 150g, makes 1500ml

Extraction buffer, indirect

ACC 00980/0025	powder, 25g, makes 1000 ml
ACC 00980/0050	powder, 50g, makes 2000 ml
ACC 00980/0250	powder, 250g, makes 10000 ml
ACC 00985	10X concentrate, 60 ml

Extraction buffer, PEB1

ACC 00890/0050	10X concentrate, 50ml
ACC 00890/0100	10X concentrate, 100ml
ACC 00890/0500	10X concentrate, 500ml

Extraction buffer, SEB1

ACC 00936/0025	Mesh bag containing SEB1, 3ml 25 Bags
ACC 00936/0250	Mesh bag containing SEB1, 3ml 250 Bags
ACC 00936/0050	Powdered buffer, 50 g, makes 2 liters
ACC 00936/0250	Powdered buffer, 250 g, makes 10 liters

Extraction buffer, SEB4

ACC 00958/0025	Mesh bag containing SEB4, 3ml 25 Bags
ACC 00958/0250	Mesh bag containing SEB4, 3ml 250 Bags
ACC 01958/005.7	Powdered buffer, 5.7 g, makes 1 liters
ACC 01958/0050	Powdered buffer, 50 g, makes 8.7 liters
ACC 01958/0250	Powdered buffer, 250g, makes 43.8 liters

MRS component

ACC 00522/0003	5X concentrate, 3ml
ACC 00522/006.5	5X concentrate, 6.5ml
ACC 00522/0011	5X concentrate, 11ml
ACC 00522/0021	5X concentrate, 21ml
ACC 00522/0102	5X concentrate, 102ml

MRS-2 component

ACC 00524/0003	5X concentrate, 3ml
ACC 00524/006.5	5X concentrate, 6.5ml
ACC 00524/0011	5X concentrate, 11ml

Nonfat dried milk

ACC 00420/0005	5g
ACC 00420/0025	25g
ACC 00420/0115	115g

OPD substrate buffer

ACC 00400	1X, ready for use, 60ml
-----------	-------------------------

OPD substrate buffer

ACC 00401/012.5	10X concentrate, 12.5ml
ACC 00401/062.5	10X concentrate, 62.5ml

OPD substrate sticks

ACC 00410/0006	6 count
ACC 00410/0012	12 count
ACC 00410/0030	30 count

PBST wash buffer

ACC 00115	powder, 110 g, makes 10 liters
ACC 00115/0014	powder, 14 kg, makes 1400 liters
ACC 00501	20X concentrate, 50ml

PCR quality water

ACC 00230/0001	1X Ready for use, 1ml
----------------	-----------------------

PNP substrate buffer

ACC 00407/0012	5X concentrate, 12ml
ACC 00407/0025	5X concentrate, 25ml
ACC 00407/0125	5X concentrate, 125ml

PNP substrate tablets

ACC 00404/0012	5 mg each, 12 count
ACC 00404/0025	5 mg each, 25 count
ACC 00404/0100	5 mg each, 100 count
ACC 00404/0200	5 mg each, 200 count

Polyvinylpyrrolidone, PVP-40T

ACC 00600/0100	mw 40,000, powder, 100g
ACC 00600/0500	mw 40,000, powder, 500g
ACC 00600/1000	mw 40,000, powder, 1000g

PZSV EC buffer

Enzyme conjugate diluent buffer for use in PZSV test

ACC 00321	powder, 8.5g, makes 200 ml
-----------	----------------------------

TMB peroxidase substrate solution

ACC 00412/0025	1X, ready for use, 25ml
ACC 00412/0040	1X, ready for use, 40ml
ACC 00412/0060	1X, ready for use, 60ml
ACC 00412/0011	1X, ready for use, 11ml
ACC 00412/0500	1X, ready for use, 500ml

Tween 20

ACC 00528/0015	15ml
ACC 00528/0030	30ml
ACC 00528/0250	250ml

LABORATORY EQUIPMENT AND SUPPLIES

Canister Racks X

Canister holding rack for Kleco 4200/8200

ACC 00733	For Set of 8
-----------	--------------

Canister Sets 802

ACC 00731/0004	Set of 4 for Kleco 4200
ACC 00731/0008	Set of 8 for Kleco 8200

KLECO 4200

Desc OK-supplied with four canisters for processing

ACC 00726	Four channel, multi tissue processor for hard and fibrous tissue.
-----------	---

KLECO 8200

Desc OK-Supplied with eight canisters for processing

ACC 00727	Eight channel, multi tissue processor for hard and fibrous tissue.
-----------	--

Microcentrifuge tubes

ACC 00340/0500	1.5 ml, 500 count
ACC 00340/1000	1.5 ml, 1000 count

Microcentrifuge tube holder

ACC 00375	holds 80 tubes
-----------	----------------

Micro dilution tubes

Tubes for extracted sap collection and storage, 1.2ml, 1000 count

ACC 00635	1.2ml, 1000count
-----------	------------------

Micro dilution tubes, racked

Tubes for extracted sap collection and storage, 1.2ml racked, 10 racks of 96 tubes

ACC 00640	Set of 10 racks with tubes
-----------	----------------------------

Nylon membrane

ACC 00370/0001	pretreated, 96 tests
----------------	----------------------

Pipetter

ACC 00770/0010	10 ul
----------------	-------

Pipet Tips

ACC 00605	Blue, 101-1000 ul, 1000 count
-----------	-------------------------------

Pipet Tips

ACC 00610	Yellow, 1-200 ul, 1000 count
-----------	------------------------------

Pipet Tips

ACC 00620	Large orifice, 1-250 ul, 1000 count
-----------	-------------------------------------

Plate, solid

ACC 00948	High bind
-----------	-----------

Plate, strip

ACC 00950	High bind
-----------	-----------

Plate reader, automated

ACC 00755	
-----------	--

Plate Washer Bottle

1 liter bottle with 8 channel machined head for ELISA plate washing

ACC 00520	
-----------	--

Precoated strip plate

ACC 00940	
-----------	--

Roller Press

Leaf squeezer to extract sap from plant tissue

ACC 00700 for sap extraction

Sample Bag

ACC 00925 4 x 5, 1000 count

Sample bag, mesh

ACC 00930/0100 100 count

ACC 00930/0500 500 count

ACC 00930/1000 1000 count

ACC 00930/5000 5000 count

MB-B0109 1000 count

Sample bag, mesh, with buffer

ACC 00780 containing BEB1, 3ml

ACC 00935 containing GEB, 3ml

ACC 00936 containing SEB1, 3ml

ACC 00958 containing SEB4, 3ml

ACC 00304/0032 32 mesh bags containing 3ml of GEB and 32straw

Sample cup holder tray

8x12 configuration for holding sample cups

ACC 00965

Sample cups

Sample cups for retention of sap from Roller Press

ACC 00960 1200 count

Screwcap vial

ACC 00800 100 count

ACC 00800 500 count

Sealed straws

Straws for transferring sample to ELISA well

ACC 00920 1000 count

Tissue Homogenizer

Hand held device for maceration of tissue in mesh sample bag

ACC 00900 Circular bearings

PCR Primers for Viral Groups

PCR Primer를 이용하여 식물이 특정 viral groups에 감염되었는지 진단할 수 있습니다.

- Sample에서 RNA를 추출하신 뒤, RT-PCR을 통해 virus 감염여부를 알 수 있습니다.
- Sensitivity가 높기 때문에 적은 양의 virus도 detection 할 수 있습니다.

<제품구성>

- Forward and reverse primers (1 tube/enough for 25 reactions)
- Sterile PCR quality water

Cat. No.	Description
PCR 93000	Carlavirus Group
PCR 94000	Carmovirus Group
PCR 95000	Dianthovirus Group
PCR 96100	Illarvirus Group
PCR 95100	Tobamovirus Group
PCR 98100	Tospovirus Group

Carlavirus primers로 detect된 pathogen

- Blueberry scorch virus (BIScV)
- Carnation latent virus (CLV)
- Chrysanthemum virus B (CVB)
- Garlic common latent virus (GarCLV)
- Kalanchoe latent virus (KLV)
- Lily symptomless virus (LSLV)
- Pea streak virus (PeSV)
- Potato latent virus (PotLV)
- Potato virus M (PVM)
- Potato virus S (PVS)
- Shallot latent virus (SLV)

바이러스 검정 RCR Mix Kit

ELISA 방법보다 100배 이상의 감도와 빠른 시간내에 검정할 수 있는 RCR Mix Kit

- RCR Mix Kit는 DNA 수준에서 바이러스를 검정할 수 있는 Kit 입니다.
- 식물, 동물 분야 등에서 사용할 수 있습니다.
- RNA 추출 시약은 물론 Dye 등 모든 시약이 혼합되어 있어 사용하기에 매우 편리합니다.
- Kit에는 RCR Mix와 특정 바이러스에 대한 primer가 들어 있습니다.
- 원하시는 Virus에 대한 primer를 따로 주문하시면 저렴한 가격으로 첨가하여 드립니다.


RCR Mix Kit 특징

- RCR Mix은 One-step master mix로 빠른 시간에 RT-PCR 수행할 수 있습니다.
- Blue dye가 포함되어 있어 RT-PCR이 끝난 후 바로 전기영동을 통해 실험 결과를 확인할 수 있습니다.

RCR Mix Kit components

Reverse Transcriptase,
RNase Inhibitor,
Taq polymerase, dNTP, MgCl₂, Reaction Buffer,
Blue dye

바이러스에 감염된 잎에 Extraction Buffer 넣고 RNA 추출



추출된 RNA와 Primer를 RCR MIX에 첨가
Total Vol. 20ul



1st strand synthesis
amplification

Results of RT-PCR

Contents

	100prep	300prep	500prep
RCR Mix	10 tubes	30 tubes	50 tubes
RCR Buffer	0.5ml x 1	0.5ml x 2	0.5ml x 3
RCR Extraction Buffer	20ml x 1	20ml x 2	20ml x 3
Primer	Virus 종류에 따라 가능한 Primer 공급		

Ordering Information

Cat.No.	Description	Size
KS-R0001-100	RCR Mix Kit – virus 이름	100 tests
KS-R0001-300	RCR Mix Kit – virus 이름	300 tests
KS-R0001-500	RCR Mix Kit – virus 이름	500 tests


Primer 제공되는 Virus

BBWN	CyMV	MNSV	SqMV
CGMMV	CyRSV	ORSV	TMV
CMV	DaMV	PLRV	TUMV
CarLV	LMV	PMMV	WMV
CarMV	LSV	PVY	ZYMV
CarRSV			

이 외에 원하시는 Virus Primer는 별도로 주문 · 제작해 드립니다.

Polyclonal Antibody Production

- Sample을 주시면 원하시는 polyclonal antibody 를 제작해드립니다.
- 생산된 antibody는 ELISA를 수행하여, OD 값이 1.0 이상인 경우, 성공으로 간주합니다.




Antibody Conjugation Service

- Antibody 에 Alkaline Phosphatase를 conjugation 시켜드림으로 ELISA asssay를 비롯한 immunoassay 실험을 쉽게 수행할 수 있도록 해드립니다.

Service 1) Antibody Conjugation Service

보유하고 계신 Crude Serum 또는 정제된 Antibody를 보내주시면, Antibody를 purify한 후, alkaline phosphatase를 Conjugation 시킵니다. (기간 : 1주일)



Service 2) ELISA Test

Conjugation된 Antibody를 이용해서 ELISA test를 해드립니다.
(최소 10 testwell 이상/antibody 1 종류당)

Animal Disease Diagnoses



Mammalian

Miscellaneous

Fish

Anabolics

Antibiotics

Sulfur drugs



1. Mammalian Diseases detection

ELISA kits

Cat. No.	Product	tests/kit
----------	---------	-----------

Antibody Detection kit

MB-K 027	BHV-1 (IBR) threshold	96
	BHV-1 (IBR) threshold (10 kits)	10 X 96
MB-K 132	BHV-1 (IBR) seroconversion	96
MB-K 004	BVDV antibody	96
MB-K 061	BRSV	96
MB-K 062	Parainfluenza 3 (PI3)	96
MB-K 063	Adenovirus 3 (ADENO 3)	96
MB-K 162	<i>Mycoplasma bovis</i>	96
MB-K 028	Respiratory pentavalent kit (BGW B-292) (BHV-1, BVDV, BRSV, PI3, ADENO 3)	32 X 5
	Respiratory trivalent RPA NEW (BRSV, PI3, ADENO 3)	
MB-K 146	BHV-4	96
MB-K 072	Abortion kit(BHV-1, BVDV, BHV-4)	48 X 3
MB-K 137	<i>Neospora caninum</i>	96
MB-K 064	<i>Fasciola Hepatica</i>	96
MB-K 126	Rotavirus (sero blocking)	192
MB-K 127	Coronavirus (sero blocking)	192
MB-K 128	<i>E. coli F5(k99)</i> (sero blocking)	192

Antigen Detection kit

MB-K 067	Rotavirus	96
MB-K 068	Coronavirus	96
MB-K 069	<i>E. coli F5(k99)</i>	96
MB-K 070	<i>Cryptosporidium parvum</i>	96
MB-K0 71	Digestive kit (BGW B-290) (Rota, Corona, F5(k99), Crypto)	24 X 4
	Duo Digestive (BGW B-242) (Rota, Corona)	
MB-K 083		48 X 2

Clostridium perfringens antigen detection

MB-K 084	alpha toxin	96
----------	-------------	----

MB-K 094	beta toxin	96
MB-K 085	epsilon toxin	96
MB-K 086	<i>Clostridium perfringens</i>	96
MB-K 095	Enterotoxaemia (toxigenicity $\alpha, \beta, \epsilon, C$. perf.)	24

Antigen Detection

MB-K 189	Pulmotest BHV-1	96
MB-K 190	Pulmotest BVDV	96
MB-K 187	Pulmotest PI3	96
MB-K 188	Pulmotest BRSV	96
MB-K 121	Pulmotest <i>Mycoplasma bovis</i>	96
MB-K 157	Pulmotest ADENO 3	96
MB-K 191	Pulmotest Tetra(BHV-1, BVDV, BRSV, PI3)	24 X 4

Quick test kits

Strips (quick tests)

MB-K 152	ROTAVIRUS	10
MB-K 153	CORONAVIRUS	10
MB-K 154	E. COLI F5(k99)	10
MB-K 155	CRYPTOSPORIDIUM PARVUM	10
MB-K 156	Tetrastrips (Rota, Corona, F5(k99), Crypto)	4 X 5

Clostridium perfringens Strips (quick tests)

MB-K 170	Alpha Toxin	10
MB-K 175	Beta Toxin	10
MB-K 176	Epsilon Toxin	10

Antibodies (dilution 1/20)

Monoclonal antibodies fluorescein conjugated

MB 039	FITC anti-ADENO 3	0.5mL
MB 026	FITC anti-BHV-1 (IBR)	"
MB 286	FITC anti-BHV-4	"
MB 032	FITC anti-BRSV	"
MB 316	FITC anti-BVDV	"
MB 023	FITC anti-Coronavirus	"
MB 073	FITC anti- <i>Cryptosporidium parvum</i>	"
MB 272	FITC anti-CSFV (Hog Cholera) BFAV B 360	"
MB 062	FITC anti-E. coli F5(k99)	"

Antibodies (dilution 1/20)		
Monoclonal antibodies fluorescein conjugated		
MB 245	FITC anti-EAV (Equine Arteritis Virus)	"
MB 274	FITC anti- <i>Mycoplasma bovis</i>	"
MB 030	FITC anti-PI3	"
MB 268	FITC anti-PRRSV	"
MB 265	FITC anti-Pseudorabies (Aujeszky) SHV-1	"
MB 020	FITC anti-Rotavirus	"
MB 318	FITC anti-Neospora caninum	"
MB 314	FITC anti-RHD	"
MB 325	FITC anti-Lawsonia intracellularis	"
Unconjugated monoclonal antibodies (+/- 100 µg/mL)		
MB 292	anti-ADENO 3	0.5mL
MB 289	anti-BHV-1 (IBR)	"
MB 279	anti-BHV-4	"
MB 296	anti-Bovine IgG1	"
MB 300	anti-Bovine IgM	"
MB 031	anti-BRSV	"
MB 295	anti-BVDV (NS3)	"
MB 288	anti-Coronavirus	"
MB 294	anti-Cryptosporidium parvum	"
MB 275	anti-CSFV BFAV/KSP/D13/01	"
MB 293	anti-E. coli F5(k99)	"
MB 278	anti-EAV	"
MB 297	anti- <i>Mycoplasma bovis</i>	"
MB 290	anti-PI3	"
MB 276	anti-PRRSV	"
MB 277	anti-Pseudorabies (Aujeszky) SHV-1	"
MB 287	anti-Rotavirus	"
MB 319	anti-Neospora caninum	"
MB 315	anti-RHD	"
MB 323	anti-Lawsonia intracellularis	"

2. Miscellaneous		
Cat. No.	Product	tests/kit
ELISA Kit		
MB-K 093	Bovine Gamma Interferon	80 / 160
MB-K 165	Dosage of bovine Immunoglobulin	96 x 2
MB-K 082	BVDV Antigen (NS3 sur leucocytes)	192
MB-K 158	Bovine Lactoferrin	96 X 2

3. Fish Diseases detection		
Cat. No.	Product	tests/kit
Immunofluorescence assay - FLUO kit		
MB-K 006	VHS (Viral Haemorrhagic Septicaemia)	24 X 2
MB-K 008	IHN (Infectious Haematopoietic Necrosis)	24 X 2
MB-K 010	IPN (Infectious Pancreatic Necrosis)	24 X 2
MB-K 012	SVC (Spring Viraemia of Carp)	24 X 2
Immunoperoxidase assay - PEROX kit		
MB-K 007	VHS	24 X 2
MB-K 009	IHN	24 X 2
MB-K 011	IPN	24 X 2
MB-K 013	SVC	24 X 2
Unconjugated monoclonal antibodies +/- 100 µg		
MB 282	anti-VHS (IP5B11) anti-N protein	0.5mL
MB 285	anti-IHN (136-3) anti-N protein	"
MB 312	anti-IHN (34D3) anti-N protein	"
MB 050	IPN (Infectious Pancreatic Necrosis)	"
MB 053	SVC (Spring Viraemia of Carp)	"
Monoclonal antibody fluorescein conjugated		
MB 281	FITC anti-VHS (IP5B11)	0.5 mL
MB 284	FITC anti-IHN (136-3)	"
MB 311	FITC anti-IHN (34D3)	"
ELISA Antigen Detection kit		
MB-K 159	<i>Aeromonas salmonicida</i>	96

4. ANABOLICS detection kits

Cat. No.	Product	tests/kit
Radioimmuno assay (Plasma, urine, feces)		
MB-K 014	3H-Ethynylestradiol	500
MB-K 015	3H-Methyltestosterone	500
MB-K 016	3H-Nortestosterone	500
MB-K 017	3H-Zeranol	500
MB-K 018	3H-Clenbuterol	500
MB-K 019	3H-β2 Agonist	500
MB-K 023	3H-17β-Estradiol	500
MB-K 025	3H-Progesterone	500
MB-K 080	3H-testosterone	500
MB-K 163	3H-Chloramphenicol	500
ELISA kit (serum)		
MB-K 022	19-Nortestosterone	96
MB-K 074	Methyltestosterone	96
MB-K 075	Clenbuterol	96
MB-K 134	β2 Agonist overnight	96
MB-K 076/1	β2 Agonist (quantitative)	96
MB-K 076/2	β2 Agonist (qualitative)	96
MB-K 076/3	β2 Agonist (tube)	25
MB-K 107	α-Trenbolone	96
MB-K 108	β-Trenbolone	96
MB-K 089	Zeranol	96
MB-K 090	Isoxsuprine	96
MB-K 091	Diethylstilbestrol (DES)	96
MB-K 097	Clostebol	96
MB-K 098	Generic Corticosteroids	96
MB-K 099	Dexamethasone (on request)	96
MB-K 109	Triamcinolone acetonide	96
MB-K 110	Generic Promazine	96
MB-K 111	Azaperone - Azaperol	96
MB-K 112	Carazolol	96
MB-K 149	Chloramphenicol	96
MB-K 102	Stanozolol	96

5. ANABOLICS

Cat. No.	Product	tests/kit
Immunoaffinity Chromatography gels		
MB 261	for Chloramphenicol	1 mL
MB 257	for Clenbuterol	1 mL
MB 259	for Dexamethasone	1 mL
MB 246	for Diethylstilbestrol	1 mL
MB 247	for Ethynylestradiol	1 mL
MB 258	for Generic Corticosteroids	1 mL
MB 256	for Generic β2 Agonists	1 mL
MB 248	for Methyltestosterone	1 mL
MB 249	for Methyltestosterone (high capacity)	1 mL
MB 250	for Nortestosterone	1 mL
MB 251	for Stanozolol	1 mL
MB 255	for Steroids multiresidues (on request)	3 mL
MB 252	for testosterone	1 mL
MB 253	for Trenbolone	1 mL
MB 260	for Triamcinolone acetonide	1 mL
MB 254	for Zeranol	1 mL

6. ANTIMBTICS and SULFUR DRUGS

Detection kits		
Cat. No.	Product	tests/kit
MB-K 033/1	Total Antibiotic (plate) milk & meat	288
MB-K 033/2	Total Antibiotic (tube) milk & meat	96
MB-K 033/3	Total Antibiotic & Sulphonamide (plate) milk	288
MB-K 033/4	Total Antibiotic & Sulphonamide (tube) milk	96